


# Low Profile/ Trans-Cuff SLAP Repair

with the JuggerKnot® Soft Anchor-1.4/1.5 mm  
with Percutaneous Instrumentation

Surgical Technique





# Table of Contents

- Patient Positioning ..... 2
- Prepare the Surface ..... 3
- Placement of Percutaneous Guide ..... 3
- Pilot Hole and Anchor Placement ..... 5
- Anchor Deployment..... 6
- Suture Retrieval ..... 6
- Knot Tying ..... 8
- Ordering Information ..... 9
- Indications for Use ..... 10
- Contraindications ..... 10

## Surgical Technique



Figure 1




Figure 2

### Patient Positioning and Portal Placement

Place the patient in either the beach chair or lateral decubitus position depending upon surgeon preference (Figure 1). Insert a 30° arthroscope through the posterior portal. Carefully assess all anatomy including any associated chondral, rotator cuff and biceps abnormality. Place a 7 mm AquaLoc® cannula appropriately through the rotator interval either from an outside in or inside out technique.

Identify the labral tear. The tear will either be posterior only, anterior only, or anterior to posterior (Figure 2).




Figure 3




Figure 4

## Prepare the Surface

Utilize either a rasp or a curved shaver through the anterior portal to prepare the bone bed until bleeding cortical bone is identified at the repair site.

## Placement of Percutaneous Guide

To repair the superior labral tear, create a transcutaneous portal and incision to place the JuggerKnot percutaneous guide. To identify the appropriate angle, utilize a spinal needle through an anterolateral percutaneous area, adjacent to the anterolateral acromion (Figure 3). Visualize the needle intra-articularly to allow a trans-cuff angle into the appropriate glenoid bed. Place a small stab incision, smaller than a portal incision, through the skin at the location of the spinal needle.

Insert the rigid trocar for the percutaneous guide through the skin incision at the same angle as previously identified by the spinal needle (Figure 4).




Figure 5




Figure 6

### Placement of Percutaneous Guide (cont.)

Insert the percutaneous guide over the rigid trocar. Utilize a rotating motion back and forth to penetrate through both subcutaneous tissue and particularly rotator cuff tendon into the intra-articular space (Figure 5). Remove the rigid trocar.

Place the guide on the prepared superior glenoid area (Figure 6).




Figure 7




Figure 8

## Pilot Hole and Anchor Placement

Chuck the drill bit to the proximal laser etch line. Advance the drill under power until the drill bottoms out on the top of the guide (Figure 7). Remove the drill. Utilize a 1.5 mm drill if a JuggerKnot 1.5 mm is selected.

ⓘ **Note:** Be sure to maintain the precise position of the guide, including the angle, over the drill hole.

Insert the JuggerKnot 1.4 mm or 1.5 mm anchor through the guide until seated at the entry of the pilot hole. Utilize a small mallet to insert the JuggerKnot anchor in the bone. Advance the inserter until it bottoms out on the top of the guide handle (Figure 8).

ⓘ **Note:** If utilizing a JuggerKnot 1.5 mm anchor, a very firm/tight fit will be experienced in passing the anchor through the percutaneous guide.




Figure 9

Figure 9a




Figure 10

## Anchor Deployment

Unscrew the white luer-lock cap to release the sutures from the inserter. Remove the inserter by pulling straight back on the inserter handle. Remove the guide and seat the anchor by lightly pulling back on both sutures exiting through the skin (Figure 9 & 9a). Check to verify the sutures slide within the anchor.

## Suture Retrieval

Through the posterior portal utilize a suture grasper to retrieve both limbs of the suture (Figure 10).






Figure 11



Figure 12




Figure 13

## Suture Retrieval (cont.)

Next, insert an angled 90° Up suture retriever through the posterior portal. Penetrate the superior portion of the labrum downward and advance the Nitinol wire into the joint (Figure 11).

Continue to advance as much of the Nitinol wire into the joint as possible. Slowly remove the suture retriever from the joint while simultaneously continuing to advance the Nitinol wire. Take care to ensure that the opening of the wire does not retreat through the tissue\* (Figure 12).

Through the posterior portal retrieve the wire using a suture grasper. Pass one strand of the MaxBraid™ suture through the wire loop and pull the Nitinol wire to pass the suture (Figure 13).

ⓘ **Note:** Due to this tear being posterior, the scope has been moved to the anterior portal to allow for easier knot tying through the posterior portal.

\* A two portal suture shuttle technique can also be used with the creation of an additional portal.




Figure 14

## Knot Tying

### Posterior/Superior Repair

For posterior/superior labral repairs, retrieve the sutures either anteriorly or posteriorly. If the decision is made to tie a knot posteriorly, insert the arthroscope through the cannula anteriorly and place a separate cannula from an outside in technique using the same standard posterior portal. Retrieve both sutures through either cannula with the opposite cannula acting as the looking portal. Use a standard knot-tying technique to secure the labrum. Cut the remaining suture tails with a Slotted MaxCutter™. Fixation is complete.

### Anterior/Superior Repair

Utilize a standard knot-tying technique with the cannula above the biceps to place the appropriate knot. Cut the remaining suture tails with a Slotted MaxCutter (Figure 14). Fixation is complete.

## Ordering Information

### Implants

<b>Part Number</b>	<b>Description</b>
912030	1.4 mm JuggerKnot Single Loaded
912010	1.4 mm JuggerKnot Package of 10
912031	1.5 mm JuggerKnot Single Loaded
912015	1.5 mm JuggerKnot Package of 10

### Instrumentation

<b>Part Number</b>	<b>Description</b>
912140C	1.4 mm JuggerKnot Curved Guide Disposable Kit with Centering Sleeve
912141C	1.5 mm JuggerKnot Curved Guide Disposable Kit with Centering Sleeve
912040P	1.4/1.5 mm JuggerKnot Percutaneous Kit

## **INDICATIONS FOR USE**

The JuggerKnot Soft Anchors are intended to be used for soft tissue to bone fixation for the following indications:

### **Shoulder**

Bankart lesion repair  
SLAP lesion repair  
Acromio-clavicular repair  
Capsular shift/capsulolabral reconstruction  
Deltoid repair  
Rotator cuff tear repair  
Biceps tenodesis

### **Foot and Ankle**

Medial/lateral repair and reconstruction  
Mid-and forefoot repair  
Hallux valgus reconstruction  
Metatarsal ligament/tendon repair or reconstruction  
Achilles tendon repair

### **Elbow**

Ulnar or radial collateral ligament reconstruction  
Lateral epicondylitis repair  
Biceps tendon reattachment

### **Knee**

Extra-capsular repair: MCL, LCL, and posterior oblique ligament  
Iliotibial band tenodesis  
Patellar tendon repair  
VMO advancement  
Joint capsule closure

### **Hand and Wrist**

Collateral ligament repair  
Scapholunate ligament reconstruction  
Tendon transfers in phalanx  
Volar plate reconstruction

### **Hip**

Acetabular labral repair

## **CONTRAINDICATIONS**

1. Infection.
2. Patient conditions including blood supply limitations and insufficient quantity or quality of bone or soft tissue.
3. Patients with mental or neurologic conditions who are unwilling or incapable of following postoperative care instructions or patients who are otherwise unwilling or incapable of doing so.
4. Foreign body sensitivity. Where material sensitivity is suspected, testing is to be completed prior to implantation of the device.








This material is intended for healthcare professionals and the Zimmer Biomet sales force. Distribution to any other recipient is prohibited. For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert or contact your local representative; visit [www.zimmerbiomet.com](http://www.zimmerbiomet.com) for additional product information. All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

Check for country product clearances and reference product specific instructions for use. Zimmer Biomet does not practice medicine. This technique was developed in conjunction with health care professional. This document is intended for surgeons and is not intended for laypersons. Each surgeon should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient, and this information does not purport to replace the comprehensive training surgeons have received. As with all surgical procedures, the technique used in each case will depend on the surgeon's medical judgment as the best treatment for each patient. Results will vary based on health, weight, activity and other variables. Not all patients are candidates for this product and/or procedure. Caution: Federal (USA) law restricts this device to sale by or on the order of a surgeon. Rx only.

©2016, 2022 Zimmer Biomet



0394.2-GLBL-en-Issue Date 2022-01-21

 **Legal Manufacturer**  
Biomet Sports Medicine  
P.O. Box 587  
56 E. Bell Drive  
Warsaw, Indiana 46581-0587  
USA

[www.zimmerbiomet.com](http://www.zimmerbiomet.com)

 2797

CE mark on a surgical technique is not valid unless there is a CE mark on the product label.