


Synovasure® RISC™ Panel





Causes of TKA Revision Surgery¹

THE BURDEN OF REVISION: 3% OF KNEE REPLACEMENTS ARE REVISED²

Total knee arthroplasty has become one of the most common orthopedic procedures for treatment of knee, with over 1 million total knee replacements performed each year in the United States, and is expected to grow to nearly 3.4 million by 2040³. Though a highly successful surgical procedure, nearly 3% of knee replacements are revised within the first 5 years for reasons such as infection, loosening, pain and instability². Ultimately, revisions place a burden on patients and the healthcare system through the potential increase in costs, readmissions and morbidity². If patients with a higher risk of post-operative complications are identified before surgery, proper treatment plans may be able to be developed to mitigate the risks.

A SINGLE ASPIRATION FOR A MORE INFORMED TREATMENT PLAN

Getting an accurate arthritis diagnosis is critical in optimizing a patient treatment plan to maximize satisfaction and minimize postoperative risks. Historically, individual laboratory tests have been utilized to diagnosis or rule out specific arthritis conditions through serum and joint fluid samples.

The Synovasure® Relative Inflammatory Status Classification (RISC) Panel is the first of its kind panel test offering a one stop, comprehensive panel of synovial fluid tests to aid in diagnosing different arthritis conditions that may correspond to relative level of inflammation and post-operative complications. The panel includes biomarker tests for:

- Osteoarthritis (OA)
- Rheumatoid Arthritis (RA)
- Crystalline Arthritis
- Native Septic Arthritis (NSA)

An accurate arthritis diagnosis helps identify which treatment options are best for the patient, including implant selection, venue, and pre- and post-operative care plans.

KEY FEATURES



Single, synovial fluid test panel for common arthritis types



Streamlined submission process



Every specimen tested for integrity



Easy-to-read results report



First and only test panel for diagnosis of OA



Quick turnaround within 72 hours

OSTEOARTHRITIS

First and Only Laboratory Test to Aid in OA Diagnosis

Historically, laboratory tests were used to exclude OA from a diagnosis. The Synovasure® RISC™ Panel includes the first test specifically designed for OA diagnosis.

Dual Biomarker Evaluation

Differentiate between primary and secondary OA through assessing levels of:

- Cartilage Oligomeric Matrix Protein (COMP) – Cartilage breakdown biomarker^{1,4}
- Interleukin-8 (IL-8) Inflammatory biomarker^{5,6}

Highly Accurate for OA Diagnosis

COMP/IL-8 ratio demonstrates high sensitivity and specificity⁷





RHEUMATOID ARTHRITIS

Testing at the Source

Synovial fluid biomarker testing at the pain source using established serum biomarkers

- Anti-Cyclic Citrullinated Peptide (Anti-CCP) Biomarker associated with presence of RA⁸
- Rheumatoid Factor (RF) Biomarker associated with presence of an autoimmune disease⁹

Result Algorithm for Diagnosis Recommendation

Results for Anti-CCP and RF combined to determine if pain potentially caused by RA or other auto-immune condition.

Performance Backed by Data⁷

- Capable of discriminating between RA and Non-RA patients
- Effective "rule-in" test 100% Specificity

CRYSTALLINE ARTHRITIS

Crystal analysis is performed as part of the Synovasure RISC Panel to detect monosodium urate (MSU) and calcium pyrophosphate dihydrate (CPPD) crystals using high-powered microscopy for gout and CPPD disease.

NATIVE SEPTIC ARTHRITIS

Evaluation of White Blood Cells (WBC) in Fluid

- Automated, high-performance cell count for total cells and breakdown of cell types present.
- A 2020 study shows nearly 4.4% of automated cells counts in native knee specimens resulted in a false positive¹⁰; therefore, all cell counts >3,000 are confirmed with a manual count.

If final reported WBC count is >3,000 cells/ μ L or neutrophil percentage is >70%, specimens are reflexed to Synovasure Alpha Defensin for NSA (alpha defensin and lactate), Synovasure Microbial Identification and Fluid Culture.

Synovasure Alpha Defensin for NSA

Laboratory-developed test (LDT) for alpha defensin biomarker and lactate. Alpha defensin demonstrates high level of accuracy in identifying infection^{11*} and is unaffected by:

- Prior Antibiotic Administration^{11,12}
- History of Inflammatory Conditions^{11,13}
- Type/Virulence of Organism^{11,14}

Synovasure Microbial Identification

A novel, antigen bead-based test for 24-hour identification of potential infection causing organism including:

	SENSITIVITY ¹⁵	SPECIFICITY ¹⁵	
Staphylococcus species	94%	99%	
Enterococcus species	97%	99%	
Candida species	90%	99%	
C. acnes	Due to difficulty culturing C. acnes	Due to difficulty culturing C. acnes, performance data not available	

The Synovasure Microbial Identification test is highly accurate and identifies more than 54% of culture-negative PJI specimens.¹⁶

Fluid Culture

Aerobic and anaerobic culture performed to aid in the diagnosis and treatment of infection. All fluid cultures are held for seven (7) days.


^{*}Validation and performance of alpha defensin is based on cases of PJI and has not been defined for NSA due to no comparative "Gold Standard" criterion.

STREAMLINED PROCESS

Submitting synovial fluid to CD Laboratories uses a simple, streamlined process from aspiration to result reporting. Results are provided within 72 hours of receipt at CD Laboratories.



ASPIRATE JOINT




TRANSFER FLUID TO PROVIDED TUBES



PACKAGE TUBES AND REQUISITION IN SHIPPER



OVERNIGHT SHIPMENT TO CD LABORATORIES



TESTS PERFORMED AND RESULTS SENT TO PROVIDER

MANAGING POST-OPERATIVE TOTAL KNEE ARTHROPLASTY COMPLICATIONS

Zimmer Biomet is your partner for assessing and preventing the potential for complications and detecting and treating complications that may arise.





Patient Specific Solutions



Limb Salvage



Re-implantation

References

- Bellucci F., et al. Synovial Fluid Levels of Bradykinin Correlate with Biochemical Markers for Cartilage Degradation and Inflammation in Knee Osteoarthritis. Osteoarthritis and Cartilage. 21(11) 1774-1780, 2013.
- Bhandari M, et al. Clinical and Economic Burden of Revision Knee Arthroplasty. Clin Med Insights Arthritis Musculoskelet Disord. 5:89-94, 2012.
- Singh J, et al. Rates of Total Joint Replacement in the United States: Future Projections to 2020–2040 Using the National Inpatient Sample. Journal of Rheumatology. doi:10.3899/jrheum.170990, 2019
- Arellano R., et al. Cartilage Oligomeric Matrix Protein Levels in Synovial Fluid in Patients with Primary Knee Osteoarthritis and Healthy Controls: A Preliminary Comparative Analysis with Serum Cartilage Oligomeric Matrix Protein. Archives of Rheumatology. 32(3):189-196, 2017
- Partsch G., et al. Highly Increased Levels of Tumor Necrosis Factoralpha and Other Proinflammatory Cytokines in Psoriatic Arthritis Synovial Fluid. The Journal of Rheumatology. 24(3) 518-523, 1997.
- Harada A., et al. Essential Involvement of Interleukin-8 (IL-8) in Acute Inflammation. Journal of Leukocyte Biology, 56(5) 559-564, 1994
- 12514-009 Technical Report: Performance of Synovasure® Differential Diagnosis Arthritis Panel Diagnostic Accuracy Test, 2020.
- 8. Freeman, J. RA and Anti-CCP: What is the Purpose of an Anti-CCP Test? https://www.rheumatoidarthritis.org/ra/diagnosis/anticcp/. Accessed 1 Sep 2021.

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for health care professionals. Distribution to any other recipient is prohibited.

For indications, contraindications, warnings, precautions, potential adverse effects and patient counselling information, see the package insert or contact your local representative; visit www.zimmerbiomet.com for additional product information.

Zimmer Biomet does not practice medicine. Each physician should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient, and this information does not purport to replace the comprehensive training physicians have received.

This test has been developed for use with synovial fluid only. The use of this test kit with any other specimen type may lead to inaccurate test results.

The Synovasure RISC Panel is not available in the state of New York.

© 2022 Zimmer Biomet

- Street, T. Rheumatoid Factor. https://emedicine.medscape.com/ article/2087091-overview. Accessed 1 Sep 2021.
- Deirmengian C., et al. False-Positive Automated Synovial Fluid White Blood Cell Counting is a Concern for Both Hip and Knee Arthroplasty Aspirations. Journal of Arthroplasty. 35(6S):S304-S307, 2020.
- Deirmengian C., et al. Combined Measurement of Synovial Fluid a-Defensin and C-Reactive Protein Levels: Highly Accurate for Diagnosing Periprosthetic Joint Infection. Journal of Bone and Joint Surgery. 96(17):1439-45, 2014.
- Shahi A., et al. The Alpha-Defensin Test for Periprosthetic Joint Infection is Not Affected by Prior Antibiotic Administration. Clinical Orthopedics and Related Research. 474(7):1610-5, 2016.
- Miyamae Y., et al. Diagnostic Accuracy of the Alpha-Defensin Test for Periprosthetic Joint Infection in Patients with Inflammatory Diseases. Journal of Arthroplasty. 34(8):1767-1771, 2019
- Deirmengian C., et al. The C-reactive Protein May Not Detect Infections Cause by Less-Virulent Organisms. Journal of Arthroplasty. 31(9 Suppl):152-5, 2016.
- M40022B Synovasure Microbial Identification Instructions for Use (IFU), 2019.
- Deirmengian C., et al. Diagnostic Performance of the Synovasure MID Test [abstract]. In: Proceedings of the MSIS Annual Meeting; 2019 Aug 2-3; New York, NY, USA.

