



Surgical Technique

Avitus® ArchiMIS™

MIS Precision
Autograft Delivery



Indications and Contraindications

Indications For Use

The Avitus® ArchiMIS™ Precision Autograft Delivery system is intended:

- 1.** To be used for the delivery of autograft or hydrated allograft bone graft material to an orthopaedic surgical site.


Contraindications


- | | |
|--|--|
| <p>1. Patients who are not a candidate for bone graft delivery.</p> | <p>2. If using non-autologous bone graft, refer to packaging insert of the specific bone graft material being used.</p> |
| <p>3. The device is intended to be used by a physician familiar with the possible side effects, typical findings, limitations, indications, and contraindications of bone graft delivery.</p> | |



Avitus[®] ArchiMIS[™]


MIS Precision Autograft Delivery





**A percutaneous controlled delivery
experience for autologous
gold-standard bone grafting.**

Avitus ArchiMIS System Contents



Ordering Information

Part#	Description	Volume / Size
AD-1CH-4.0-120	Avitus ArchiMIS Autograft Delivery	1cc of Graft
		ø3.8mm x 106mm - Straight Delivery Cannula
		ø4.0mm x 120mm - Angled Delivery Cannula

Technique

1. How to get a proper scoop of autograft


Use the Bone Scoop to obtain a level scoop of graft (fig. 1).



2. Remove excess graft

Use the Pack Tool to remove any excess bone graft until the scoop is sub-flush (fig. 2a, 2b).


WARNING: Do not use more than 1 sub-flush scoop per graft cartridge



Technique


3. Packing the Graft Cartridge

Delivery Channel for
packing bone graft




Distribute graft along channel.

Distribute the bone graft along the length of the delivery channel from the scoop while aligning along the delivery channel with the pack tool (fig. 3b).



Compress the graft into channel.


Use the pack tool to compress the bone graft down into the delivery channel with the pack tool (fig. 3c).



Technique


Check for proper pack.

Ensure that the graft is compressed into each channel without residual material above the delivery channel (fig. 3d).



4. Close the Graft Cartridge


Properly packed autograft cartridge.



Technique

5. Assembling the cartridge


Connect the graft cartridge to the desired delivery cannula (Fig. 5a, 5b)



Technique


6. To prime the cannula

Push the delivery plunger into the cartridge until the threads meet (Fig. 7a). Device is “primed” right before the threads begin to engage. DO NOT rotate Delivery Plunger unless ready to deploy bone graft (Fig. 7b).



7. Advance the autograft

Turn the threaded delivery plunger clockwise to advance autograft into the intended delivery site (Fig. 8).



8. Reloading the cartridge

Repeat steps 3-6 with the replacement graft cartridge and reconnect the desired delivery cannula.

Press the sideclips to disconnect the empty graft cartridge and reconnect a newly packed graft cartridge.

TROUBLESHOOTING STEPS:

1. Purging the cannula between reloads to get every last bit of bone out of the cannula.
2. If doing an OCD delivery, slide the angled cannula opening along the rim of the defect to dislodge any remaining bone chunks before removing device.



Harvest. Debride.
Deliver.

Avitus®

DragonWing™ | Bone Harvester | Bone Debrider | ArchiMIS™



 **ZIMMER BIOMET**
Moving You Forward.™



Harvest. Debride.
Deliver.

Avitus®

DragonWing™ | Bone Harvester | Bone Debrider | ArchiMIS™



This material is intended for health care professionals. Distribution to any other recipient is prohibited. For indications, contraindications, warnings, precautions, potential adverse effects and patient counselling information, see the instructions for use or contact your local representative; visit www.zimmerbiomet.com for additional product information.

Zimmer Biomet does not practice medicine. This technique was developed in conjunction with a health care professional. This document is intended for surgeons and is not intended for laypersons. Each surgeon should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient, and this information does not purport to replace the comprehensive training surgeons have received. As with all surgical procedures, the product/technique used in each case will depend on the surgeon's medical judgment as the best treatment for each patient. Results will vary based on health, weight, activity, and other variables. Not all patients are candidates for this product and/or procedure. Caution: Federal (USA) law restricts this device to sale by or on the order of a surgeon. Rx only.

All content herein is protected by copyright, trademarks, and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated, or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

©2025 Zimmer Biomet. All rights reserved.



ST-0808-03-EN-Issue Date YYYY-MM-DD
VV-10989



Legal Manufacturer

Biomet Orthopedics
56 East Bell Drive
PO Box 587
Warsaw, IN 46581
www.zimmerbiomet.com