

**Surgical Technique**

# Avitus® Bone Debrider





# Table of Contents

---

<b>System Overview/Features</b> .....	<b>03</b>
<b>Device Setup/Contents</b> .....	<b>04</b>
General Debridement Technique .....	07
Example Debridement Applications .....	08
Estimating Specimen Volume.....	09
Specimen Retrieval.....	10
Device Re-Assembly Steps.....	12
<b>Troubleshooting/Declogging</b> .....	<b>14</b>
<b>Ordering Information</b> .....	<b>16</b>

---

# System Overview

---

## Indications For Use

The Avitus Bone Debrider is intended:

**1.** To harvest cancellous bone and marrow

**2.** To debride and capture infected, necrotic or diseased cancellous bone (e.g. osteomyelitis, cancellous bone tumors)

## Contraindications

**1.** Osteoporosis, or other disorder that diminishes the quality of bone tissue.

**2.** Active infection in or around the donor site (unless using device to remove bone infection).

**3.** Bone anatomy in which the device cannot fit.

**4.** Previous donor site harvest.

**5.** Patients who are not a candidate for autograft.

## Summary

The Avitus® Bone Debrider is a closed-capture suction curettage system that can be used for the debridement and capture of infected, necrotic, or diseased cancellous bone. Specimen collects inside the Filter Insert for easy retrieval. It is a plug-and-play system that connects to standard operating room suction sources.

## System Overview



## Device Setup - Kit Content



## Device Setup - Connect to Suction

<p>1. Connect Avitus® Bone Debrider to the suction source with a <b>dedicated suction line that is not shared by other suction devices in the room.</b></p>	<p>2. Keep Avitus® Bone Debrider easily accessible so that it is readily available for debridement.</p>
<p>3. Turn suction source to HIGH/MAX.</p>	<p>4. Prepare specimen cups on the back table (not included).</p>

# General Debridement Technique - Pearls

## Pearls

### Scrape the Canals

Drag the angled cutting tip blade along surfaces to scrape. Limited to canals 13mm in diameter or greater (Fig. 1).


### Work the Clock

For thorough debridement, scrape the length of the canal at each circumferential increment (Fig 2).

E.g. Scrape along the length of the canal with the angled blade facing the 12 o'clock position around the canal. Rotate the handle to align the blade to the 2 o'clock position and scrape again. Repeat until satisfied with debridement.


### Empty the Device when FULL

Suction power may be reduced when device reaches capacity. Empty contents and reconnect to suction to resume debridement (Fig. 3).




## Example Debridement Applications

---



**TKA – Femur & Tibia  
Canal Debridement**




**THA – Hip Canal  
Debridement**



## Example Debridement Applications


---



## Estimating Specimen Volume

---

When the Filter Insert is removed from the Handle of the Debrider, specimen volume can be estimated in the Bone Receptacle as follows:



## Specimen Retrieval

---

### Remove the Cap

Hold the Debrider horizontally over the Collection Cup and unscrew the Cap. Clean off any contents stuck to the Cap.


**NOTE:** Take care to stop the Filter Insert from sliding out of the handle at this point.



### Remove the Filter Insert

Ensure that the handle is held horizontally with the logo facing upwards. Carefully remove the Filter Insert.

**NOTE:** If the logo is not facing upwards, the specimen will fall out of the Filter Insert as it is removed from the handle.



# Specimen Retrieval

---

## Retrieve Bone from the Filter Insert Receptacle



Use the rounded end of the Plunger accessory (included) to remove the specimen from the Filter Insert receptacle into the specimen cup. Take care to retrieve every chunk.

# Specimen Retrieval

## Retrieve Remaining Debridement Contents from the Debrider Handle


Use the Spatula accessory (Fig. 1) to remove any remaining debridement contents from the Debrider handle.

Insert the Spatula into the Debrider handle and press gently against the wall while rotating the Debrider handle to scrape any remaining cancellous graft into the Spatula (Fig. 2).



**NOTE:** Spatula accessory included in kit.

Applying gentle pressure against the flat side of the Debrider handle, drag the Spatula out of the handle to retrieve the remaining cancellous bone (Fig.3).





# Device Re-Assembly

## Re-Insert the Filter Insert

### Orientation

Hold the Filter Insert with the rectangular opening of the receptacle facing upward. Similarly, hold the Debrider with the Avitus logo facing upwards.

### Insertion

Insert the Filter Insert into the handle of the Debrider until the Filter Insert sits flush at the bottom of the handle.

**NOTE:** The Filter Insert will not sit flush against the bottom of the Debrider handle if the orientation shown below is not followed.



# Device Re-Assembly

---

## Screw on the Cap & Attach Suction

### Screw On Cap

Screw on the Cap and tighten until the flat sides of the Cap and Handle align.

### Connect Suction Tube

Attach a standard suction tube (6-10mm) to the barbed nozzle of the Cap. Ensure that the suction tube is connected to a waste canister and suction source (HIGH/MAX setting is recommended). Check the cutting tip of the metal shaft to confirm that suction is active.



# Troubleshooting

---

## De-clogging

If the cutting tip of the Debrider gets clogged with cancellous bone, follow these steps:

1. Confirm that the suction tube forms an air-tight connection with the barbed nozzle and that the suction source is powered on high.
2. Perform several additional cutting strokes to allow the Debrider to de-clog itself.

If the cutting tip is still clogged, use the Plunger accessory as shown above. Holding the Plunger by the thumb holder, hook the L of the Plunger into the Debrider opening and follow the curve of the Debrider to de-clog and resume debridement.

---





## Ordering Information

---



Part #	Description	Size
BH-330	<b>Avitus® Bone Debrider</b> - Ø8mm, 255mm length. Sterile, disposable surgical device for debridement of bone such as osteomyelitis, bone infection, prosthetic joint infection etc. Closed capture system for specimen collection of contaminated tissue.	8mm (Shaft Diameter) 255mm (Shaft Working Length)



Harvest. Debride.  
Deliver.

# Avitus®

DragonWing™ | Bone Harvester | Bone Debrider | ArchiMIS™



 **ZIMMER BIOMET**  
Moving You Forward.™

This material is intended for health care professionals. Distribution to any other recipient is prohibited. For indications, contraindications, warnings, precautions, potential adverse effects and patient counselling information, see the instructions for use or contact your local representative; visit [www.zimmerbiomet.com](http://www.zimmerbiomet.com) for additional product information.

Zimmer Biomet does not practice medicine. This technique was developed in conjunction with a health care professional. This document is intended for surgeons and is not intended for laypersons. Each surgeon should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient, and this information does not purport to replace the comprehensive training surgeons have received. As with all surgical procedures, the product/technique used in each case will depend on the surgeon's medical judgment as the best treatment for each patient. Results will vary based on health, weight, activity, and other variables. Not all patients are candidates for this product and/or procedure. Caution: Federal (USA) law restricts this device to sale by or on the order of a surgeon. Rx only.

All content herein is protected by copyright, trademarks, and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated, or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

©2025 Zimmer Biomet. All rights reserved.



ST-0808-01-EN-Issue Date YYYY-MM-DD  
VV-12062



**Legal Manufacturer**

Biomet Orthopedics  
56 East Bell Drive  
PO Box 587  
Warsaw, IN 46581  
[www.zimmerbiomet.com](http://www.zimmerbiomet.com)