

2-Year Outcomes of the Treatment of Defects from Bone Marrow Lesions with Subchondroplasty¹

Zimmer Biomet Summary of Published Article

66 Procedures Performed by Steven B. Cohen, MD, of The Rothman Institute, Philadelphia

Bone marrow lesions (BML) are MRI-visible reactions to stress injuries or fractures of the subchondral bone.¹ These defects occur when bone remodeling fails due to ongoing or increased stress and/or reduced healing capacity in the region.² Patients with BML defects typically have substantial chronic, aching pain, primarily on weight-bearing, and may have a worsening prognosis when left untreated.^{3,4}

In patients with BML defects of the knee, current treatment guidelines recommend an initial course of conservative care, which may improve symptoms and delay further treatment.^{5,6} Patients with these defects, however, are likely to progress to requiring total knee replacement (TKR) or osteotomy,^{6,7} with an incidence reported in one study as 9 times more likely than patients without marrow involvement.⁶ Although generally successful, TKR is associated with prolonged recovery and ongoing functional limitations. Less invasive treatment options with more rapid recovery would be desirable, especially for younger or more-active patients.

A promising treatment for subchondral bone defects associated with BML is Subchondroplasty[®] (SCP[®]), developed by Peter Sharkey, MD and Charles Leinberry, MD.⁸ The SCP Procedure fills cancellous bone defects with AccuFill[®], an injectable, flowable, engineered calcium phosphate bone substitute.^{9,10} AccuFill crystallizes and hardens in an endothermic reaction at 37° C to form a nanocrystalline, macroporous scaffold in the subchondral bone, while also promoting cell-mediated remodeling.

Methods

Consecutive patients treated with SCP between May 2008 and May 2012 were contacted at least 2 years following SCP to determine post-operative state. Patients were evaluated for pain (VAS score) and length of delay to TKR. 2-year status was captured for 60 of 66 patients (91%).

Patient Demographics

- BML on MRI determined to be primary driver of patient symptoms
 - Pain isolated to compartment with BML
- Advanced DJD in compartment with BML, Mid-stage DJD in other 2 compartments
- All patients had been scheduled for TKR
- Non-responsive to other conservative or minimally-invasive therapy

	Age	Height	Weight (lbs)	BMI	Alignment	Symptom Length (months)	ICRS Grade in Treated Compartment
Min	33	5'0"	115	20.3	-10° Valgus	2	0
Avg	55.9	5'7"	195	30.1	1.9° Varus	22.4	3.6
Max	76	6'2"	350	53.2	8° Varus	180	4


Summary of Results

- Substantial, consistent reduction in pain
- 70% of patients delayed TKR at least 2 years
- Reduced rehabilitation time when compared to TKR

Summary of Results


1) SCP can help a patient delay a total knee replacement:

70% of patients who were scheduled to receive TKR were able to delay arthroplasty by at least 2 years.



2) SCP provided significant pain relief:

Pain improved by > 4 points.



3) SCP reduced rehabilitation and recovery as compared to TKR

Study postoperative protocol:

- Strict pain management 48-72 hours postop
- Weight bearing as tolerated with crutches as needed
- Standard post-arthroscopy PT regimen
- Return to Work and Activity as tolerated

Discussion

- SCP was offered as a targeted treatment of BML defect, as an alternative to total knee replacement
- Patients demonstrated and maintained a strong response to SCP treatment
- Reduced rehabilitation and recovery time when compared to TKR
- SCP demonstrated efficacy in a population with advanced DJD

References

1. Eriksen EF, Ringe JD. Bone marrow lesions: a universal bone response to injury? *Rheumatol Int.* 2012;32(3):575-584.
2. Roemer FW, Neogi T, Nevitt MC, *et al.* Subchondral bone marrow lesions are highly associated with, and predict subchondral bone attrition longitudinally: the MOST study. *Osteoarthritis Cartilage.* 2010;18(1):47-53.
3. Lo GH, Hunter DJ, Zhang Y, *et al.* Bone marrow lesions in the knee are associated with increased local bone density. *Arthritis Rheum.* 2005;52(9):2814-2821.
4. Felson DT, Chaisson CE, Hill CL, *et al.* The association of bone marrow lesions with pain in knee osteoarthritis. *Ann Intern Med.* 2001;134(7):541-549.
5. Jevsevar DS, Brown GA, Jones DL, *et al.* The American Academy of Orthopaedic Surgeons evidence-based guideline on: treatment of osteoarthritis of the knee, 2nd edition. *J Bone Joint Surg Am.* 2013;95(20):1885-1886.
6. Scher C, Craig J, Nelson F. Bone marrow edema in the knee in osteoarthrosis and association with total knee arthroplasty within a three-year follow-up. *Skeletal Radiol.* 2008;37(7):609-617.
7. Tanamas SK, Wluka AE, Pelletier JP, *et al.* Bone marrow lesions in people with knee osteoarthritis predict progression of disease and joint replacement: a longitudinal study. *Rheumatology (Oxford).* 2010;49(12):2413-2419.
8. Sharkey PF, Cohen SB, Leinberry CF, Parvizi J. Subchondral bone marrow lesions associated with knee osteoarthritis. *Am J Orthop (Belle Mead NJ).* 2012;41(9):413-417.
9. Cohen SB, Sharkey PF. Surgical Treatment of Osteoarthritis Pain Related to Subchondral Bone Defects or Bone Marrow Lesions: Subchondroplasty. *Techniques in Knee Surgery.* 2012;11(4):170-175. Available at: <http://accellorthopedics.com/wp-content/uploads/2013/01/Cohen-2012-Knee-Creations.pdf>.
10. Farr J,II, Cohen SB. Expanding Applications of the Subchondroplasty Procedure for the Treatment of Bone Marrow Lesions Observed on Magnetic Resonance Imaging. *Operative Techniques in Sports Medicine.* ;21(2):138-143. Available at: [http://www.optechsportsmed.com/article/S1060-1872\(13\)00029-4/abstract](http://www.optechsportsmed.com/article/S1060-1872(13)00029-4/abstract).


All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for health care professionals and the Zimmer Biomet sales force. Distribution to any other recipient is prohibited.

AccuPort® Delivery Cannula and AccuMix® Mixing System are Manufactured by: Zimmer Knee Creations, 56 East Bell Drive, P.O.Box. 587, Warsaw, IN 46581 US. AccuFill® BSM is Manufactured by: ETEX Corporation, 55 Messina Drive, Braintree, MA 02184, USA.

For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert or contact your local representative; visit www.zimmerbiomet.com for additional product information.

©2022, 2023 Zimmer Biomet

 **Legal Manufacturer**
Zimmer Knee Creations
56 East Bell Drive
P.O. Box 587
Warsaw, IN 46581
USA

 **Legal Manufacturer**
AccuFill BSM Manufactured By:
Etex Corporation
55 Messina Drive,
Braintree, MA 02184
USA

