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BI-METRIC XR SERIES

126.5° RPP

A close-up view of the proximal femoral component of a hip implant. The head is covered in a porous, textured coating. The neck is polished metal with a locking mechanism. The neck is angled at 126.5 degrees relative to the proximal femoral shaft.

131.5° RPP

A close-up view of the proximal femoral component of a hip implant. The head is covered in a porous, textured coating. The neck is polished metal with a locking mechanism. The neck is angled at 131.5 degrees relative to the proximal femoral shaft.

131.5° SPP

A close-up view of the proximal femoral component of a hip implant. The head is covered in a porous, textured coating. The neck is polished metal with a locking mechanism. The neck is angled at 131.5 degrees relative to the proximal femoral shaft.

131.5° CDH

A close-up view of the proximal femoral component of a hip implant. The head is covered in a porous, textured coating. The neck is polished metal with a locking mechanism. The neck is angled at 131.5 degrees relative to the proximal femoral shaft.

The Bi-Metric XR Hip Series implants offer four distinctive design strategies for restoration of lateral offset.

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
ANATOMIC OFFSETS
ENHANCE JOINT STABILITY
AND HELP RESTORE HIP
BIOMECHANICS BY
PROVIDING THE
OPPORTUNITY TO TIGHTEN
SOFT TISSUE WITHOUT
CREATING LEG LENGTH
DISCREPANCIES

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POROUS COATING VS. OSTEOLYSIS

Author	Reference	Hip Implant System	Years Followed	Osteolysis
Evans and DeLee	Submitted for Publication	Bi-Metric (Biomet)	5-10 years	0.0%
Mauerhan, et al.	J. Arthroplasty, 1997	Integral (Biomet)	5-8 years	0.0%
McLaughlin	JBJS	Taperloc (Biomet)	8-12.5 years	6.0%
Head, et al.	Orthopedics, 1999	Mallory-Head (Biomet)	11 years avg.	0.0%
Rothman	Orthopedics, 1994	Taperloc (Biomet)	7 years	3.0%
Bourne, et al.	Hip Society (March), 2001	Mallory-Head (Biomet)	10-13 years	0.0%
Multi-Center Study	Biomet Clinical Report, 1994	Taperloc (Biomet)	5 years	0.4%
Meding, et al.	AAOS (Fed.), 2001	Bi-Metric, Integral (Biomet)	10-12 years	0.0%
Capello, McClain	Trans. Int'l Sym., 1992	Omnifit (Osteonics)	2-6 years	44.7%
Heekin, et al.	JBJS, 1993	PCA (Howmedica)	5-7 years	18.0%
Woolson, Maloney	J. Arthroplasty, 1992	Harris/Galante (Zimmer)	3.5 yrs. avg. follow-up	22.0%
Kim, et al.	Orthop. Trans, 1992-3	PCA (Howmedica)	2-7 years	37.0%
Kim, et al.	Orthop. Trans, 1992-3	AML (DePuy)	2-7 years	55.8%
Smith, Harris	CORR, 1995	Harris/Galante (Zimmer)	4.5 yrs. avg. follow-up	31.0%
Engl	Presentation, 1992	AML (DePuy)	7.5 yrs. avg. follow-up	28.0%
Kim, et al.	CORR, 1999	AML (DePuy)	11.3 yrs. follow-up	17.0%

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BIOMET INNOVATIONS:

- TITANIUM
- BI-PLANAR TAPER GEOMETRY
- CIRCUMFERENTIAL PLASMA SPRAY POROUS COATING

UNCOMPROMIZED TRADITION AND TECHNOLOGY

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**Primary Results at 10-12 year follow-up/
100% survivorship in 105 hips¹**

- No Osteolysis
- HHS Increased from an Average of 26 Preoperatively to 92 Postoperatively
- No Revisions

8-11 year follow-up/ 100% Aseptic Survivorship in 67 hips²

- 3% Thigh Pain
- 1.5% Osteolysis
- No Revisions for Aseptic Loosening

5.15 year follow-up of 118 hips³

- No Osteolysis
- No Thigh Pain.
- 100% Survivorship

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REFERENCES & ADDITIONAL SUPPORT MATERIAL

- Anthony, P.; *et al.*: "Localized Endosteal Bone Lysis in Relation to the Femoral Components of Cemented Total Hip Arthroplasties." *J. Bone Joint Surg.*, B: 971-979, November, 1990.
- Bourne, R.B.; *et al.*: "Ingrowth Surfaces: Plasma Spray Coating to Titanium Alloy Hip Replacements." *CORR*, 298: 37-46, 1994.
- "Clinical Evaluation of Titanium Alloy Cementless Total Hip Replacement: A1-5 Year Multi-Center Study." Biomet, Inc., Clinical Report, 1994.
- Davey, J.R.; Tozakoglou, E.: "The Role of Lateral Offset Stems." *Orthop. Trans.*, 22(1): 273, 1999.
- DeHeer, *et al.*: "Differential Activation of Macrophages by Implant Wear Debris." *Trans. Implant Retrieval Symposium of the Society for Biomaterials*, 85, 1992.
- Evans, J.: "Outcome of a Tapered, Titanium, Proximal Load-Bearing Non-Cemented Femoral THA Component: A Minimum 5-Year Follow-Up Study." Presented at AAOS, New Orleans, LA, March 19-23, 1998.
- Head, W.C.: "Mallory-Head Porous Press-Fit Primary Hip Replacement." Presented at the Tenth Annual International Symposium: New Developments in Total Joint Reconstruction, Lake Tahoe, Nevada, June 14-16, 1993.
- Head, W.C.; Mallory, T.H.; Emerson Jr., R.H.: "The Proximal Porous Coating Alternative for Primary Total Hip Arthroplasty." *Orthop.*, 22: 813, 1999.
- Horowitz, S.; *et al.*: "Microphage Exposure to Polymethyl Methacrylate Leads to Mediator Release and Injury." *J. Orthop. Research*, 9(3): 406-413, 1991.
- Jiranek, W.: "The Bi-Metric Component at 8-11 Years." Presented at the '98 Harvard Hip Course.
- Keisu, K.; Orozco, F.; Sharkey, P.; Hozack, W.; Rothman, R.: "Primary Cementless Total Hip Arthroplasty in Octogenarians." *J. Bone Joint Surg.*, 83-A: 359, 2001.
- Mallory, T.H.: "Minimum 10-Year Results of a Tapered Cementless Femoral Component in Total Hip Arthroplasty." Presented at Festschrift Celebration, May 2001—to be published in *J. Arthroplasty*.
- Mallory, T.H.; *et al.*: "Clinical and Radiographic Outcome of a Cementless Titanium Plasma-Spray Coated Total Hip Arthroplasty Femoral Component: Justification for Continuance of Use." Presented at the Annual Meeting of the AAOS, Orlando, FL, February, 1995.
- Maloney, W.; *et al.*: "Fibroblastic Response to Particulate Metallic Debris." *Trans. Implant Retrieval Symposium of the Society of Biomaterials*, 34, 1992.
- Mauerhan, D.R.; Mesa, J.; Gregory, A.; Mokris, J.: "Integral Porous Femoral Stem 5 to 8 Year Follow-up Study." *J. of Arthroplasty*, 12(3): 250-255, 1997.
- Meding, J.B.: "Minimum Ten-Year Follow-Up of a Straight-Stemmed, PlasmaSprayed, Titanium-Alloy, Uncemented Femoral Component." Presented at AAOS, San Francisco, CA, February 28 - March 4, 2001.
- Symposium: "Porous Coating Methods: The Pro's and Cons." *Contemporary Orthop.*, 27(3): 269-296, 1993.
- Tanzer, M.; *et al.*: "The Progression of Femoral Cortical Osteolysis in Association with Total Hip Arthroplasty Without Cements." *J. Bone Joint Surg.*, 74-A: March, 1992.

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