
E1 Bi-Polar

BIOMET

E1 RINGLOC BI-POLAR

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References

- Kurtz, S., et al. The UHMWPE Handbook: Ultra High Molecular Weight Polyethylene in Total Joint Replacement (2nd ed.). Elsevier Academic Press. San Diego, CA. 2004, pp 333-5.
- Tower, et al. Rim Cracking of the Cross-Linked Longevity Polyethylene Acetabular Liner After Total Hip Arthroplasty. The Journal of Bone and Joint Surgery. 89:2212-7, 2007.
- Wang, A., et al. Effect of Thermal Treatment of the Wear of Radiation-Crosslinked UHMWPE with and without Vitamin E. 5th UHMWPE International Meeting, 2011.
- Oral, E., et al. The Effects of High Dose Irradiation on the Cross-linking of Vitamin E-Blended Ultrahigh Molecular Weight Polyethylene. Biomaterials. 29(26): 3557-60, 2008.
- Currier, B., et al. Evaluation of Oxidation and Fatigue Damage of Retrieved Crossfire Polyethylene Acetabular Cups. The Journal of Bone and Joint Surgery. 89:2023-2029, 2007.
- Reinitz, et al. Early indications of Oxidative Degradation in Retrieved Annealed UHMWPE Bearings Resembles Gamma Sterilized Materials. ORS. Paper 0307, 2012.
- Currier, B., et al. In Vivo Oxidation in Remelted Highly Cross-Linked Retrievals. The Journal of Bone & Joint Surgery. 92:2409-18, 2010.
- Lombardi, A., et al. Mid-Term Results of a Polyethylene-Free Metal-on-Metal Articulation. Journal of Arthroplasty. 19(7 Suppl. 2): 42-7, 2004.
- McLaughlin, J.R. and Lee, K.R. Cementless Total Hip Replacement Using Second-generation Components. Journal of Bone and Joint Surgery (Br). 92(12):1636-41, 2010.
- Nabar, S., et al. Comparison of Second Generation Highly Crosslinked Polyethylenes Under Adverse Aging Conditions. ORS. 2008. Poster No. 1684.
- Freedman, J. & Schroeder, D. Oxidative Stability of Vitamin E Blended and Subsequently Crosslinked UHMWPE. ORS. Paper No. 0236, 2012
- Rowell, S., et al. Comparative Oxidative Stability of -Tocopherol Blended and Diffused UHMWPEs at 3 Years of Real-Time Aging. Journal of Orthopedic Research. 29:773-780, 2011. (This study was funded by Zimmer, Inc., Biomet, Inc. and DePuy, Inc.)
- Audrey K. et al. Five year experience of Vitamin E Stabilized Irradiated Ultrahigh Molecular Weight Polyethylene wear and stability of Regenerex Acetabular Shells and Femoral Components using Radiostereometric Analysis 2014 AAOS Annual Meeting.
- Oral, E., and Muratoglu, O. Vitamin E Diffused, Highly Crosslinked UHMWPE: A Review. International Orthopaedics (SICOT). Online publication. December 2010.
- Kosnir, et al. Lipid Doping and Aging of Various UHMWPEs. ORS. 2012 Paper No. 0311.

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