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
# References

---

1. Amis, A. A.: The strength of artificial ligament anchorages. A comparative experimental study. *J Bone Joint Surg Br*, 70(3): 397–403, 1988.
2. Bailey, S. B.; Grover, D. M.; Howell, S. M.; and Hull, M. L.: Foam-reinforced elderly human tibia approximates young human tibia better than porcine tibia: A study of the structural properties of three soft-tissue fixation devices. *Am J Sports Med*, 2, In Press.
3. Chen, C. H.; Chen, W. J.; Shih, C. H.; Yang, C. Y.; Liu, S. J.; and Lin, P. Y.: Enveloping the tendon graft with periosteum to enhance tendon-bone healing in a bone tunnel: A biomechanical and histologic study in rabbits. *Arthroscopy*, 19(3): 290–6, 2003.
4. Clatworthy, M.; Buelow, J. U.; Pinczewski, L.; Howell, S. M.; Fowler, P.; and Amendola, A.: Tunnel widening in hamstring ACL reconstruction: A prospective clinical and radiographic evaluation of four different techniques. In *ISAKOS*. Edited, Montreux, Switzerland, 1999.
5. Clatworthy, M. G.; Annear, P.; Bulow, J. U.; and Bartlett, R. J.: Tunnel widening in anterior cruciate ligament reconstruction: a prospective evaluation of hamstring and patella tendon grafts. *Knee Surg Sports Traumatol Arthrosc*, 7(3): 138–45, 1999.
6. Corry, I. S.; Webb, J. M.; Clingeffer, A. J.; and Pinczewski, L. A.: Arthroscopic reconstruction of the anterior cruciate ligament: A comparison of patellar tendon autograft and four-strand hamstring tendon autograft. *Am J Sports Med*, 27: 444–454, 1999.
7. Feller, J. A.; Webster, K. E.; and Gavin, B.: Early post-operative morbidity following anterior cruciate ligament reconstruction: patellar tendon versus hamstring graft. *Knee Surg Sports Traumatol Arthrosc*, 9(5): 260–6, 2001.
8. Giurea, M.; Zorilla, P.; Amis, A. A.; and Aichroth, P.: Comparative pull-out and cyclic-loading strength tests of anchorage of hamstring tendon grafts in anterior cruciate ligament reconstruction. *Am J Sports Med*, 27(5): 621–5, 1999.
9. Grassman, S. R.; McDonald, D. B.; Thornton, G. M.; Shrive, N. G.; and Frank, C. B.: Early healing processes of free tendon grafts within bone tunnels is bone-specific: a morphological study in a rabbit model. *Knee*, 9(1): 21–6, 2002.
10. Greis, P. E.; Burks, R. T.; Bachus, K.; and Luker, M. G.: The influence of tendon length and fit on the strength of a tendon-bone tunnel complex. A biomechanical and histologic study in the dog. *Am J Sports Med*, 29(4): 493–7, 2001.
11. Howell, S. M.: Compaction of an autogenous dowel of cancellous bone minimizes tunnel widening at 4 months. *Am J Sports Med*, Submitted.
12. Kyung, H. S.; Kim, S. Y.; Oh, C. W.; and Kim, S. J.: Tendon-to-bone tunnel healing in a rabbit model: the effect of periosteum augmentation at the tendon-to-bone interface. *Knee Surg Sports Traumatol Arthrosc*, 11(1): 9–15, 2003.
13. Magen, H. E.; Howell, S. M.; and Hull, M. L.: Structural properties of six tibial fixation methods for anterior cruciate ligament soft tissue grafts. *Am J Sports Med*, 27(1): 35–43, 1999.
14. Rodeo, S. A.; Suzuki, K.; Deng, X. H.; Wozney, J.; and Warren, R. F.: Use of recombinant human bone morphogenetic protein-2 in enhance tendon healing in a bone tunnel. *American Journal of Sports Medicine*, 27(4): 476–488, 1999.
15. Roos, P.; Howell, S. M.; and Hull, M. L.: Compaction of an autogenous dowel of cancellous bone improves strength and stiffness of tibial fixation. *Am J Sports Med*, Submitted.
16. Singhatat, W.; Lawhorn, K. W.; Howell, S. M.; and Hull, M. L.: How four weeks of implantation affect the strength and stiffness of a tendon graft in a bone tunnel: a study of two fixation devices in an extraarticular model in ovine. *Am J Sports Med*, 30(4): 506–13, 2002.
17. To, J. T.; Howell, S. M.; and Hull, M. L.: Contributions of femoral fixation methods to the stiffness of anterior cruciate ligament replacements at implantation. *Arthroscopy*, 15(4): 379–87, 1999.
18. Tomita, F.; Yasuda, K.; Mikami, S.; Sakai, T.; Yamazaki, S.; and Tohyama, H.: Comparisons of intraosseous graft healing between the doubled flexor tendon graft and the bone-patellar tendon-bone graft in anterior cruciate ligament reconstruction. *Arthroscopy*, 17(5): 461–76, 2001.
19. Wallace, M. P.; Howell, S. M.; and Hull, M. L.: In vivo tensile behavior of a four-bundle hamstring graft as a replacement for the anterior cruciate ligament. *J Orthop Res*, 15(4): 539–45, 1997.
20. Weiler, A.; Hoffmann, R. F.; Bail, H. J.; Rehm, O.; and Sudkamp, N. P.: Tendon healing in a bone tunnel. Part II: Histologic analysis after biodegradable interference fit fixation in a model of anterior cruciate ligament reconstruction in sheep. *Arthroscopy*, 18(2): 124–35, 2002.

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