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
References

1. Giron, Francesco: Femoral Attachment of the Anterior Cruciate Ligament. *Knee Surg Sports Traumatol Arthrosc* (2006) 14: 250-256.
2. Zavras, T. D.; Race, A.; and Amis, A. A.: The Effect of Femoral Attachment Location on Anterior Cruciate Ligament Reconstruction: Graft Tension Patterns and Restoration of Normal Anterior-Posterior Laxity Patterns. *Knee Surg Sports Traumatol Arthrosc*, 13(2): 92-100, 2005.
3. Simmons, R.; Howell, S. M.; and Hull, M. L.: Effect of the Angle of the Femoral and Tibial Tunnels in the Coronal Plane and Incremental Excision of the Posterior Cruciate Ligament on Tension of an Anterior Cruciate Ligament Graft: An *in vitro* Study. *J Bone Joint Surg Am*, 85-A(6): 1018-29, 2003.
4. Sasaki, Norihiro; Ishibashi, Yasuyuki ; Tsuda, Eiichi: Femoral Insertion of the Anterior Cruciate Ligament: Discrepancy Between Macroscopic and Histological Observations. *Arthroscopy: The Journal of Arthroscopic and Related Surgery*, Vol 28, No 8 (August), 2012: pp 1135-1146.
5. Howell, S.M.; Gittins, M.E.; Gottlieb, J.E.; Traina, S.M.; and Zoellner, T.M.: The Relationship between the Angle of the Tibial Tunnel in the Coronal Plane and Loss of Flexion and Anterior Laxity after Anterior Cruciate Ligament Reconstruction. *American Journal of Sports Medicine*, 29(5); 2001, pp 567-74.
6. Markolf, K. L.; Park, S.; Jackson, S. R.; and McAllister, D. R.: Anterior-Posterior and Rotatory Stability of Single and Double-Bundle Anterior Cruciate Ligament Reconstructions. *J Bone Joint Surg Am*, 91(1): 107-18, 2009.
7. Markolf, K. L.; Park, S.; Jackson, S. R.; and McAllister, D. R.: Contributions of the Posterolateral Bundle of the Anterior Cruciate Ligament to Anterior-Posterior Knee Laxity and Ligament Forces. *Arthroscopy. The Journal of Arthroscopic and Related Surgery*, VOL 24, No 7 (July), 2008: pp 805-809.
8. Markolf, K. L.; Park, S.; Jackson, S. R.; and McAllister, D. R.: A Comparison of 11 O'clock Versus Oblique Femoral Tunnels in the Anterior Cruciate Ligament-Reconstructed Knee. *The American Journal of Sports Medicine*, Vol. 38, No. 5.
9. Rahr-Wagner, L. R.: Increased Risk of Revision After Anteromedial Compared With Transtibial Drilling of the Femoral Tunnel During Primary Anterior Cruciate Ligament Reconstruction. Results from the Danish Knee Ligament Reconstruction Register. *Arthroscopy: The Journal of Arthroscopic and Related Surgery*, 29(1); (January), 2013: pp 98-105.
10. Morgan, Joseph A.: Femoral Tunnel Malposition in ACL Revision Reconstruction. *J Knee Surg*. 2012 November; 25(5): 361-368.

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