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INDICATIONS: The Nextremity Solutions Nextra® Hammertoe Correction Implant is indicated for small bone reconstruction limited to inter-phalangeal repair and fusion of the lesser toes.
CONTRAINDICATIONS: The Nextremity Solutions Nextra® Hammertoe Correction Implant is NOT intended for use in procedures involving the great toe. In addition, the device is contraindicated in the following: (1) Pathological conditions of bone, such as cystic changes or severe osteopenia, osteoporosis, bony deficiencies, or comminuted bone surface which would adversely affect the Nextra® Hammertoe Correction Implant; (2) Physical conditions that would eliminate, or tend to eliminate, adequate implant support or retard healing, including inadequate soft tissue coverage; (3) Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period; (4) Sensitivity or allergy to the metal implant. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation; (5) Blood supply limitations and previous infections that may retard healing; (6) Presence of an active infection; (7) Pediatric patients with open epiphyseal plates.

This material is intended for health care professionals. Distribution to any other recipient is prohibited. For product information, including indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert. Check for country product clearances and reference product-specific instructions for use.

Zimmer Biomet and Nextremity Solutions do not practice medicine. This technique was developed in conjunction with health care professionals. This document is intended for surgeons and is not intended for laypersons. Each surgeon should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient, and this information does not purport to replace the comprehensive training surgeons have received. As with all surgical procedures, the technique used in each case will depend on the surgeon's medical judgment as the best treatment for each patient. Results will vary based on health, weight, activity and other variables. Not all patients are candidates for this product and/or procedure. Caution: Federal (USA) law restricts this device to sale by or on the order of a surgeon. Rx only.

Nextremity Solutions, Nextra, and RevLock are trademarks of Nextremity Solutions, Inc. Zimmer Biomet is the exclusive distributor of the Nextra® Hammertoe Correction System.

The Nextra Hammertoe Correction system is manufactured using 316L Stainless Steel.

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©2018 Nextremity Solutions LIT-1480 Rev.08
Zimmer, Inc.
Warsaw, IN 46580 U.S.A.