



THE **ABCs** OF LOP



Have any of your patients complained of post-operative tourniquet pain or transient nerve pain? If so, are you using a standard 300 mm Hg for thigh pressures & 250 mm Hg for upper arm?

Make no mistake; Limb Occlusion Pressure (LOP) is not just any pressure that stops blood flow to the limb as the competition may want to portray.

We know that patients all have different systolic blood pressure, limb shape & circumference – so why would we want to use unnecessarily high tourniquet pressure on all of your patients?

The Zimmer Biomet A.T.S.[®] 4000 measures the minimum cuff pressure to stop arterial blood flow and reports this back to the healthcare professional. Concurrently, it also displays the Recommended Tourniquet Pressure (RTP) for that specific patient, which accounts for any fluctuations in blood pressure throughout the case.

Let's discuss the ABCs of LOP

The background features three large, stylized letters: a dark blue 'A' on the left, a dark blue 'B' on the right, and a dark blue 'C' at the bottom. The letters are semi-transparent, allowing the lighter blue background to show through. The text is overlaid on these letters.

A IS FOR “ABSOLUTE MINIMUM”

B IS FOR “BLOODLESS FIELD”

C IS FOR “CONFIDENCE”

Absolute Minimum

The true definition of LOP as defined by McEwen *et al*, is the minimum cuff pressure that stops arterial blood flow distal to the cuff.⁶ This definition is recently well supported by Tuncali *et al*, as being called arterial occlusion pressure (AOP).⁴ The term AOP and LOP are synonyms. Using the minimum effective pressure (i.e. LOP) added to a pressure margin has been shown to allow the use of lower tourniquet cuff pressures when compared to recommended standards without affecting limb hemostasis.²

Bloodless Field

It has been well reported and studied that high tourniquet pressure can cause injury to the tissue under and near the tourniquet cuff, nerve damage, and could cause weakening, partial or complete muscle paralysis.^{1,2,7} Tourniquet cuff pressures that are too low could result in bleeding distal to the cuff which can result in undesirable blood in the surgical field and blood loss. Insufficient tourniquet cuff pressures have been reported to impair venous return while allowing arterial blood flow to continue.³ Hence, setting tourniquet pressures too low can cause pooling of blood in the field.



Confidence

Researchers have shown that using a Doppler stethoscope (the gold standard) in the determination of the minimum effective tourniquet pressure to be quite useful and effective.⁵ A similar study concluded that there was no significant difference between the LOP obtained by the plethysmograph-equipped tourniquet system [A.T.S.® 3000] and the LOP obtained with a Doppler stethoscope.⁶ Moreover, the study found that using a standard width tourniquet cuff and using the plethysmograph technique to sense LOP and adding a pressure margin, the patient experiences tourniquet cuff pressure of 12 to 25% less when compared to current practice and tourniquet hemostasis was successfully obtained in 93% of the cases.⁶ The A.T.S. 4000, similar to the A.T.S. 3000, does this entire process in approximately 30 seconds. Reid *et al*, showed that the LOP plus a pressure margin gave good hemostasis in 93.2% of the upper extremity procedures and 97.5% in lower extremity procedures.² The A.T.S. 4000 reports the LOP and RTP and gives the healthcare professional an option in setting the tourniquet cuff pressures.

With the minimal time required, ease of use and well documented benefits of lower tourniquet pressures, the A.T.S. 4000 Automatic Tourniquet System offers the healthcare professional an option in setting the tourniquet cuff pressures while always keeping the final setting of pressure in the healthcare professional's hands. Minimizing tourniquet pressure and pressure gradients with LOP helps minimize the risk of post-op tourniquet pain and nerve-related tourniquet injuries.⁸

References

1. Tredwell S.J., Wilmink M., Inkpen K., McEwen J.A., "Pediatric tourniquets: Analysis of cuff and limb interface, current practice, and guidelines for use.", *J Pediatr Orthop* 21:671-676, 2001.
2. Reid H.S., Camp R.A., and Jacob W.H., "Tourniquet hemostasis: a clinical study", *Clin Orthop* 177:230-234, 1983
3. Massey K.A., Blakeslee C., Martin W, Pitkow H.S., "Pneumatic ankle tourniquets: physiological factors related to minimal arterial occlusion pressure", *J Foot Ankle Surg* 38:256-263, 306-307, 1999
4. Tuncali B, Karci A, Tuncali BE, Mavioglu O, Ozkan M, Bacakoglu AK, Baydur H, Ekin A, Elar Z. "A new method for estimating arterial occlusion pressure in optimizing pneumatic tourniquet inflation pressure.", *Anesth Analg*. 2006 Jun; 102(6):1752-7
5. Diamond E.L., Sherman M., Lenet M., "A quantitative method of determining the pneumatic ankle tourniquet setting", *J Foot Surg* 24:330-334, 1985
6. McEwen J.A., Inkpen K.B., Younger A., "Thigh tourniquet safety: Limb occlusion pressure measurement and a wide contoured cuff allow lower cuff pressure" *Surg Tech* 34:8-18, 2002
7. McEwen J.A., Kelly D.L., Jardanowski T., Inkpen K., "Tourniquet safety in lower leg applications", *Orthop Nurs* 21:55-62, 2002
8. Noordin S, McEwen JA, Kragh Jr, Eisen A, Masri BA. "Surgical Tourniquets in Orthopaedics" *The Journal of Bone & Joint Surgery* 2009 Dec; 91(12):2958-67

This publication and all content, is protected by copyright, trademarks, and other intellectual property rights owned by, or licensed to, Zimmer Biomet or its affiliates unless otherwise indicated. This material is intended for physicians and the Zimmer Biomet sales force only. Distribution to any other recipient is prohibited. This publication must not be used, copied, or reproduced in whole or in part without the express written consent of Zimmer Biomet or its authorized representatives. For product information see the package insert and zimmerbiomet.com. Zimmer Biomet does not practice medicine, does not recommend a particular orthopaedic implant or surgical technique, and is not responsible for the kind of treatment selected for a specific patient. The surgeon who performs any implant procedure is responsible for determining and utilizing the appropriate techniques for implanting a prosthesis in each individual patient. ©2017 Zimmer Biomet



ZIMMER BIOMET

Your progress. Our promise.™