



DRAGONWING™

LARGE VOLUME **AUTOGRAFT** DELIVERY

STREAMLINE LARGE VOLUME AUTOLOGOUS BONE GRAFT DELIVERY

Rapidly advance reloadable graft cartridge

A SMOOTH MINIMALLY INVASIVE EXPERIENCE

Ø5.2mm Delivery Cannula | Deliver Autograft without Jamming



DRAGONWING™


LARGE VOLUME AUTOGRAFT DELIVERY

Autograft is often used for:

Large Cysts, Lesions, AVN¹⁻³



A large volume, minimally invasive autograft delivery experience




Arthroscopic and Minimally Invasive Fusions and Nonunions^{4,5}

DRAGONWING™

LARGE VOLUME AUTOGRAFT DELIVERY

CLINICAL APPLICATION EXAMPLES

FUSIONS & NONUNIONS^{4,5}
ARTHROSCOPIC/MIS GRAFTING⁴
CYSTS, BONE VOIDS, AVN, LESIONS¹⁻³




The device is intended to be used by a physician familiar with the possible side effects, typical findings, limitations, indications, and contraindications of bone graft delivery

DRAGONWING™

LARGE VOLUME AUTOGRAFT DELIVERY

SYSTEM CONTENTS

1 x Delivery Plunger



1 x Pack Tool

1 x 1 cc Bone Scoop

2 x 5cc Reloadable Graft Cartridges

ø5.2mm x 120mm
Straight Delivery Cannula

REF CODE	DESCRIPTION	RECOMMENDED FOR:
AD-5ch-5.2-120	Avitus® DragonWing™ Large Volume Auograft Delivery. ø5.2mm Delivery Cannula, 120mm Length	-Filling 3cc-20cc Defects -Delivering morcelized autogenous bone graft

Refer to device specific instructions for use for information on indications for use, contraindications, warnings, and precautions. APM042.B 2023.12 (DCO-0499)

REFERENCES

- 1 Singh, H., Agarwal, K. K., Tyagi, S., Rampurwala, A., Singh, A. K., Bhrambhath, P., Unjia, R., Agarwal, N., & Rampurwala, A. J. (2023). An evaluation of core decompression and cancellous bone grafting for Early-Stage avascular necrosis of the femoral head. *Cureus*. <http://doi.org/10.7759/cureus.37878>
- 2 Alisi, M. S., Hammad, Y., Khanfar, A., Samarah, O. Q., & Hassan, F. A. (2022). Percutaneous curettage and local autologous cancellous bone graft: a simple and efficient method of treatment for benign bone cysts. *PubMed*, 10(1), 104-111. <https://doi.org/10.22038/abjs.2021.55189.2747>
- 3 Lui, T. H. (2017). Arthroscopic curettage and bone grafting of bone cysts of the talar body. *Arthroscopy Techniques*, 6(1), e7-e13. <https://doi.org/10.1016/j.eats.2016.08.029>
- 4 Iturregui, J. M., Moses, A. M., Shi, G. G., & Haupt, E. T. (2023). Contemporary Review: Autograft Bone use in foot and ankle surgery. *Foot & Ankle Orthopaedics*, 8(1), 247301142311531. <https://doi.org/10.1177/24730114231153153>.
- 5 Roberts, T. T., & Rosenbaum, A. J. (2012). Bone grafts, bone substitutes and orthobiologics. *Organogenesis*, 8(4), 114-124. <https://doi.org/10.4161/org.23306>

Intended for Healthcare Professionals. All content herein is protected by copyright, trademarks, and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated, or disclosed, in whole or in part, without the express written consent of Zimmer Biomet. ©2023 Zimmer Biomet.

