

Convenient Closing for Shoulder Arthroplasty

SUBSCAPULARIS REPAIR KITS

Streamline your subscapularis repair process with the convenience of **Subscapularis Repair Kits**.


The sterilized kits contain a comprehensive suture, tape and drill bit arrangement to facilitate soft tissue repair.

Flexible configurations equip you with everything needed for this procedure in one concise kit.



BroadBand[®] Tape

- Provides a broader footprint over tissue¹
- Easy-to-handle material: Smooth, silk-like feel
- Coreless, single weave braid allows for low profile knots



MaxBraid[®] Sutures

- Proprietary braid allows suture to lie flat when tied
- Made from 100% polyethylene
- Diverse selection of suture sizes, colors and needle options

Convenient

- Complete, sterilized kit
- Intuitive packaging provides accessible suture identification and application throughout procedure
- Minimal materials to manage

Gentle

- BroadBand® Tape provides a larger footprint over tissue¹
- Non-abrasive, gentle on tissue and gloves

Versatile

- Assorted needle types to help facilitate suture/tape passage
- Multiple kit configurations for maximizing repair effort
- Distinct colors provide optimal visibility for suture management and repair


Ordering Information

Basic Kit for Subscapularis Repair

| Part Number | Description | Qty |
|-------------|---|-----|
| 110046602 | Tagging Suture Pack – MaxBraid USP#2 Suture w/ MO-4 Needle | 2 |
| | Shuttle Suture Pack – 5 in. MaxBraid USP#2 Suture Loop, USP#2 Suture Tail w/ CP-1 Needle and CCS Needle | 2 |
| | Repair Suture Pack – 1.5 mm BroadBand Suture Tape w/ MO-4 Needle | 4 |
| | Closure Pack – 1.5 mm BroadBand Suture Tape w/ MO-4 Needle | 1 |
| | Drill Bit – 2.0 x 100 mm | 1 |

Ultra Kit for Subscapularis Repair

| Part Number | Description | Qty |
|-------------|--|-----|
| 110046455 | Tagging Suture Pack – MaxBraid USP#2 Suture w/ MO-4 Needle | 2 |
| | Shuttle Suture Pack – 5 in. MaxBraid USP#2 Suture Loop, 3 in. USP#2 Suture Tail w/ Needles | 3 |
| | 36 mm Reverse Cutting Needle | – |
| | 40 mm Taper Needle | – |
| | 48 mm Conventional Cutting Needle | – |
| | Repair Suture Pack – 1.5 mm BroadBand Dual-Strand Suture Tape w/ CT Needle | 3 |
| | Closure Pack – 1.5 mm BroadBand Suture Tape w/ MO-4 Needle | 2 |
| | Drill Bit – 2.0 x 100 mm | 1 |



material is intended for health care professionals. Distribution to any other recipient is prohibited. For indications, contraindications, warnings, precautions, potential adverse effects and patient counselling information, see the package insert or contact your local representative; visit www.zimmerbiomet.com for additional product information. ©2022 Zimmer Biomet