



ACL RECONSTRUCTION

# ToggleLoc<sup>®</sup> Flip Inline with ZipLoop<sup>®</sup> Technology

**TACTILE FLIPPING.  
CONTROLLED TENSIONING.**

ToggleLoc Flip Inline with ZipLoop Technology features an adjustable button loop that allows for one size fits all surgeon-controlled graft tensioning for ACL reconstruction. The strong titanium suspensory button fixation device is designed with a flipping suture for reproducible button deployment.

## DYNAMIC

- Self-adjusting loop covers 15 to 60 mm ACL loop repair lengths.
- Additional suture holes allow surgeon-controlled suture flipping and tactile deployment of the button.
- Designed with a single, long protective card to reduce coiling of the loop and suture.

## STRONG

- Strong titanium suspensory button fixation device with a high suture pullout tensile strength of 900 Newtons.<sup>1</sup>
- Distinct 3-point locking ZipLoop design and braid construction allow for easy placement and are designed to help minimize suture creep.

## FLEXIBLE

- Saddle design provides smooth control for ACL graft tensioning.
- Adjustable loop can tension and close as low as 5 mm to the end of the button, maximizing graft biology in bone.
- Plate button style fits flush in extender or washer if ACL socket is fully reamed out.



### PRODUCT NAME

ToggleLoc Flip Inline with ZipLoop Technology

### PART NUMBER

OBAL80

References:

1. Riverpoint Internal Test Report [Product Bench Testing Button Loop, 2022] VP-02N-TR on file.

ToggleLoc Flip Inline with ZipLoop Technology is an OrthoButton AL product.

OrthoButton AL is a trademark of Riverpoint Medical.

This material is intended for healthcare professionals. For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert or contact your local representative; visit [www.zimmerbiomet.com](http://www.zimmerbiomet.com) for additional product information. All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

©2023 Zimmer Biomet



**Legal Manufacturer**  
 Riverpoint Medical  
 825 NE 25th Ave.  
 Portland, OR 97232

[www.zimmerbiomet.com](http://www.zimmerbiomet.com)