

XPERIENCE™

Advanced Surgical Irrigation

No-Rinse Wound Lavage





What is XPERIENCE?

XPERIENCE™ is a **no-rinse** wound lavage.

It is designed to rinse away debris and microorganisms that may cause post-operative complications.

Features

- 5 hours of ongoing protection against bacteria¹
- No secondary rinse out required, which may reduce surgery time
- Safety tested showing no irritation²
- Broad spectrum efficacy within the solution¹
- Compatible with most commonly used implants and closure methods*,³
- 99.9% bacterial elimination⁴

XPERIENCE[®]
Advanced Surgical Irrigation








Warning: NOT FOR IV USE

**Warm to Body Temperature
Prior to Use**
(Do not exceed 40°C for more than 6 months)

REF SL-0005 **QTY** 1 (500 mL)

CONTAINS Sodium Citrate,
Citric Acid, Sodium Lauryl
Sulfate, and Purified Water.

15°C 25°C  

  **STERILE** 


Rx Only

Manufactured for:
Next Science, LLC
10550 Deerwood Park Blvd
Suite 300
Jacksonville, FL 32256 USA

*Does not include hydroxyapatite (HA) or fibrin sealants

Efficacy and Safety


Bacteria Reduction Over Time



Laboratory results are not necessarily indicative of clinical performance.

Proven Safety

XPERIENCE has been safety tested showing no irritation.²



Ordering Information

Description	Item Number
XPERIENCE Wound Lavage 4pk	0088750400
XPERIENCE Wound Lavage 10pk	0088751000

For more information, please contact your Zimmer Biomet sales representative or visit zimmerbiomet.com/XPERIENCE

References

1. Bashyal RK, Mathew M, et. al. A Novel Irrigant to Eliminate Planktonic Bacteria and Eradicate Biofilm Superstructure with Persistent Effect During Total Hip Arthroplasty. *Journal of Arthroplasty*. 2022.
2. Data on file. Next Science Testing Report CD-0040
3. Data on file. Next Science Testing Report TR-06-19-012
4. Data on file. Next Science Testing Report TR-06-19-007

This material is intended for health care professionals. Distribution to any other recipient is prohibited.

For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert or contact your local representative; visit www.zimmerbiomet.com for additional product information.

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

Check for country-specific product clearances. Not for distribution in France.

©2022 Zimmer Biomet



3875.1-US-en-Issue Date 2022-06



Legal Manufacturer
Next Science, LLC
10550 Deerwood Park Blvd
Suite 300
Jacksonville, FL 32256 USA

zimmerbiomet.com
1-800-253-6190