

Acetabular Bone Screw Combinations

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	Bone Screw, Self-tapping, 4.5mm 00-6250-045-xx	Bone Screw, Self-Tapping, 6.5mm 00-6250-065-15/70	Bone Screw, Self-Tapping, 6.5mm 00-6250-065-xx	Cancellous Bone Screw 02.03147,020/110 P-100	Cancellous Screw 02.06000,015/070	Cancellous Screw Universal 14.00.015/060	Countersunk Cancellous Bone Screw 01.00118,015/070 P-64	Countersunk Cancellous Bone Screw 42.19.1570 P-100	Countersunk Cancellous Bone Screw 4301-07-015/080	HGP II, 4.5mm Bone Screws 00-6624-045-xx	HGP II, 6.5mm Bone Screws 00-6624-065-xx	Novafit Cancellous Screw 41.14.27/40	Peripheral Bone Screw 5.0 00-7250-050-20/65	Self-Tapping Cancellous Bone Screw 6301-07-015/080	Stabilizing Nails 1194/1197	Biomet Low Profile Self-Tapping Bone Screws 1035-30/39-00
Acetabular Cups (Shells):																
ACA cup									✓							
Continuum Acetabular Cup, Elliptical Porous Coated													✓			
ALLOFIT-S IT Alloclassic Shell (8755)		✓	✓						✓							
ALLOFIT-S Alloclassic Shell		✓	✓		✓				✓					✓		
APR Acetabular Shell									✓					✓		
Cedior Press-fit Cup						✓										
Continuum Acetabular Cup (8757)	✓	✓	✓													✓
Continuum Hedrocel Acetabular Revision Shell											✓					
Converge Porous Shell									✓					✓		
FITMORE shell					✓		✓	✓								
G7 PPS Limited Hole Acetabular Shell		✓														
G7 PPS Bonemaster Limited Hole Acetabular Shell		✓														
G7 PPS Finned Limited Hole Acetabular Shell		✓														
G7 PPS Bonemaster Finned Limited Hole Acetabular Shell		✓														
G7 Osseo Ti Limited Hole Acetabular Shell		✓														
G7 Osseo Ti Multi Hole Acetabular Shell		✓														
G7 Bispherical Limited Hole Acetabular Shell		✓														
G7 Bispherical HA Limited Hole Acetabular Shell		✓														
G7 Bispherical Non-Plugged Limited Hole Acetabular Shell		✓														
HGP shell										✓	✓					
LOGOS shell						✓										
LOR shell					✓		✓	✓	✓					✓		
Marburg METASUL cup															✓	
Marburg PE Cup															✓	
Novafit PE Cup												✓				
SL Titanium shell							✓	✓								
St. Nabor shell				✓			✓	✓								
STANDARD cup					✓		✓	✓								
TM Modular Shell (00-6202-XXX-XX)	✓	✓	✓													✓
TM Monoblock Cup													✓			
TM Revision Shell (00-7000-XXX-20)											✓					
TM Revision Shell, Flat Ring (00-7000-XXX-70)											✓					
Trilogy AB Shell	✓	✓	✓													
Trilogy IT Acetabular Cup (8753)	✓	✓	✓													✓
Trilogy Shell	✓	✓	✓													✓

Acetabular Reconstruction Cages/Rings

Burch-Schneider Cage				✓			✓	✓	✓						✓	
Zimmer Reinforcement Ring with Hook*				✓			✓	✓	✓						✓	
Original M.E. Müller Roof Reinforcement Ring				✓			✓	✓							✓	
TM Revision Cage		✓	✓				✓	✓			✓					
Wagner Reinforcement Cage							✓	✓								
ZCA Cage	✓	✓	✓													
ZCA Roof Ring	✓	✓	✓													

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.