

Reference Combination Manual/Automated Cleaning and Disinfection Instructions (Section H) of the Instructions for Care, Cleaning, Maintenance and Sterilization Manual (97-5000-170-00). Do not utilize the optional ultrasonic cleaning note in Section H.



## Persona® Shims

42-5279-001-00

42-5279-003-00


42-5279-005-00

42-5279-007-00

47-5279-009-00

### Cleaning Instructions

Do not ultrasonic clean



97-5026-065-00 Rev. 2 Printed in USA ©2015 Zimmer, Inc.