

Articulating Combinations: Metal on Metal

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

Metal femoral heads:

Continuum Bearing Head, Hi Carbon
DUROM Femoral component *
 METASUL eccentric head, 28 & 32mm
 Metasul LDH (Large Diameter Head)
 METASUL standard head, 28 & 32mm
 METASUL standard head, 36 & 40mm

Metal Inserts (Liners):

ALLOCLASSIC CSF METASUL insert
 ALPHA METASUL insert
 APR METASUL insert
 BETA METASUL insert
 CLS SPOTORNO METASUL insert
 CLW METASUL insert
 Epsilon METASUL insert
 GAMMA METASUL insert
 Kappa METASUL insert
 LAMELLA METASUL insert
 Metasul Taper Liner (8770)
 St. Nabor METASUL insert
 STANDARD METASUL insert

		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
				✓	✓
		✓		✓	
		✓		✓	

Metal-lined Articulating Cups:

DUROM Acetabular Component *
 Low Profile METASUL cup
 Marburg METASUL cup
 Morscher Press-fit METASUL cup
 Weber METASUL cup
 Zimmer MMC Cup *

	✓		✓		
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
	✓		✓		

*** NOTE: The Durom Femoral Component is not approved for use with the Durom Acetabular cup and Zimmer MMC Cup for total hip surface replacement Arthroplasty in the U.S.A.**