

**Head and Stem Combinations:
Stainless Steel Titanium Femoral Heads**

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	CPT stainless steel head (6 degree)	PROTASUL-S30 head (12/14 taper)	PROTASUL-S30 head (8/10 taper)	Testa Jupe head (12/14 taper)	Tete a Jupe Acier Inox Stainless Steel 316L Cone 12/14 D28	Tete a Jupe Inox HTA Cone 8/10 D28	TRIBOSUL head (12/14 taper)	TRIBOSUL head (14/16 taper)	Zimtron stainless steel femoral head (12/14 taper)	Zimtron stainless steel femoral head (6 degree)
ACA stem	✓									
Acora cemented stem		✓				✓				
ALLOCLASSIC SLA and SL-64 stem†	✓						✓			
ALLOCLASSIC SLL-Revision stem	✓						✓			
ALLOCLASSIC stem	✓						✓			
ALLOCLASSIC VARIALL SLV stem	✓						✓			
ALLOCOR AN stem	✓						✓			
Anatomic stem										
Anatomique uncemented HA stem	✓			✓			✓			
Apollo Hip System	✓									
APR Hip System	✓									
APS Natural-Hip System										
Avenir Müller stem	✓						✓			
Avenir Cemented stem	✓								✓	
Avenir Complete stem	✓						✓			
BEO stem	✓						✓			
CDH stem PROTASUL-S30 (monoblock)										
CLS Brevius Stem with Kinectiv Technology	Refer to Kinectiv Technology Page									
CLS SPOTORNO stem	✓						✓			
Concept cemented stem		✓					✓			
Continuum Cobrex stem										
Continuum F-115 Hip stem										
Continuum F-117 Collarless Hip stem										
Continuum J-Taper stem										
Continuum Proxilock stem (porous and HA coated)										
CPS uncemented stem	✓		✓				✓			
CPT (6 Deg, CoCr) stem										
CPT (6 Deg, SST) stem	✓									✓
CPT 12/14 CoCr stem										
DSP / TPP / PDC	✓						✓			
Emeraude stem	✓						✓			
Epoch stem										
ESKA 8/10 stem										
ESKA 12/14 stem										
Etoile cemented stem		✓					✓			
EXAFIT Primary stem		✓					✓			
EXAFIT Revision stem		✓					✓			

†ALLOCLASSIC SLA andALLOCLASSIC SL-64 stem cannot be used with VerSys 12/14 femoral heads 32,36, and 40 mm.

Continued on Page 2

**Head and Stem Combinations:
Stainless Steel Titanium Femoral Heads**

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	CPT stainless steel head (6 degree)	PROTASUL-S30 head (12/14 taper)	PROTASUL-S30 head (8/10 taper)	Testa Jupe head (12/14 taper)	Tete a Jupe Acier Inox Stainless Steel 316L Cone 12/14 D28	Tete a Jupe Inox HTA Cone 8/10 D28	TRIBOSUL head (12/14 taper)	TRIBOSUL head (14/16 taper)	Zimtron stainless steel femoral head (12/14 taper)	Zimtron stainless steel femoral head (6 degree)
Fitmore Hip Stem	✓									
FracSure Hip system	✓									
Heritage stem (XX-8095-XXX-XX)										
Intermedia stem	✓			✓			✓			
Lamella stem	✓						✓			
M/L Taper stem										
M/L Taper Stem with Kinectiv Technology										
MARBURG stem	✓						✓			
Mayo 12/14 stem										
Mayo 6 Deg stem										
Mercure/F23 cementless stem										
Mercure/F23 cemented stem									✓	
METABLOC stem	✓						✓			
Minima cemented stem		✓					✓			
Modular Austin Moore stem										
Modular Calcar Replacement stem										
Moore monoblock stem										
MOST Options system	✓									
MS-30 stem	✓						✓			
Muller stem (XX-8050-051/056-XX)									✓	
Müller Straight stem PROTASUL-10	✓						✓			
Müller Straight stem PROTASUL-S30	✓						✓			
Müller Straight stem PROTASUL-S30 (monoblock)										
MultiLock stem										
Natural-Hip LD stem	✓									
Natural-Hip Short Narrow Neck stem	✓									
Natural-Hip System	✓									
Omnicares stem	✓			✓			✓			
OPTAN stem	✓						✓			
Orphee uncemented HA stem	✓			✓			✓			
PF Primary stem	✓			✓			✓			
PF Revision stem	✓			✓			✓			
PFM-Revision stem	✓						✓			
Precedent Revision Hip system	✓						✓			
PYRAMID stem	✓						✓			
REVITAN PRIMARY stem	✓						✓			
REVITAN straight/curved stem *	✓						✓			
Segmental Proximal Femur										
Segmental TM Proximal Femur										

* REVITAN straight stems in 14mm diameter cannot be used with femoral heads longer than +4mm.

Continued on Page 3

Head and Stem Combinations: Stainless Steel Titanium Femoral Heads

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	CPT stainless steel head (6 degree)	PROTASUL-S30 head (12/14 taper)	PROTASUL-S30 head (8/10 taper)	Testa Jupe head (12/14 taper)	Tete a Jupe Acier Inox Stainless Steel 316L Cone 12/14 D28	Tete a Jupe Inox HTA Cone 8/10 D28	TRIBOSUL head (12/14 taper)	TRIBOSUL head (14/16 taper)	Zimtron stainless steel femoral head (12/14 taper)	Zimtron stainless steel femoral head (6 degree)
SF stem										
SKT stem										
SL Müller stem (cemented)		✓					✓			
STELCOR stem										
TCM stem		✓					✓			
Thompson monoblock stem										
Trabecular Metal primary hip stem										
VARIALL stem (except ALLOCLASSIC VARIALL SLV stem)		✓					✓			
VerSys Epoch RNL Stem										
VerSys Epoch Stem										
VerSys System **										
VIRTEC stem		✓					✓			
VIRTEC, Long stem		✓					✓			
Wagner Cone prosthesis stem		✓					✓			
Wagner SL Revision stem ***		✓					✓			
Weber/Stühmer and Weber Curved stem		✓					✓			
ZMR (Body)										

** VerSys Reduced Metaphyseal (RM) stems cannot be used with femoral heads longer than +7mm.

*** Wagner SL Revision Stem in 14mm and 15mm diameters cannot be used with femoral heads longer than +4mm.