

**Acetabular Reconstruction: Trabecular Metal Augments with Acetabular Shells and Cages**

- Indicates "Combination has been assessed as functionally compatible by Zimmer Biomet"
- Indicates "Combination has either not been assessed or is not functionally compatible"

	Trabecular Metal™ Augments and Restrictors	Trabecular Metal™ Buttress and Shim Augments
Continuum® Acetabular System	✓	✓
G7® OsseoTi Shells	✓	✓
Trabecular Metal™ Modular Acetabular System	✓	✓
Trabecular Metal™ Revision Shell	✓	✓

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.