

- Indicates "Combination has been assessed as functionally compatible by Zimmer Biomet"
- Indicates "Combination has either not been assessed or is not functionally compatible"

CoCr Femoral Heads

All heads listed are compatible with the correspondingly sized liners with no limitations.

		Longevity IT XLPE Liners, Neutral & Vivacit-E Liners, Neutral	Longevity IT XLPE Liners, 7mm Offset, Oblique & Vivacit-E Liners, Elevated	Trilogy Longevity Liners, Neutral	Trilogy Longevity Liners, 10° Elevated	Longevity IT XLPE Liners, Elevated	Trilogy Longevity Liners, 7mm Offset	Trilogy PE Liners, Neutral, 10° Elevated, 7mm Offset	Trilogy PE 20° Elevated & Trilogy Longevity 20° Elevated	Trilogy PE, Trilogy Longevity, Longevity IT XLPE Liners, Constrained	Durasul PE Alpha Inserts (01.00013.XXX)	Sulene PE Alpha Inserts (01.00010.XXX)	Original M.E. Müller Low Profile Cup Sulene PE (63.XX.XX)	Original M.E. Müller Low Profile Cup Durasul XLPE (01.00XXX.XXX)	Brunswick Cup Sulene PE (01.00134.10X)	ZCA All-Poly Acetabular Cup, PE (00-8005-XXX-XX)	ZCA All-Poly Acetabular Cup, Longevity (00-8065-XXX-XX)	TMARS Longevity Revision Shell Liner, neutral (00-7105-0XXX-XX)	TMARS Longevity Revision Shell Liner, 10° oblique (00-7110-0XXX-XX)	TMARS Cemented Constrained Liner (00-7115-0XXX-XX)	
010001010	COCR FEM HD 28MM TYPE 1 +12MM																				
010001011	COCR FEM HD 28MM TYPE 1 +9MM																				
010001012	COCR FEM HD 28MM TYPE 1 +6MM																				
010001013	COCR FEM HD 28MM TYPE 1 +3MM																				
010001014	COCR FEM HD 28MM TYPE 1 STD																				
010001015	COCR FEM HD 28MM TYPE 1 -3MM																				
010001016	COCR FEM HD 28MM TYPE 1 -6MM																				
010001017	COCR FEM HD 32MM TYPE 1 +12MM																				
010001018	COCR FEM HD 32MM TYPE 1 +9MM																				
010001019	COCR FEM HD 32MM TYPE 1 +6MM																				
010001020	COCR FEM HD 32MM TYPE 1 +3MM																				
010001021	COCR FEM HD 32MM TYPE 1 STD																				
010001022	COCR FEM HD 32MM TYPE 1 -3MM																				
010001023	COCR FEM HD 32MM TYPE 1 -6MM																				
010001024	COCR FEM HD 36MM TYPE 1 +12MM																				
010001025	COCR FEM HD 36MM TYPE 1 +9MM																				
010001026	COCR FEM HD 36MM TYPE 1 +6MM																				
010001027	COCR FEM HD 36MM TYPE 1 +3MM																				
010001028	COCR FEM HD 36MM TYPE 1 STD																				
010001029	COCR FEM HD 36MM TYPE 1 -3MM																				
010001030	COCR FEM HD 36MM TYPE 1 -6MM																				
010001031	TYPE 1 COCR FEM HEAD 40MM +12MM, G7																				
010001032	TYPE 1 COCR FEM HEAD 40MM +9MM, G7	✓		✓																	
010001033	TYPE 1 COCR FEM HEAD 40MM +6MM, G7	✓		✓																	
010001034	TYPE 1 COCR FEM HEAD 40MM +3MM, G7	✓		✓																	
010001035	TYPE 1 COCR FEM HEAD 40MM STD, G7	✓		✓																	
010001036	TYPE 1 COCR FEM HEAD 40MM -3MM, G7	✓		✓																	
010001037	TYPE 1 COCR FEM HEAD 40MM -6MM, G7	✓		✓																	
010001048	COCR FEM HD 28MM TP 12/14+10.5										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
010001049	COCR FEM HD 28MM TP 12/14 +7MM										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
010001050	COCR FEM HD 28MM TP 12/14 +3.5										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
010001051	COCR FEM HD 28MM TP 12/14 STD										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
010001052	COCR FEM HD 28MM TP 12/14 -3.5										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
010001053	COCR FEM HD 32MM TP 12/14 +12										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
010001054	COCR FEM HD 32MM TP 12/14 +8MM										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

CoCr Femoral Heads (Cont.)

All heads listed are compatible with the correspondingly sized liners with no limitations.

		Longevity IT XLPE Liners, Neutral & Vivacit-E Liners, Neutral	Longevity IT XLPE Liners, 7mm Offset, Oblique & Vivacit-E Liners, Elevated	Trilogy Longevity Liners, Neutral	Trilogy Longevity Liners, 10° Elevated	Longevity IT XLPE Liners, Elevated	Trilogy Longevity Liners, 7mm Offset	Trilogy PE Liners, Neutral, 10° Elevated, 7mm Offset	Trilogy PE 20° Elevated & Trilogy Longevity 20° Elevated	Trilogy PE, Trilogy Longevity, Longevity IT XLPE Liners, Constrained	Durasul PE Alpha Inserts (01.00013.XXX)	Sulene PE Alpha Inserts (01.00010.XXX)	Original M.E. Müller Low Profile Cup Sulene PE (63.XX.XX)	Original M.E. Müller Low Profile Cup Durasul XLPE (01.000XX.XXX)	Brunswick Cup Sulene PE (01.00134.10X)	ZCA All-Poly Acetabular Cup, PE (00-8005-XXX-XX)	ZCA All-Poly Acetabular Cup, Longevity (00-8065-XXX-XX)	TMARS Longevity Revision Shell Liner, neutral (00-7105-0XXX-XX)	TMARS Longevity Revision Shell Liner, 10° oblique (00-7110-0XXX-XX)	TMARS Cemented Constrained Liner (00-7115-0XXX-XX)
010001055	COCR FEM HD 32MM TP 12/14 +4MM										✓	✓	✓	✓						
010001056	COCR FEM HD 32MM TP 12/14 STD										✓	✓	✓	✓						
010001057	COCR FEM HD 32MM TP 12/14 -4MM										✓	✓	✓	✓						
010001058	COCR FEM HD 36MM TP 12/14 +12										✓	✓	✓	✓						
010001059	COCR FEM HD 36MM TP 12/14 +8MM										✓	✓	✓	✓						
010001060	COCR FEM HD 36MM TP 12/14 +4MM										✓	✓	✓	✓						
010001061	COCR FEM HD 36MM TP 12/14 STD										✓	✓	✓	✓						
010001062	COCR FEM HD 36MM TP 12/14 -4MM										✓	✓	✓	✓						
010001063	12/14 COCR FEM HEAD 40MM +12MM, G7	✓		✓															✓	✓
010001064	12/14 COCR FEM HEAD 40MM +8MM, G7	✓		✓															✓	✓
010001065	12/14 COCR FEM HEAD 40MM +4MM, G7	✓		✓															✓	✓
010001066	12/14 COCR FEM HEAD 40MM STD, G7	✓		✓															✓	✓
010001067	12/14 COCR FEM HEAD 40MM -4MM, G7	✓		✓															✓	✓
11-363660	TYPE 1 COCR MOD HEAD 36MM -6MM	✓	✓	✓	✓	✓					✓	✓	✓	✓					✓	✓
11-363661	TYPE 1 COCR MOD HEAD 36MM -3MM	✓	✓	✓	✓	✓					✓	✓	✓	✓					✓	✓
11-363662	TYPE 1 COCR MOD HEAD 36MM STD	✓	✓	✓	✓	✓					✓	✓	✓	✓					✓	✓
11-363663	TYPE 1 COCR MOD HEAD 36MM +3MM	✓	✓	✓	✓	✓					✓	✓	✓	✓					✓	✓
11-363664	TYPE 1 COCR MOD HEAD 36MM +6MM	✓	✓	✓	✓	✓					✓	✓	✓	✓					✓	✓
11-363665	TYPE 1 COCR MOD HEAD 36MM +9MM	✓	✓	✓	✓	✓					✓	✓	✓	✓					✓	✓
11-363666	TYPE 1 COCR MOD HEAD 36MM +12MM	✓	✓	✓	✓	✓					✓	✓	✓	✓					✓	✓
163638	TYPE 1 COCR MOD HEAD NO SKIRT 28MM +6MM	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓
163651	TYPE 1 COCR MOD HEAD 22.2 MM STD				✓	✓	✓	✓	✓										✓	✓
163652	TYPE 1 COCR MOD HEAD 22.2 MM -5MM				✓	✓	✓	✓	✓										✓	✓
163653	TYPE 1 COCR MOD HEAD 22.2 MM -3MM				✓	✓	✓	✓	✓										✓	✓
163660	TYPE 1 COCR MOD HEAD 28 MM -6MM	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓					✓	✓
163661	TYPE 1 COCR MOD HEAD 28 MM -3MM	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓					✓	✓
163662	TYPE 1 COCR MOD HEAD 28 MM STD	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓					✓	✓
163663	TYPE 1 COCR MOD HEAD 28 MM +3MM	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓					✓	✓
163664	TYPE 1 COCR MOD HEAD 28 MM +6MM	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓					✓	✓
163665	TYPE 1 COCR MOD HEAD 28 MM +9MM	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓					✓	✓
163666	TYPE 1 COCR MOD HEAD 28 MM +12MM	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓					✓	✓

CoCr Femoral Heads (Cont.)

All heads listed are compatible with the correspondingly sized liners with no limitations.

		Longevity IT XLPE Liners, Neutral & Vivacit-E Liners, Neutral	Longevity IT XLPE Liners, 7mm Offset, Oblique & Vivacit-E Liners, Elevated	Trilogy Longevity Liners, Neutral	Trilogy Longevity Liners, 10° Elevated	Longevity IT XLPE Liners, Elevated	Trilogy Longevity Liners, 7mm Offset	Trilogy PE Liners, Neutral, 10° Elevated, 7mm Offset	Trilogy PE 20° Elevated & Trilogy Longevity 20° Elevated	Trilogy PE, Trilogy Longevity, Longevity IT XLPE Liners, Constrained	Durasul PE Alpha Inserts (01.00013.XXX)	Sulene PE Alpha Inserts (01.00010.XXX)	Original M.E. Müller Low Profile Cup Sulene PE (63.XX.XX)	Original M.E. Müller Low Profile Cup Durasul XLPE (01.000XX.XXX)	Brunswick Cup Sulene PE (01.00134.10X)	ZCA All-Poly Acetabular Cup, PE (00-8005-XXX-XX)	ZCA All-Poly Acetabular Cup, Longevity (00-8065-XXX-XX)	TMARS Longevity Revision Shell Liner, neutral (00-7105-0XXX-XX)	TMARS Longevity Revision Shell Liner, 10° oblique (00-7110-0XXX-XX)	TMARS Cemented Constrained Liner (00-7115-0XXX-XX)
163667	TYPE 1 COCR MOD HEAD 32 MM -6MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
163668	TYPE 1 COCR MOD HEAD 32 MM -3MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
163669	TYPE 1 COCR MOD HEAD 32 MM STD	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
163670	TYPE 1 COCR MOD HEAD 32 MM +3MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
163671	TYPE 1 COCR MOD HEAD 32 MM +6MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
163672	TYPE 1 COCR MOD HEAD 32 MM +9MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
163673	TYPE 1 COCR MOD HEAD 32 MM +12MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
163674	TYPE 1 COCR MOD HEAD NO SKIRT 32MM +6MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
164131	12/14 COCR 5DEG 42TAPER MOD HEAD 22.2MM +2MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
164132	12/14 COCR 5DEG 42TAPER MOD HEAD 22.2MM 0MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
164133	12/14 COCR 5DEG 42TAPER MOD HEAD 22.2MM -2MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0839	TYPE 1 COCR BIOMET FEM HEAD 36MM -6MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0840	TYPE 1 COCR BIOMET FEM HEAD 36MM -3MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0841	TYPE 1 COCR BIOMET FEM HEAD 36MM STD, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0842	TYPE 1 COCR BIOMET FEM HEAD 36MM +3MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0843	TYPE 1 COCR BIOMET FEM HEAD 36MM +6MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0863	TYPE 1 COCR BIOMET FEM HEAD 28MM -6MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0864	TYPE 1 COCR BIOMET FEM HEAD 28MM -3MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0865	TYPE 1 COCR BIOMET HEAD 28MM STD, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0866	TYPE 1 COCR BIOMET FEM HEAD 28MM +3MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0867	TYPE 1 COCR BIOMET FEM HEAD 28MM +6MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0870	TYPE 1 COCR BIOMET FEM HEAD 32MM -6MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0871	TYPE 1 COCR BIOMET FEM HEAD 32MM -3MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0872	TYPE 1 COCR BIOMET FEM HEAD 32MM STD, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0873	TYPE 1 COCR BIOMET FEM HEAD 32MM +3MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0874	TYPE 1 COCR BIOMET FEM HEAD 32MM +6MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0877	12/14 COCR TPR FEM HEAD 28MM -3.5MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0878	12/14 COCR TPR HEAD 28MM STD, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0879	12/14 COCR TPR HEAD 28MM +3.5MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0882	12/14 COCR TPR HEAD 32MM -4MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0883	12/14 COCR TPR HEAD 32MM STD, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0884	12/14 COCR TPR HEAD 32MM +4MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0887	12/14 COCR TPR HEAD 36MM -4MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0888	12/14 COCR TPR HEAD 36MM STD, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0889	12/14 COCR TPR HEAD 36MM +4MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0890	12/14 COCR TPR HEAD 36MM +8MM, EXCEED	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1020	12/14 COCR MOD HEAD 40MM -4MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1021	12/14 COCR MOD HEAD 40MM 0MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1022	12/14 COCR MOD HEAD 40MM +4MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1023	12/14 COCR MOD HEAD 40MM +8MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1024	12/14 COCR MOD HEAD 40MM +12MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1030	TYPE 1 COCR MOD HEAD 40MM -6MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

CoCr Femoral Heads (Cont.)

All heads listed are compatible with the correspondingly sized liners with no limitations.

		Longevity IT XLPE Liners, Neutral & Vivacit-E Liners, Neutral	Longevity IT XLPE Liners, 7mm Offset, Oblique & Vivacit-E Liners, Elevated	Trilogy Longevity Liners, Neutral	Trilogy Longevity Liners, 10° Elevated	Longevity IT XLPE Liners, Elevated	Trilogy Longevity Liners, 7mm Offset	Trilogy PE Liners, Neutral, 10° Elevated, 7mm Offset	Trilogy PE 20° Elevated & Trilogy Longevity 20° Elevated	Trilogy PE, Trilogy Longevity, Longevity IT XLPE Liners, Constrained	Durasul PE Alpha Inserts (01.00013.XXX)	Sulene PE Alpha Inserts (01.00010.XXX)	Original M.E. Müller Low Profile Cup Sulene PE (63.XX.XX)	Original M.E. Müller Low Profile Cup Durasul XLPE (01.00XXXX.XXX)	Brunswick Cup Sulene PE (01.00134.10X)	ZCA All-Poly Acetabular Cup, PE (00-8005-XXX-XX)	ZCA All-Poly Acetabular Cup, Longevity (00-8065-XXX-XX)	TMARS Longevity Revision Shell Liner, neutral (00-7105-0XXX-XX)	TMARS Longevity Revision Shell Liner, 10° oblique (00-7110-0XXX-XX)	TMARS Cemented Constrained Liner (00-7115-0XXX-XX)
650-1031	TYPE 1 COCR MOD HEAD 40MM -3MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
650-1032	TYPE 1 COCR MOD HEAD 40MM 0MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
650-1033	TYPE 1 COCR MOD HEAD 40MM +3MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
650-1034	TYPE 1 COCR MOD HEAD 40MM +6MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
650-1035	TYPE 1 COCR MOD HEAD 40MM +9MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
650-1036	TYPE 1 COCR MOD HEAD 40MM +12MM, EXCEED ABT M2A-40	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
BH400140	12/14 COCR MOD HEAD 40MM -6MM	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
BH400240	12/14 COCR MOD HEAD 40MM -3MM	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
BH400340	12/14 COCR MOD HEAD 40MM STD	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
BH400440	12/14 COCR MOD HEAD 40MM +3MM	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
BH400540	12/14 COCR MOD HEAD 40MM +6MM	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
BH400640	12/14 COCR MOD HEAD 40MM +9MM	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
P0206228	12/14 COCR MOD FEM HEAD 28MM +10.5MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
P0206232	12/14 COCR MOD FEM HEAD 32MM +12MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
P0206C28	12/14 COCR HEAD SHT NK 28MM -3.5MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
P0206C32	12/14 COCR HEAD SHT NK 32MM -4MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
P0206E28	12/14 COCR HEAD XL NK 28MM +7MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
P0206E32	12/14 COCR HEAD XL NK 32MM +8MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
P0206L28	12/14 COCR HEAD LNG NK 28MM +3.5MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
P0206L32	12/14 COCR HEAD LNG NK 32MM +4MM	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
P0206M28	12/14 COCR HEAD MD NK 28MM STD	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
P0206M32	12/14 COCR HEAD MD NK 32MM STD	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
P0206028	CrCo FEMORAL HEAD 5°42' Ø 28/-5/12-14 EXTRA SHORT NECK	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
P0206032	CrCo FEMORAL HEAD 5°42' Ø 32/-6/12-14 EXTRA SHORT NECK	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
S001140	TYPE 1 COCR MOD HEAD 40MM STD, SELEX/MAGNUM	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
S031140	TYPE 1 COCR MOD HEAD 40MM -3MM, SELEX/MAGNUM	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
S061140	TYPE 1 COCR MOD HEAD 40MM -6MM, SELEX/MAGNUM	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
S331140	TYPE 1 COCR MOD HEAD 40MM +3MM, SELEX/MAGNUM	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
S661140	TYPE 1 COCR MOD HEAD 40MM +6MM, SELEX/MAGNUM	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
S991140	TYPE 1 COCR MOD HEAD 40MM +9MM, SELEX/MAGNUM	✓	✓	✓	✓	✓	✓	✓	✓									✓	✓	
P0206C22	12/14 COCR HEAD SHORT NECK 22.2MM -2MM		✓	✓	✓	✓	✓	✓	✓							✓	✓			
P0206L22	12/14 COCR HEAD LONG NECK 22.2MM -2MM		✓	✓	✓	✓	✓	✓	✓							✓	✓			
P0206M22	12/14 COCR HEAD LONG NECK 22.2MM 0MM		✓	✓	✓	✓	✓	✓	✓							✓	✓			
163620	26MM DIA COCR MOD HD -6MM NK NO SKIRT		✓	✓	✓	✓	✓	✓	✓							✓	✓			
163621	26MM DIA COCR MOD HD -3MM NK NO SKIRT		✓	✓	✓	✓	✓	✓	✓							✓	✓			
163622	26MM DIA COCR MOD HD STD NK NO SKIRT		✓	✓	✓	✓	✓	✓	✓							✓	✓			
163623	26MM DIA COCR MOD HD +3MM NK NO SKIRT		✓	✓	✓	✓	✓	✓	✓							✓	✓			
163624	26MM DIA CO CR MOD HEAD +6		✓	✓	✓	✓	✓	✓	✓							✓	✓			
163625	26MM DIA COCR MOD HD +9MM NK SKIRTED		✓	✓	✓	✓	✓	✓	✓							✓	✓			
163626	26MM DIA COCR MOD HD +12MM NK SKIRTED		✓	✓	✓	✓	✓	✓	✓							✓	✓			
163627	26MM DIA COCR MOD HD +6MM NK NO SKIRT		✓	✓	✓	✓	✓	✓	✓							✓	✓			

Ceramic Femoral Heads

All heads listed are compatible with the correspondingly sized liners with no limitations.

		Longevity IT XLPE Liners, Neutral & Vivacit-E Liners, Neutral	Longevity IT XLPE Liners, 7mm Offset, Oblique & Vivacit-E Liners, Elevated	Trilogy Longevity Liners, Neutral	Trilogy Longevity Liners, 10° Elevated	Longevity IT XLPE Liners, Elevated	Trilogy Longevity Liners, 7mm Offset	Trilogy PE Liners, Neutral, 10° Elevated, 7mm Offset	Trilogy PE 20° Elevated & Trilogy Longevity 20° Elevated	Trilogy PE, Trilogy Longevity, Longevity IT XLPE Liners, Constrained	Durasul PE Alpha Inserts (01.00013.XXX)	Sulene PE Alpha Inserts (01.00010.XXX)	Original M.E. Müller Low Profile Cup Sulene PE (63.XX.XX)	Original M.E. Müller Low Profile Cup Durasul XLPE (01.000XX.XXX)	Brunswick Cup Sulene PE (01.00134.10X)	ZCA All-Poly Acetabular Cup, PE (00-8005-XXX-XX)	ZCA All-Poly Acetabular Cup, Longevity (00-8065-XXX-XX)	TMARS Longevity Revision Shell Liner, neutral (00-7105-0XXX-XX)	TMARS Longevity Revision Shell Liner, 10° oblique (00-7110-0XXX-XX)	TMARS Cemented Constrained Liner* (00-7115-0XXX-XX)
12-115109	TYPE 1 CERAMIC FEM HEAD 28MM -3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
12-115110	TYPE 1 CERAMIC FEM HEAD 28MM -0M, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
12-115111	TYPE 1 CERAMIC FEM HEAD 28MM +3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
12-115114	TYPE 1 CERAMIC FEM HEAD 32MM -3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
12-115115	TYPE 1 CERAMIC FEM HEAD 32MM -0M, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
12-115116	TYPE 1 CERAMIC FEM HEAD 32MM +3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
12-115117	TYPE 1 CERAMIC FEM HEAD 32M +6MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
12-115120	TYPE 1 CERAMIC FEM HEAD 36MM -3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
12-115121	TYPE 1 CERAMIC FEM HEAD 36MM +0MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
12-115122	TYPE 1 CERAMIC FEM HEAD 36MM +3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
12-115123	TYPE 1 CERAMIC FEM HEAD 36MM +6MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
164135	TYPE 1 CERAMIC FEM HEAD 28MM -3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
164136	TYPE 1 CERAMIC FEM HEAD 28MM +0MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
164137	TYPE 1 CERAMIC FEM HEAD 28MM +0MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
164185	TYPE 1 CERAMIC FEM HEAD 32MM -3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
164186	TYPE 1 CERAMIC FEM HEAD 32MM +0MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
164187	TYPE 1 CERAMIC FEM HEAD 32MM +3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0660	TYPE 1 CERAMIC FEM HEAD 36MM -3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
650-0661	TYPE 1 CERAMIC FEM HEAD 36MM +0MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
650-0662	TYPE 1 CERAMIC FEM HEAD 36MM +3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
650-0663	TYPE 1 CERAMIC FEM HEAD 36MM +6MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
650-0667	12/14 CERAMIC FEM HEAD 36MM +8MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓								
650-0830	12/14 CERAMIC FEM HEAD 28MM -3.5MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0831	12/14 CERAMIC FEM HEAD 28MM +0MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0832	12/14 CERAMIC FEM HEAD 28MM +3.5MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0833	12/14 CERAMIC FEM HEAD 32MM -4.0MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0834	12/14 CERAMIC FEM HEAD 32MM +0MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0835	12/14 CERAMIC FEM HEAD 32MM +4MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0836	12/14 CERAMIC FEM HEAD 36MM -4.0MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0837	12/14 CERAMIC FEM HEAD 36MM +0MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-0838	12/14 CERAMIC FEM HEAD 36MM +4MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1055	OPTION CERAMIC HEAD 28MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1056	OPTION CERAMIC HEAD 32MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1057	OPTION CERAMIC HEAD 36MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1058	OPTION CERAMIC HEAD 40MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1157	TYPE1 CERAMIC FEM HEAD 28MM +3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1158	TYPE1 CERAMIC FEM HEAD 28MM +0MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1159	TYPE1 CERAMIC FEM HEAD 28MM -3MM, DELTA	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1160	DELTA CER FEM HD 32/+6MM T1	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1161	DELTA CER FEM HD 32/+3MM T1	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1162	DELTA CER FEM HD 32/0MM T1	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
650-1163	DELTA CER FEM HD 32/-3MM T1	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

*TMARS Cemented Constrained Liners are not available for the use with ceramic heads in the US

Stainless Steel Femoral Heads All heads listed are compatible with the correspondingly sized liners with no limitations.		Longevity IT XLPE Liners, Neutral & Vivacit-E Liners, Neutral	Longevity IT XLPE Liners, 7mm Offset, Oblique & Vivacit-E Liners, Elevated	Trilogy Longevity Liners, Neutral	Trilogy Longevity Liners, 10° Elevated	Longevity IT XLPE Liners, Elevated	Trilogy Longevity Liners, 7mm Offset	Trilogy PE Liners, Neutral, 10° Elevated, 7mm Offset	Trilogy PE 20° Elevated & Trilogy Longevity 20° Elevated	Trilogy PE, Trilogy Longevity, Longevity IT XLPE Liners, Constrained	Durasul PE Alpha Inserts (01.00013.XXX)	Sulene PE Alpha Inserts (01.00010.XXX)	Original M.E. Müller Low Profile Cup Sulene PE (63.XX.XX)	Original M.E. Müller Low Profile Cup Durasul XLPE (01.000XX.XXX)	Brunswick Cup Sulene PE (01.00134.10X)	ZCA All-Poly Acetabular Cup, PE (00-8005-XXX-XX)	ZCA All-Poly Acetabular Cup, Longevity (00-8065-XXX-XX)	TMARS Longevity Revision Shell Liner, neutral (00-7105-0XXX-XX)	TMARS Longevity Revision Shell Liner, 10° oblique (00-7110-0XXX-XX)	TMARS Cemented Constrained Liner (00-7115-0XXX-XX)
		<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<
78-620-22952	12/14 STSL FEM HEAD 22.2 M +0MM	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<
78-620-22962	12/14 STSL FEM HEAD 22.2 M +4MM	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<
P0201C28	12/14 STSTL FEM HEAD 28MM -3.5MM	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<
P0201C32	12/14 STSTL FEM HEAD 32MM -4MM	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<
P0201E28	12/14 STSTL FEM HEAD 28MM +7MM	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<
P0201E32	12/14 STSL FEM HEAD 32MM/NL +8MM	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<
P0201L28	12/14 STSTL FEM HEAD 28MM +3.5MM	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<
P0201L32	12/14 STSTL FEM HEAD 32MM +4MM	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<
P0201M28	12/14 STSTL FEM HEAD 28MM 0MM	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<
P0201M32	12/14 STSTL FEM HEAD 32MM 0MM	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<	<