

- ✔ Indicates that the combination CAN be used
- ✘ Indicates that the combination CANNOT be used

System Designation	ZNN	ZNN	ZNN	ZNN	ZNN	ZNN	ZNN	ZNN	ZNN	ZNN	ZNN	ZNN	ZNN	ZNN	ZNN	ZNN
Note: This table represents functional compatibility and does not indicate the regulatory approval status. xx or xxx = These are placeholders that represent the item numbers for the range of device sizes.	4.0mm Cortical Screw, Partially Threaded 47/00-2483-xxx-40	4.0mm Cortical Screw, Fixed Angle 47/00-2484-xxx-40	5.0mm Cortical Screw, Partially Threaded, 20-60mm 47/00-2483-0xx-50	5.0mm Cortical Screw, Partially Threaded, 65-100mm 47/00-2483-xxx-50	5.0mm Cortical Screw, Fixed Angle 20-60mm 47/00-2484-0xx-50	5.0mm Cortical Screw, Fixed Angle 65-100mm 47/00-2484-xxx-50	6.0mm Cancellous Screw, Partially Threaded 47/00-2483-xxx-60	6.0mm Cancellous Screw, Fixed Angle, 30-45mm 47/00-2484-xxx-60	6.0mm Cancellous Screw, Fixed Angle, 50-110mm 47/00-2484-xxx-60	6.0mm Cancellous Screw, Fixed Angle, 115-140mm 47/00-2484-xxx-60	Antegrade Femoral Nail Cap, 0-15mm Height 47-2487-001-xx	Retrograde Femoral Nail Cap, 0-15mm Height 47-2487-008-xx	Tibial Nail Cap, 0-15mm Height 47-2487-005-xx	Cephalomedullary Nail Cap 0-15mm Height 00/47-2487-002-xx	10.5mm Lag Screw, 70-130mm Length 00/47-2485-xxx-10	Cephalomedullary Nail Set Screw 47-2493-000-00
Zimmer Natural Nail (ZNN) System Piriformis Fossa Femoral Nail-Green, 9.3-15mm Diameter, 24-48cm Length 47-2491-xx0-09/10/11/12/13/14/15	✘	✘	✔	✔	✔	✔	✔	✘	✔	✔	✔	✘	✘	✘	✘	✘
Zimmer Natural Nail (ZNN) System Greater Trochanter Femoral Nail-Green, Right, 8.3mm Diameter, 24-40cm Length 47-2492-xx0-08	✔	✔	✔	✔	✔	✔	✔	✘	✔	✔	✔	✘	✘	✘	✘	✘
Zimmer Natural Nail (ZNN) System Greater Trochanter Femoral Nail-Green, Right, 9.3-14mm Diameter, 24-48cm Length 47-2492-xx0-09/10/11/12/13/14	✘	✘	✔	✔	✔	✔	✔	✘	✔	✔	✔	✘	✘	✘	✘	✘
Zimmer Natural Nail (ZNN) System Greater Trochanter Femoral Nail-Green, Left, 8.3mm Diameter, 24-40cm Length 47-2492-xx1-08	✔	✔	✔	✔	✔	✔	✔	✘	✔	✔	✔	✘	✘	✘	✘	✘
Zimmer Natural Nail (ZNN) System Greater Trochanter Femoral Nail-Green, Left, 9.3-14mm Diameter, 24-48cm Length 47-2492-xx1-09/10/11/12/13/14	✘	✘	✔	✔	✔	✔	✔	✘	✔	✔	✔	✘	✘	✘	✘	✘
Zimmer Natural Nail (ZNN) System Retrograde Femoral Nail Short-Purple 10mm Diameter, 16-44cm Length 47-2494-xx0-10	✘	✘	✔	✘	✔	✘	✘	✔	✔	✘	✘	✔	✘	✘	✘	✘
Zimmer Natural Nail (ZNN) System Retrograde Femoral Nail Long-Purple 11.5-13mm Diameter, 16-44cm Length 47-2494-xx0-11/13	✘	✘	✔	✘	✔	✘	✘	✔	✔	✘	✘	✔	✘	✘	✘	✘
Zimmer Natural Nail (ZNN) System Tibial Nail-Yellow 8.3mm Diameter 24-44cm Length 47-2495-xx0-08	✔	✔	✔	✔	✔	✔	✘	✘	✘	✘	✘	✘	✔	✘	✘	✘
Zimmer Natural Nail (ZNN) System Tibial Nail-Yellow 9.3-14mm Diameter 24-44cm Length 47-2495-xx0-09/10/11/12/13/14	✘	✘	✔	✔	✔	✔	✘	✘	✘	✘	✘	✘	✔	✘	✘	✘
Zimmer Natural Nail (ZNN) System Cephalomedullary Femoral Nail-Short, Blue, 125-135 Deg, Left/Right, 10-14.5mm Diameter 21.5cm Length 47-2493-2xx-10/11/13/14	✘	✘	✔	✔	✔	✔	✘	✘	✘	✘	✘	✘	✘	✔	✔	✔
Zimmer Natural Nail (ZNN) System Cephalomedullary Femoral Nail-Asia Short, Blue, 125-130 Deg, Left/Right, 9.3-14.5mm Diameter 18cm Length 47-2493-18x-09/10/11/13/14	✘	✘	✔	✔	✔	✔	✘	✘	✘	✘	✘	✘	✘	✔	✔	✔
Zimmer Natural Nail (ZNN) System Cephalomedullary Femoral Nail-Long, Blue, 125-135 Deg, Left/Right, 10-13mm Diameter 30-38cm Length 47-2493-3xx-10/11/13	✘	✘	✔	✔	✔	✔	✘	✘	✘	✘	✘	✘	✘	✔	✔	✔
Zimmer Natural Nail (ZNN) System Cephalomedullary Femoral Nail-Long, Blue, 125-135 Deg, Left/Right, 10-13mm Diameter 40-48cm Length 47-2493-4xx-10/11/13	✘	✘	✔	✔	✔	✔	✘	✘	✘	✘	✘	✘	✘	✔	✔	✔

	ZNN	ZNN
Cortical Nut* 47-2488-000-01	✔	✔
Cortical Washer* 47/00-2488-000-02	✔	✔

*Cortical Nut and Washer should only be used with the 6.0mm Fixed Angle screws when also used with the retrograde femoral nails.

Zimmer Natural Nail System 6.0 mm Cancellous Screw, Fixed Angle, 30-110mm, 47/00-2484-xxx-60	✔	✔
--	---	---