

Zimmer Universal Locking (ULS) System Compatibility

- █ Indicates that the combination CAN be used
- █ Indicates that the combination CANNOT be used

Material		Titanium					Stainless Steel					Titanium				Stainless Steel			
		Screw Head Size		2.7mm heads			Titanium		3.5mm heads			4mm heads							
System Designation		ULS	ULS	ULS	ZPS	ZPS	ULS	ULS	ZPLP	ZPLP	Peri	ZPS	ZPS	ULS	ULS	ZPS	ZPS		
Note: This table represents functional compatibility and does not indicate the regulatory approval status. Be sure to check for proper Regulatory approvals in your country prior to using any products found on this website. Some devices may not be currently licensed with Health Canada. Some device compatibilities may not be approved for use by Health Canada. xx or xxx = These are placeholders that represent the item numbers for the range of device sizes.		2.7mm Locking Screw 47/00-2369-xxx-27	2.7mm Cortical Screw, Self-Tapping 47/00-4627-xxx-01	2.7mm Locking Screw 47/00-4828-xxx-02	2.7mm Cortical Screw 47/00-4827-xxx-00	2.7mm Cortical Screw, Self-Tapping 47/00-4827-xxx-01	3.5mm Locking Screw 47/00-2369-xxx-35	3.5mm Cortical Screw, Self-Tapping 47/00-4935-xxx-01	3.5mm Locking Screw 47/00-2359-xxx-35	3.5mm Cannulated Locking Screw 47/00-2359-xxx-36	3.5mm Periarthral Screw, Self-Tapping 47/00-2347-020/120-xx	3.5mm Cortical Screw 47/00-4835-xxx-00	3.5mm Cortical Screw, Self-Tapping 47/00-4835-xxx-01	4.0mm Cancellous Screw, Partially Threaded 47-00-4640-xxx-00	4.0mm Cancellous Screw, Fully Threaded 47/00-4640-xxx-01	4.0mm Cancellous Partially Threaded Screw 47-00-4840-xxx-01	4.0mm Cancellous Fully Threaded Screw 47/00-4840-xxx-00		
STAINLESS STEEL	2.7mm plates	Universal Locking System 2.7mm Locking Dual Compression Plates 47/00-4928-xxx-07			█	█													
		Universal Locking System 2.7mm Locking T-Plates 47/00-4928-xxx-02			█	█	█												
		Universal Locking System 2.7mm Locking L-Plates, and L-Plates Oblique (Left/Right) 47/00-4928-xxx-02			█	█	█												
		Universal Locking System 2.7mm Locking Reconstruction Plates 47/00-4928-xxx-13			█	█	█												
		Universal Locking System 2.7mm Locking 1/4 Tubular Plates 47/00-4928-xxx-03			█	█	█												
	3.5mm plates	Universal Locking System 3.5mm Locking Dual Compression Plates 47/00-4936-xxx-07								█	█	█	█	█			█	█	
		Universal Locking System 3.5mm Locking T-Plates 47/00-4936-035/046-04								█	█	█	█	█			█	█	
		Universal Locking System 3.5mm Locking T-Plates Oblique (Left/Right) 47/00-4936-xxx-05/06								█	█	█	█	█			█	█	
		Universal Locking System 3.5mm Locking 1/3 Tubular Plates 47/00-4936-xxx-03								█	█	█	█	█			█	█	
		Universal Locking System 3.5mm Locking Reconstruction Plates 47/00-4936-xxx-13								█	█	█	█	█			█	█	
TITANIUM	2.7mm plates	Universal Locking System 2.7mm Locking Dual Compression Plates 47/00-4938-xxx-07	█	█															
		Universal Locking System 2.7mm Locking T-Plates 47/00-4938-xxx-02	█	█															
		Universal Locking System 2.7mm Locking L-Plates/L-Plates Oblique (Left/Right) 47/00-4938-xxx-02	█	█															
		Universal Locking System 2.7mm Locking Reconstruction Plates 47/00-4938-xxx-13	█	█															
		Universal Locking System 2.7mm Locking 1/4 Tubular Plates 47/00-4938-xxx-03	█	█															
	3.5mm plates	Universal Locking System 3.5mm Locking Dual Compression Plates 47/00-4946-xxx-07						█	█						█	█			
		Universal Locking System 3.5mm Locking T-Plates 47/00-4946-035/046-04						█	█						█	█			
		Universal Locking System 3.5mm Locking T-Plates Oblique (Left/Right) 47/00-4946-xxx-05/06						█	█						█	█			
		Universal Locking System 3.5mm Locking 1/3 Tubular Plates 47/00-4946-xxx-03						█	█						█	█			
		Universal Locking System 3.5mm Locking Reconstruction Plates 47/00-4946-xxx-13						█	█						█	█			