

## Acetabular Reconstruction: Polyethylene Liners/Cups and Metal Cages/Rings

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	TM Revision Cage	Burch-Schneider Cage	Zimmer Reinforcement Ring with Hook*	Original M.E.Müller Roof Reinforcement Ring	Wagner Reinforcement Cage	ZCA Cage	ZCA Roof Ring
Acora cup (cemented)		✓	✓	✓			
Apollo All-poly (PE) Cemented cup		✓	✓	✓			
BETA DURASUL insert					✓		
BETA METASUL insert					✓		
BETA PE insert					✓		
Brunswick PE cup		✓	✓	✓			
Continuum Hydrocel Acetabular Revision Shell Liner							
Full Profile PE cup		✓	✓	✓			
HG II PE Liner							
HG II XLPE Liner							
Longevity IT Liner							
Low Profile DURASUL cup		✓	✓	✓			
Low Profile METASUL cup		✓	✓	✓			
Low Profile PE cup		✓	✓	✓			
TM Revision Cemented Constrained Liner	✓						
TM Revision Liner	✓						
TM Revision Shell (Liner/Cup/Cage Construct)	✓						
Trilogy AB Alumina Insert							
Trilogy Longevity Liner							
Trilogy PE Liner							
ZCA Cup						✓	✓

\*Formerly commercialized under the name Ganz Ring