

ZIMMER BIOMET

(Please note that only actively sold device pairings are listed)

"X" Indicates "Functional compatibility is approved by Zimmer Biomet"
 "Greyed out" columns Indicates "Combination is not approved by Zimmer Biomet"

Type 1 Femoral Heads	
Item # / Scope Group ID	Brand / Product Family
Freedom Type 1	Freedom Head Type 1
Biomet CoCr Type 1	CoCr Head Type 1
Biolox Delta	BIOLOX® Delta Ceramic Fem Head
ZB Ceramic T1, 22.2*	ZB Ceramic Femoral Heads 22.2mm *
ZB Ceramic T1	ZB Ceramic Femoral Heads 28mm-44mm
ZB Select Ceramic	ZB Select Ceramic Head Implants
Ringloc Heads	RingLoc® Type 1 CoCr Heads
Endo II **	Endo II CoCr Head**
Biolox Delta Option	BIOLOX® Delta Option Head

Type 1 Adaptors	
Item # / Scope Group ID	Brand / Product Family
ZB Select Adaptor T1	ZB Type 1 Select Adaptor
BioloX Option Adaptor	BioloX® Option Adaptor (Type 1)
Endo II Adaptor	Endo II Taper Insert

X	X	X
X	X	X
X		X
X		X
X	X	X
X		X
X		X
X		X
X		X
X		X
X	X	X
X	X	X
X	X	X
X	X	X
X	X	X

SCOPE

Type 1 Stems		
Group ID	Item Number	Description
Arcos One		Arcos One Piece Revision
Arcos One	11-3030XX (00_16,18,20,22)	Arcos Broach Body STD One Piece Revision Femur
Arcos One	11-3031XX (00_16,18,20,22)	Arcos Broach Body STD One Piece Revision Femur
Arcos One	11-3034XX (00_16,18,20,22)	Arcos Calcar Body STD One Piece Revision Femur
Arcos One	11-3040XX (00_16,18,20,22)	Arcos Broach Body HI One Piece Revision Femur
Arcos One	11-3041XX (00_16,18,20,22)	Arcos Broach Body HI One Piece Revision Femur
Arcos One	11-304XXX (200_326)	Arcos Broach Body HI One Piece Revision Femur
Arcos One	11-3044XX (00_16,18,20,22)	Arcos Calcar Body HI One Piece Revision Femur
Arcos Mod		Arcos Revision Bodies
Arcos Mod	11-3011XX	Arcos Femoral Revision Broached Proximal Body Std/Hi Offset
Arcos Mod	11-30120X	Arcos +0 Calcar Proximal Body Standard Offset 60mm
Arcos Mod	11-30121X	Arcos +0 Calcar Proximal Body High Offset 60mm
Arcos Mod	11-30122X	Arcos +10 Calcar Proximal Body Standard Offset 60mm
Arcos Mod	11-30123X	Arcos +10 Calcar Proximal Body High Offset 60mm
Arcos Mod	11-30124X	Arcos +20 Calcar Proximal Body Standard Offset 60mm
Arcos Mod	11-30125X	Arcos +20 Calcar Proximal Body High Offset 60mm
Arcos Mod	11-30130X	Arcos Cone Proximal Body Standard Offset 50/60mm
Arcos Mod	11-30131X	Arcos Cone Proximal Body High Offset 50/60mmmm
Arcos Mod	11-30132X	Arcos Cone Proximal Body Standard Offset 70mm
Arcos Mod	11-30133X	Arcos Cone Proximal Body High Offset 70mm
Arcos Mod	11-30134X	Arcos Cone Proximal Body Standard Offset 80mm
Arcos Mod	11-30135X	Arcos Cone Proximal Body High Offset 80mm
BiMetric		Bi-Metric CoCr Type 1 Stems
BiMetric	CP1524XX (05,06,44,55,57)	BiMetric Stem
BiMetric	CP155XXX (919_920)	BiMetric Interlok
BiMetric	CP160XXX (400_401)	BiMetric Meta
Compress		Compress Proximal Femorals
Compress	1787XX	Compress Proximal Femorals
Echo BiMetric		Echo Bi-Metric PF Stem Type 1
Echo BiMetric	1920XX	Echo Bi-metric Full Proximal Profile - Standard 135°
Echo BiMetric	1921XX	Echo Bi-metric Full Proximal Profile - Lateralized 130°
Echo BiMetric	1924XX	Echo Bi-metric Reduced Proximal Profile - Standard 135°
Echo BiMetric	1925XX (08_21)	Echo Bi-metric Reduced Proximal Profile - Lateralized 130°

Compatibility Matrix (Intended for countries outside the EU)

Echo BiMetric	192807_8	Echo Bi-Metric Microplasty Reduced Proximal Profile - Standard 135°
Echo BiMetric	192907_8	Echo Bi-Metric Microplasty Reduced Proximal Profile - High Offset 130°
Echo BiMetric	1930XX (07_21)	Echo Bi-Metric Microplasty Full Proximal Profile - Standard 135°
Echo BiMetric	1931XX (07_21)	Echo Bi-Metric Microplasty Full Proximal Profile -High Offset 130°
Echo FX		Echo Hip Fracture Cemented Stem Type 1
Echo FX	12-1513XX	Echo FX Cobalt Chromium Femoral
Echo FX	12-1514XX	Echo FX Cobalt Chromium Femoral, Lateralized
OSS Low Profile		OSS Low Profile ProXimal Stem
OSS Low Profile	151820_21	OSS Proximal Femorals
OSS Elliptical		OSS Elliptical Stem
OSS Elliptical	150497_8	OSS Proximal Femorals
OSS Modular		OSS Modular Stem
OSS Modular	150457_60	OSS Proximal Femorals
Taperloc Legacy Uncemented		Taperloc Legacy Primary Uncemented Stems Type 1
Taperloc Legacy Uncemented	16440X	Taperloc Primary Stem Porous Ti Type 1 Taper Cementless
Taperloc Legacy Cemented		Taperloc Legacy Primary Cemented Stems
Taperloc Legacy Cemented	650-0325_30	Taperloc Cemented Primary Stem - Standard - Type 1
Taperloc Legacy Cemented	650-0331_36	Taperloc Cemented Primary Stem - Lateralised - Type 1
Taperloc Complete BM		Taperloc Complete Bonemaster Type 1 Stems
Taperloc Complete BM	51-111XX0 (05_18,20,22,24)	Taperloc Complete Bm Type 1
Taperloc Complete BM	51-113XX0 (09_18,20,22,24)	Taperloc Complete Bm Type 1
Taperloc Complete BM	51-114XX0 (09_18,20,22,24)	Taperloc Complete Bm Type 1
Taperloc Complete BM	51-117XX0 (09_18,20,22,24)	Tprlc Bm High O/S Mp Type 1
Taperloc Complete BM	51-116XX0 (09_18,20,22,24)	Tprlc Bm Std O/S Mp Type 1
Taperloc Complete BM	51-1180X0 (4_8)	Tprlc Bm Std O/S Mp Type 1
Taperloc Complete BM	51-1190X0 (5_8)	Tprlc Bm High O/S Mp Type 1
Taperloc Complete PPS		Taperloc Complete PPS Full Length Type 1 Stem
Taperloc Complete PPS	51-1020X0 (4_8)	Taperloc FP Type 1 PPS
Taperloc Complete PPS	51-101XX0 (05_20,22,24)	Taperloc FP Type 1 PPS HO
Taperloc Complete PPS	51-105XX0 (09_19)	Taperloc Type 1 PPS
Taperloc Complete PPS	51-100XX0 (04_20,22,24)	Taperloc Complete PPS Type 1 Stem
Taperloc Complete PPS	51-103XX0 (09_20,22,24)	Taperloc Complete PPS Type 1 Stem
Taperloc Complete PPS	51-104XX0 (09_20,22,24)	Taperloc Complete PPS Type 1 Stem
Taperloc Complete PPS	51-106XX0 (09_20,22,24)	Taperloc Complete PPS Std Offset
Taperloc Complete PPS	51-107XX0 (09_20,22,24)	Taperloc Complete PPS Hi Offset
Taperloc Complete PPS	51-1080X0 (4_8)	Taperloc Complete PPS Std Offset
Taperloc Complete PPS	51-110XX0 (04_18,20,22,24)	Taperloc Complete PPS Hi Offset
Taperloc Complete PPS	51-1090X0 (5_8)	Taperloc Complete PPS Std Offset

Compatibility Matrix (Intended for countries outside the EU)

Taperloc Complete PPS	51-1490X0 (4_8)	Taperloc Complete PPS Coxa Vara
Taperloc Complete Microplasty		Taperloc Complete Microplasty Stem
Taperloc Complete Microplasty	51-145XX0 (09_20,22,24)	Micro Taperloc Complete PPS Type1
Taperloc Legacy Lateralized		Taperloc Legacy Full Distal Lateralized Offset
Taperloc Legacy Lateralized	11-1032XX(11-103200_216)	Taperloc POR FMRL
Taperloc Legacy Standard		Taperloc Legacy Full Distal Standard Offset
Taperloc Legacy Standard	1032XX (103200_216)	Taperloc POR FMRL

Type 1 Heads		
Group ID	Item Number	Description
Freedom Type 1	110025130_134	FREEDOM CONSTR HD, D32 with neck length -6/-3/0/+3/+6
Freedom Type 1	11-107016_107021	FREEDOM CONSTR HD, D36 with neck length -6/-3/0/+3/+6/+9
Biomet CoCr Type 1	010001031_44	CoCr Head Type 1, D40/44 with neck length -6/-3/0/+3/+6/+9/+12
BioloX Delta	650-1157_1159	BioloX Delta Ceramic Fem Head, D28 with neck length -3/0/+3
BioloX Delta	650-1160_1163	BioloX Delta Ceramic Fem Head, D32 with neck length -3/0/+3/+6
BioloX Delta	650-0660_0663	BioloX Delta Ceramic Fem Head, D36 with neck length -3/0/+3/+6
ZB Ceramic T1, 22.2*	802802201_03*	ZB Ceramic Femoral Heads*, D22.2 with neck length -3/0/+3*
ZB Ceramic T1	802802801_03	ZB Ceramic Femoral Heads, D28 with neck length -3/0/+3
ZB Ceramic T1	802803201_03	ZB Ceramic Femoral Heads, D32 with neck length -3/0/+3
ZB Ceramic T1	802803601_4404	ZB Ceramic Femoral Heads, D36/40/44 with neck length -3/0/+3/+6
ZB Select Ceramic	803003201_803004401	ZB Select Ceramic Head Implants D32/36/40/44
Ringloc Heads	163627	COCR MOD HD NO SKIRT, D26 with neck length +6
Ringloc Heads	163638	COCR MOD HD NO SKIRT, D28 with neck length +6
Ringloc Heads	163674	COCR MOD HD NO SKIRT, D32 with neck length +6
Ringloc Heads	163651_53	COCR MOD HD STD NK, D22.2 with neck length -5/-3/0
Ringloc Heads	163620_26	DIA COCR MOD HD, D26 with neck length -6/-3/0/+3/+6/+9/+12
Ringloc Heads	163660_666	DIA COCR MOD HD, D28 with neck length -6/-3/0/+3/+6/+9/+12
Ringloc Heads	163667_673	DIA COCR MOD HD, D32 with neck length -6/-3/0/+3/+6/+9/+12
Ringloc Heads	11-363660_666	DIA COCR MOD HD, D36 with neck length -6/-3/0/+3/+6/+9/+12
Endo II**	12-139006_46**	Endo II CoCr Head**
BioloX Delta Option	650-1055_1059	BioloX Delta Option Head, D28/32/36/40

Type 1 Adaptors		
Group ID	Item Number	Description
ZB Select Adaptor T1	803000101_105	ZB Type 1 Select Adaptor
BioloX Option Adaptor	650-1064_1068	BIOLOX Option Adaptor (Type 1)
Endo II Adaptor	139245_139249	Endo II Taper Insert

Abbreviations used in this document:

CoCr = Cobalt Chromium

Notes, references:

*Size 22.2mm heads are contraindicated for use with Stainless Steel and Cobalt Chromium stems

** Endo II CoCr Head must be used with Endo II Taper Insert for Type 1 Stems

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for health care professionals and Zimmer Biomet Sales Staff. Distribution to any other recipient is prohibited. For product information, including indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the Instructions for Use (IFU) and www.zimmerbiomet.com.

An electronic copy of the Instructions for Use (eIFU) and this Product Compatibility Matrix can be obtained from the Zimmer Biomet Electronic Labeling Service (E-labeling) website: <https://labeling.zimmerbiomet.com>

Check for country product clearances and reference product specific Instructions for Use.

© 2025 Zimmer Biomet



CM-0001-01-EN Rev. 1 Issue Date 2025-10-16