

SCOPE

12/14 Femoral Stems		
Group ID	Item Number	Description
ALC SL		Alloclassic® Zweymüller SL Stem
ALC SL	28XX (2839-2852)	Alloclassic SL
ALC SLO		Alloclassic® Zweymüller SL Offset Stem
ALC SLO	01.00121.XXX	Alloclassic SL Offset
ALC SLL		Alloclassic® Zweymüller SLL Revision Hip Stem
ALC SLL	28XX (2883-2891)	Alloclassic SLL
Avenir Cemented		Avenir Müller Cemented Hip Stems-
Avenir Cemented	01.06010.20X	Avenir Müller Cemented Hip Stem Standard
Avenir Cemented	01.06010.30X	Avenir Müller Cemented Hip Stem Lateral
Avenir Complete, Uncem.		Avenir Complete® Stems
Avenir Complete, Uncem.	574101XXX(000-090)	Avenir Complete HA Standard Collarless
Avenir Complete, Uncem.	574102XXX(000-090)	Avenir Complete HA High Offset Collarless
Avenir Complete, Uncem.	574103XXX(000-090)	Avenir Complete HA Coxa Vara Collarless
Avenir Complete, Uncem.	574201XXX(000-090)	Avenir Complete HA Standard Collared
Avenir Complete, Uncem.	574202XXX(000-090)	Avenir Complete HA High Offset Collared
Avenir Complete, Uncem.	574203XXX(000-090)	Avenir Complete HA Coxa Vara Collared
Avenir Uncemented		Avenir Müller Uncemented Hip Stem
Avenir Uncemented	01.06010.00X	Avenir Müller Uncemented Hip Stem Standard
Avenir Uncemented	01.06010.10X	Avenir Müller Uncemented Hip Stem Lateral
CLS		CLS Spotorno Hip Stem
CLS	01.00295.0XX	CLS Spotorno 125°
CLS	29.00.09-XXX	CLS Spotorno 145°
CLS	29.00.39-XXX	CLS Spotorno 135°
CMK		CMK Primary Stems
CMK	3553.XXX(101-421)	CMK OC Stems - Stainless Steel 12/14 Taper
CMK	3553.XXX(501-521)	CMK OC Stems - Stainless Steel 12/14 Taper
CMK	3553.291_492-200_250	CMK OC Stems - Stainless Steel 12/14 Taper - Revision
CMK	3553.591-200_250	CMK OC Stems - Stainless Steel 12/14 Taper - Revision
CPT		CPT® 12/14 CoCr stem
CPT	00-8114-00X-00	CPT Standart Offset
CPT	00-8114-0X0-00 (4_5)	CPT Small
CPT	00-8114-00X-18 (2_3)	CPT Small
CPT	00-8114-00X-10 (0_5)	CPT Extended Offset
CPT	00-8114-01X-18 (2_3)	CPT Revision Long
CPT	00-8114-004-20	CPT Revision Long
CPT	00-8114-004-23	CPT Revision Long
CPT	00-8114-004-26	CPT Revision Long
CPT	00-8114-00X-30 (2_5)	CPT Extra Extended Offset
Exception, Cemented		EXception® Primary Cemented Stems (SS)

Compatibility Matrix (Intended for countries outside the EU)



Exception, Cemented	PV1250XX	Exception Vara Cemented Stem - Left
Exception, Cemented	PV1260XX	Exception Vara Cemented Stem - Right
Exception, Uncemented		EXception® Primary Uncemented Stems
Exception, Uncemented	PV125YXX	Exception Vara Stem - Cementless - Left
Exception, Uncemented	PV126YXX	Exception Vara Stem - Cementless - Right
Fitmore		Fitmore® Hip Stem
Fitmore	01.00551.1XX	Fitmore Family A
Fitmore	01.00551.2XX	Fitmore Family B
Fitmore	01.00551.3XX	Fitmore Family B, Extended Offset
Fitmore	01.00551.4XX	Fitmore Family C
M/L Taper		M/L Taper stem
M/L Taper	00-7711-0XX-00	ML Taper Standard
M/L Taper	65-7711-0XX-00	ML Taper Standard, HA/TCP
M/L Taper	00-7711-0XX-10	ML Taper, Reduced Neck Length
M/L Taper	65-7711-004-10	ML Taper, Reduced Neck Length, HA/TCP
M/L Taper	00-7711-0XX-20	ML Taper Extended offset
M/L Taper	65-7711-0XX-20	ML Taper Extended offset, HA/TCP
M/L Taper	00-7711-0XX-40	ML Taper, Extended Offset, Reduced Neck Length
M/L Taper	65-7711-004-40	ML Taper, Extended Offset, Reduced Neck Length, HA/TCP
M/L Taper, Kin.		M/L Taper stem with Kinectiv® Technology
M/L Taper, Kin.	00-7713-0XX-00	ML Taper Kinectiv
M/L Taper, Kin.	65-7713-0XX-00	ML Taper Kinectiv, HA/TCP
M/L Taper, Kin.	00-7848-0XX-00	Kinectiv Modular Necks
MS-30		MS-30® Stem
MS-30	01.00351.XXX	MS-30 Stem Lateral
MS-30	30.00.49-XXX	MS-30 Stem Standard
Müller Straight (CoCr)		Original Müller Straight stem PROTASUL-10
Müller Straight (CoCr)	12.00.29.XXX	Original Mueller Straight Stem Standard Protasul 10
Müller Straight (CoCr)	12.00.39.XXX	Original Mueller Straight Stem Lateral Protasul 10
Müller Straight (CoCr)	05.95001.0XX	Additional 8.5 Size in Protasul 10 Only
Müller Straight (SS)		Original Müller Straight stem PROTASUL-S30
Müller Straight (SS)	35.00.29.XXX	Original Mueller Straight Stem Standard Protasul 30
Müller Straight (SS)	35.00.39.XXX	Original Mueller Straight Stem Lateral Protasul 30
Müller Straight (SS)	05.95001.066	Mueller Straight Stem Lateral
Müller Straight (SS)	05.95001.067	Mueller Straight Stem Lateral
Revitan Revision		Revitan® Revision Hip Stem (straight/curved)
Revitan Revision	01.00401.0XX	Revitan spout body
Revitan Revision	01.00402.XXX	Revitan Cylindrical Body
Revitan Revision	01.00405.XXX	Revitan Straight
Revitan Revision	01.00406.XXX	Revitan Curved
Segmental		Segmental ProXimal Femur
Segmental	00-5850-0XX-XX	Zimmer Segmental Proximal Femoral Body
Taperloc L.+C., MP		Taperloc® Complete and Legacy Microplasty Stems
Taperloc L.+C., MP	51-126090	Taperloc Complete Micro RD 12/14

Compatibility Matrix (Intended for countries outside the EU)



Taperloc L.+C., MP	51-1261X0(0_9)	Taperloc Complete Micro RD 12/14
Taperloc L.+C., MP	51-1262X0(0,2,4)	Taperloc Complete Micro RD 12/14
Taperloc L.+C., MP	51-127090	Taperloc Complete Micro RD Ho TI 12/14
Taperloc L.+C., MP	51-1271X0(0_9)	Taperloc Complete Micro RD Ho TI 12/14
Taperloc L.+C., MP	51-1272X0(0,2,4)	Taperloc Complete Micro RD Ho TI 12/14
Taperloc L.+C., MP	51-1280XX(40_70,80)	Taperloc Complete Micro TI PC 12/14
Taperloc L.+C., MP	51-1290XX(50_70,80)	Taperloc Complete Micro HO TI PC 12/14
Taperloc L.+C., MP	650-09XX	Taperloc MicroPlasty Standard PC 12/14 Taper
Taperloc L.+C., Uncem.		Taperloc® Complete and Legacy Uncemented Stems
Taperloc L.+C., Uncem.	51-120XX0(06_24)	Taperloc Complete Standard Offset 12/14 Taper - PPS
Taperloc L.+C., Uncem.	51-121XX0(05_24)	Taperloc Complete High Offset 12/14 Taper - PPS
Taperloc L.+C., Uncem.	51-123090, 1X0 (0_9)	Taperloc Complete 12/14 PPS SO
Taperloc L.+C., Uncem.	51-1232X0 (0,2,4)	Taperloc Complete 12/14 PPS SO
Taperloc L.+C., Uncem.	51-124XX0 (09_18)	Taperloc Complete 12/14 PPS HO
Taperloc L.+C., Uncem.	650-0580_90	Taperloc PPS 12/14
Taperloc L.+C., Uncem.	650-026X	Taperloc PPS 12/14
Taperloc L.+C., Uncem.	650-03XX	Taperloc PPS 12/14
Taperloc L., Cem.		Taperloc® Legacy Cemented Stems
Taperloc L., Cem.	650-0337_42	Taperloc Cemented Primary Stem - Standard - 12/14
Taperloc L., Cem.	650-0343_48	Taperloc Cemented Primary Stem - Lateralised - 12/14
Taperloc C., BM, RD		Taperloc® Complete Bonemaster™ Reduced Distal 12/14
Taperloc C., BM, RD	51-133XX0(09_80)	Taperloc Complete Reduced Distal Std Offset - Bonemaster
Taperloc C., BM, RD	51-1332X0(0,2,4)	Taperloc Complete Reduced Distal Std Offset - Bonemaster
Taperloc C., BM, RD	51-134XX0(09_18,20,40)	Taperloc Complete Reduced Distal Ext Offset - Bonemaster
Taperloc C., BM, RD	51-136XX0(09_18,20,40)	Taperloc Complete Reduced Distal Std Offset MP - Bonemaster HA
Taperloc C., BM, RD	51-137XX0(09_18,20,40)	Taperloc Complete Reduced Distal High Offset MP - Bonemaster HA
Taperloc C., BM, MP		Taperloc® Complete Bonemaster™ Microplasty Stem
Taperloc C., BM, MP	51-1380X0(4_7,8)	Taperloc Complete Std Offset MP - Bonemaster HA
Taperloc C., BM, MP	51-1390X0(5_7,8)	Taperloc Complete High Offset MP - Bonemaster HA
Taperloc C., BM		Taperloc® Complete Bonemaster™ Stems
Taperloc C., BM	51-130XX0 (04_18,20,22,24)	Taperloc Complete Standard Offset 12/14 Taper - Bonemaster HA
Taperloc C., BM	51-131XX0 (04_18,20,22,24)	Taperloc Complete High Offset 12/14 Taper - Bonemaster HA
Taperloc C., XR		Taperloc® Complete XR 123 Stems
Taperloc C., XR	51-1220X0 (4_8)	Taperloc Complete 12/14 PPS XR 123
Taperloc C., XR	51-125XX0 (09_17)	Taperloc Complete 12/14 PPS XR 123
TM Primary		Trabecular Metal® Primary Hip Stem
TM Primary	00-7864-0XX-00	TM Primary
TM Primary	00-7864-0XX-20	TM Primary Extended Offset
Versys Advocate		Versys® Advocate
Versys Advocate	00-7850-0XX-00	Versys Advocate Standard Offset-Non-VLign
Versys Advocate	00-7850-0XX-20	Versys Advocate Standard Offset-Non-Vlign, Extended offset
VerSys CRC		VerSys® Cemented Revision Calcar (CRC)
VerSys CRC	00-7871-0XX-XX	VerSys Cemented Revision Calcar
VerSys Heritage		VerSys Heritage®

Compatibility Matrix (Intended for countries outside the EU)



VerSys Heritage	00-7857-0XX-00 (09_17)	Versys-Heritage CDH
VerSys Heritage	00-7857-0XX-20	Versys Heritage Primary Extended Offset
VerSys FMT		Versys® Fiber Metal Taper
VerSys FMT	00-7862-0XX-00	Versys Fiber Metal Taper
VerSys FMT	65-7862-0XX-00	Versys Fiber Metal Taper, HA/TCP
VerSys FMT	00-7862-0XX-20	Versys Fiber Metal Taper, Extended Offset
VerSys FMT	65-7862-0XX-20	Versys Fiber Metal Taper, Extended Offse, HA/TCP
VerSys FMT	00-7862-0XX-30	Versys Fiber Metal Taper, LM
VerSys FMT	65-7862-0XX-30	Versys Fiber Metal Taper, LM, HA/TCP
VerSys FMT	00-7862-0XX-50	Versys Fiber Metal Taper, LM, Extended Offset
VerSys FMT	65-7862-0XX-50	Versys Fiber Metal Taper, LM, Extended Offset, HA/TCP
VerSys FMT MC HATCP		Versys® Fiber Met MidCoat HATCP
VerSys FMT MC HATCP	65-7641-0XX-00	VerSys Fiber Metal Midcoat CL, HA
VerSys FMT MC HATCP	65-7641-0XX-30	VerSys Fiber Metal Midcoat CL, LM, HA
VerSys FMT MC HATCP	65-7645-0XX-00	VerSys Fiber Metal Midcoat, HA
VerSys FMT MC HATCP	65-7645-0XX-30	VerSys Fiber Metal Midcoat, LM HA
VerSys FMT HATCP		Versys® FM Taper HATCP
VerSys FMT HATCP	65-7662-0XX-00	VerSys Fiber Metal Taper
VerSys FMT HATCP	65-7662-0XX-30	VerSys Fiber Metal Taper, LM
Wagner Cone		Wagner Cone Prosthesis® Hip Stem
Wagner Cone	01.00561.2XX	Wagner Cone Prosthesis 125°
Wagner Cone	01.00561.3XX	Wagner Cone Prosthesis 135°
Wagner SL Revision		Wagner SL Revision® Hip Stem
Wagner SL Revision	01.00101.9XX	Wagner SL Revision Stem, 190mm
Wagner SL Revision	01.00102.2XX	Wagner SL Revision Stem, 225mm
Wagner SL Revision	01.00102.6XX	Wagner SL Revision Stem, 265mm
Wagner SL Revision	01.00103.0XX	Wagner SL Revision Stem, 305mm
Z1, Collarless		Z1® Collarless Cementless
Z1, Collarless	611101XXX	Standard Offset
Z1, Collarless	611102XXX	High Offset
Z1, Collarless	611103XXX	Coxa Vara
Z1, Collared		Z1® Collared Cementless
Z1, Collared	611201XXX	Standard
Z1, Collared	611202XXX	High Offset
Z1, Collared	611203XXX	Coxa Vara

12/14 Heads		
Group ID	Item Number	Description
G7 CoCr Freedom Heads, Biomet	110025126_129	G7® CoCr Freedom Head (Biomet 12/14), D32 with neck length -3/0/+3/+6
BioloX d., Biomet	650-0667	BioloX Delta Head Biomet 12/14, D36 with neck length +8
BioloX d., Biomet	650-0836_0838	BioloX Delta Head Biomet 12/14, D36 with neck length -4/0/+4
BioloX d., Biomet	650-0830_0832	BioloX Delta Head Biomet 12/14, D28 with neck length -3.5/0/+3.5
BioloX d., Biomet	650-0833_0835	BioloX Delta Head Biomet 12/14, D32 with neck length -4/0/+4
BioloX d., Zimmer	00-8775-028-01_03	BioloX Delta Heads Zimmer 12/14, D28 with neck length -3.5/0/+3.5

Compatibility Matrix (Intended for countries outside the EU)



Bilox d., Zimmer	00-8775-032_040-01_03	Bilox Delta Heads Zimmer 12/14, D28/32/36 with neck length -3.5/0/+3.5
Bilox d., Zimmer	00-8775-032_040-04	Bilox Delta Heads Zimmer 12/14, D28/32/36 with neck length +7
Bilox Opt., Zimmer (Head+Adapter)	00-8777-028-01_04	Bilox Option Head Zimmer 12/14, D28 with neck length -3.0/0/+3.5 (Head and Adaptor, Packed Together)
Bilox Opt., Zimmer (Head+Adapter)	00-8777-032_040-01_04	Bilox Option Head Zimmer 12/14, D32/36/40 with neck length -3.0/0/+3.5/+7.0 (Head and Adaptor, Packed Together)
Bilox Opt., Zimmer (Head+Adapter)	00-8777-044-01_04	Bilox Option Head Zimmer 12/14, D44 with neck length -3.0/0/+3.5/+7.0 (Head and Adaptor, Packed Together)
Bilox Opt., Zimmer (Head+Adapter)	00-8777-048-01_04	Bilox Option Head Zimmer 12/14, D48 with neck length -3.0/0/+3.5/+7.0 (Head and Adaptor, Packed Together)
Bilox Opt., Biomet (Head only)	650-1055_1059	Bilox Delta Option Biomet 12/14, D28/32/36/40/44 (Head only)
Durasul CoCr Head	01.01012.384_388	Durasul CoCr Ball Head, D38 with neck length -8/-4/0/+4/+8
ZB 12/14 CoCr	802202201_03	ZB 12/14 CoCr Femoral Heads Universal, D22.2 with neck length -2/0/+3
ZB 12/14 CoCr	802202801_04	ZB 12/14 CoCr Femoral Heads Universal, D28 with neck length -3.5/0/+3.5/+7
ZB 12/14 CoCr, 10.5	802202805	ZB 12/14 CoCr Femoral Heads Universal, D28 with neck length +10.5
ZB 12/14 CoCr	8022032_4401_04	ZB 12/14 CoCr Femoral Heads Universal, D32/36/40/44 with neck length -3.5/0/+3.5/+7
ZB 12/14 CoCr, 10.5	8022032_4405	ZB 12/14 CoCr Femoral Heads Universal, D32/36/40/44 with neck length +10.5
ZB 12/14 Ceramic Head*	802602202_03*	ZB 12/14 Ceramic Femoral Heads, D22 with neck length 0/+3*
ZB 12/14 Ceramic Head	802602801_03	ZB 12/14 Ceramic Femoral Heads, D28 with neck length -3.5/0/+3.5
ZB 12/14 Ceramic Head	802604001_4404	ZB 12/14 Ceramic Femoral Heads, D40/44 with neck length -3.5/0/+3.5/+7
ZB 12/14 CoCr Freedom Head	802403201_802403203	ZB 12/14 CoCr Freedom Heads, D32 with neck length -3/0/+3//+6
ZB 12/14 CoCr Freedom Head	802403601_802403605	ZB 12/14 CoCr Freedom Heads, D36 with neck length -3/0/+3/+6
ZB 12/14 CoCr Freedom Head, +9	802403606	ZB 12/14 CoCr Freedom Heads, D36 with neck length +9
Sulox Head	17.28.05_07	Sulox Head, D28 with neck length -3.5/0/+3.5
Sulox Head	17.32.05_07	Sulox Head, D32 with neck length -3.5/0/+3.5
Unipolar Head	30.38_60.09	Unipolar Head, D38/40/42/44/46/48/50/52/54/56/58/60 with neck length -4.0
ZB 12/14 Select (Head only)	8030028_4401	ZB Ceramic Select Head (Head Only)
ZB 12/14 Select (Head+Adapter)	8030032_4400	ZB Ceramic Select Head (Head and Adaptor, Packed Together)
Versys Endo (Head only)	00-7818-0xx-00	Versys Endo Femoral Heads
Protasul-S30	30.22.05_07	Protasul-S30 Head, D22 with neck length -3.5/0/+3.5
Protasul-S30	30.28.05_07	Protasul-S30 Head, D28 with neck length -4.0/0/+4.0
Protasul-S30	30.32.05_07	Protasul-S30 Head, D32 with neck length -4.0/0/+4.0

Adaptors		
Group ID	Item Number	Description
Bilox Option, Biomet Adaptor**	650-1060_1063	Bilox Option (Biomet 12/14 Taper Adaptor)
ZB 12/14 Select Adaptor***	803000000_04	ZB 12/14 Select Adaptor
Versys Endo Adaptor ⁽⁵⁾ , +7	00-7818-099-00	Versys Endo Femoral Head Adaptor +7.0
Versys Endo Adaptor ⁽⁵⁾ , +10.5/+14.0	00-7818-099-01_2	Versys Endo Femoral Head Adaptor +10.5/14.0

Abbreviations used in this document:

ZB - Zimmer Biomet

TM - Trabecular Metal

SS - Stainless Steel

CoCr - Cobalt-chrome

Notes, references:

*Size 22.2 mm heads are contraindicated for use with Stainless Steel and Cobalt Chromium stems

**Only for Biomet Biolox Option Heads in group "Biolox Opt., Biomet (Head only)"

*** Only for Zimmer Biomet Select Ceramic Head in group "ZB 12/14 Select (Head only)"

⁽⁴⁾Revitan Straight Revision Stem Ø14 and Wagner SL Revision Stems Ø14-15 mm cannot be used with head neck lengths longer than +4 mm

⁽⁵⁾Only for Versys Endofemoral Heads in group "Versys Endo (Head only)"

This material is intended for health care professionals and the Zimmer Biomet sales force. Distribution to any other recipient is prohibited. For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the instructions for use or contact your local representative; visit www.zimmerbiomet.com for additional product information. All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet. An electronic copy of the Instructions for Use (eIFU) and this Product Compatibility Matrix can be obtained from the ZimmerBiomet Electronic Labeling Service (E-labeling) website: <https://labeling.zimmerbiomet.com>

BIOLOX® Ceramic is a trademark of CeramTec GmbH.

Check for country product clearances and reference product specific Instructions for Use.

© 2025, 2026 Zimmer Biomet



CM-0001-02-EN Issue Date 2026-02-27