



# Bonus<sup>®</sup> CC Matrix

## Bone Graft System

An adaptable graft hydration system

### A formulation of demineralized cortical bone and mineralized cancellous chips.

This bone graft system offers all the benefits of an allograft while providing the surgeon a consistent and efficient method of hydration, including optional use of outputs from Zimmer Biomet Biologics autologous products as hydration media.

#### Bone Remodeling

Osteoconductive cancellous bone

Osteoinductivity-verified Demineralized Bone Matrix (DBM) powder

Hydrate with autologous blood products, bone marrow aspirate, saline, or antibiotic solution

#### Handling Characteristics

User controlled graft consistency

Hydration flexibility

Blend of cortical and cancellous bone matrices

#### Convenience

All inclusive bone grafting kit



Available in 5cc and 10cc sizes

Ease of preparation



**ZIMMER BIOMET**  
Your progress. Our promise.®



## Graft Preparation

Production	Description	Part Number
	Bonus CC Matrix Graft System 5cc	48-1805
	Bonus CC Matrix Graft System 10cc	48-1810

## Autologous Product Offerings

Production	Description	Part Number
	BioCUE® Mini Concentration System	800-0610A
	GPS® III Mini Kit	800-0505A
	Clotalyst® Kit with GPS III Mini Separator	800-0724
	Plasmax® Mini Kit	800-0516

## Accessories\*

Production	Description	Part Number
	Dual Ratio Applicator Kit	800-0251
	Mixing Connector with Spray Tip	800-0252



All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for health care professionals. Distribution to any other recipient is prohibited.

For product information, including indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert and [www.zimmerbiomet.com](http://www.zimmerbiomet.com).

© 2018 Zimmer Biomet



**Legal Manufacturer**  
 Biomet Biologics  
 P.O. Box 587  
 56 E. Bell Drive  
 Warsaw, Indiana 46581-0587  
 USA

\* **Legal Manufacturer**  
 Micromedics, Inc.  
 1270 Eagan Industrial Road  
 St. Paul, MN 55121-1385  
 USA