

iQ™ Intelligent System



BIOMET®
MICROFIXATION
One Surgeon. One Patient.™

A revolutionary device that combines high-speed drilling with cordless convenience

- Bicortical drilling capability at 10,000 RPM
- Eliminates cords, foot pedals and equipment setups
- Activates with a simple press of a button

Midface Drills

Part #	Description
72-2030	1.1mm drill w/ 5mm stop
72-2034	1.1mm drill w/ 7mm stop
72-2038	1.1mm drill w/ 9mm stop
72-2042	1.1mm drill w/ 15mm stop

Mandible Drills

Part #	Description
72-2049	1.5mm drill w/ 5mm stop
72-2052	1.5mm drill w/ 7mm stop
72-2055	1.5mm drill w/ 22mm stop
72-2057	1.6mm drill w/ 26mm stop
72-2060	1.8mm drill w/ 26mm stop

Percutaneous Drills*

Part #	Description
72-2019	1.5mm drill w/ 22mm stop
72-2021	1.6mm drill w/ 26mm stop**
72-2023	1.8mm drill w/ 26mm stop**

* These drills are intended to be used with a trocar.

**The drills have colored bands to determine the bicortical screw length needed, eliminating the need for a depth gauge (color/screw length legend shown below).

BLACK	YELLOW	PURPLE	GREEN	BROWN	RED	BLUE
6MM	8MM	10MM	12MM	14MM	16MM	18MM



Experience convenient drilling with our comprehensive plating systems.

Biomet Microfixation provides a comprehensive array of craniomaxillofacial products to meet a broad range of your procedure needs.

