



## Polaris™ 5.5 Spinal System


We've Got the Lock on  
Pedicle Screw Technology

- 5.5mm diameter titanium rod system
- Helical Flange® enhanced locking technology
- Multi-axial and reduction screws available

**BIOMET**  
SPINE

## Polaris™ 5.5 Spinal System

### Helical Flange® Technology



- Starts easily
- Minimizes cross threading and seat splay
- Forces are concentrated inward

### Multi-axial Screw



- Friction fit screw seat stays in place once aligned in desired position; easy rod placement
- Anatomical fit minimizes adjacent segment impingement
- Color-coded implants for easy determination of screw sizes
- Site hole on seat for simple rod manipulation and excellent interface with instruments

### Screw Insertion



- Screw driver that is efficient and easy to use
- Allows excellent visualization
- Control of screw trajectory

### Rod Reduction



- Very powerful and easily controlled
- User-friendly
- Locks onto seat for easy rod reduction

**BIOMET**  
SPINE

Broomfield, CO • 800.447.3625  
biometspine.com • BSP216508L 01/15

©2015 Biomet Spine, LLC. All rights reserved. All trademarks are the property of Biomet, Inc. or one of its subsidiaries unless otherwise indicated. Helical Flange® is a registered trademark of Roger P. Jackson. Rx Only.