


Reference Combination Manual/Automated Cleaning and Disinfection Instructions (Section H) of the Instructions for Care, Cleaning, Maintenance and Sterilization Manual (97-5000-170-00). Do not utilize the optional ultrasonic cleaning note in Section H.



**Persona® Shims**  
42-5279-000-01\_04  
42-5279-002-01\_04  
42-5279-004-01\_04  
42-5279-006-01\_04  
42-5279-008-01-04  
Cleaning Instructions

Do not ultrasonic clean



97-5026-055-00 Rev. 1 Printed in USA ©2015 Zimmer, Inc.