

Head and Stem Combinations: Ceramic Femoral Heads

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	Alumina 12/14 Femoral Head (XX-6418-026-XX)	Alumina 12/14 Femoral Head (XX-6418/6428-028-XX)	Alumina 12/14 Femoral Head (XX-6418/6428-032-XX)	BIOLOX-forte head (12.2832.05/07) (12/14 taper)	CERASUL head (12/14 taper)	CERASUL head (8/10 taper)	Continuum Alumina Bearing Head (XX-4102-XXX-XX) (12/14 taper)	Continuum Zirconia Bearing Head (XX-4101-XXX-XX) (12/14 taper)	ESKA BIOLOX Forte head (12/14 taper)	ESKA BIOLOX Forte head (8/10 taper)	NGK Zirconia femoral head (6 degree)	SULOX head (12/14 taper)	SULOX Head (8/10 taper)	Tete Alumine Cone 8/10 D28	Zirconia head (7676-XX-XXX) (12/14 taper)	BIOLOX delta head (00-8775-028/040-xx)(12/14 taper)	BIOLOX OPTION head system (00-8777-028/040-xx)(12/14 taper)	BIOLOX OPTION head system (00-8778-028/040-xx)(8/10 taper)
ACA stem	✓	✓	✓	✓								✓				✓	✓	✓
Acora cemented stem						✓							✓	✓				✓
ALLOCLASSIC SLA and SL-64 stem	✓	✓	✓	✓								✓			✓	✓	✓	✓
ALLOCLASSIC SLL-Revision stem	✓	✓	✓	✓								✓			✓	✓	✓	✓
ALLOCLASSIC stem	✓	✓	✓	✓								✓			✓	✓	✓	✓
ALLOCLASSIC VARIALL SLV stem	✓	✓	✓	✓								✓			✓	✓	✓	✓
ALLOCOR AN stem	✓	✓	✓	✓								✓			✓	✓	✓	✓
Anatomic stem											✓							
Anatomic II stem	✓	✓	✓															
Anatomique uncemented HA stem				✓	✓							✓				✓	✓	✓
Apollo Hip System				✓	✓							✓			✓	✓	✓	✓
APR Hip System, Non-porous				✓	✓							✓			✓	✓	✓	✓
APR Hip System, Textured		✓	✓	✓	✓							✓			✓	✓	✓	✓
APS Natural-Hip System	✓	✓	✓	✓	✓							✓				✓	✓	✓
Avenir Müller stem		✓	✓	✓	✓							✓				✓	✓	✓
Avenir Cemented Stem		✓	✓	✓	✓							✓				✓	✓	✓
Avenir Complete stem		✓	✓	✓	✓							✓				✓	✓	✓
BEO stem		✓	✓	✓	✓							✓				✓	✓	✓
CDH stem PROTASUL-S30 (monoblock)																		
CLS Brevius Stem with Kinectiv Technology																		
CLS SPOTORNO stem	✓	✓	✓	✓								✓			✓	✓	✓	✓
Concept cemented stem						✓							✓	✓				✓
Continuum Cobrex stem							✓	✓										
Continuum F-115 Hip stem							✓	✓										
Continuum F-117 Collarless Hip stem							✓	✓										
Continuum J-Taper stem																		
Continuum Proxilock stem (HA coated)		✓					✓	✓										
Continuum Proxilock stem (porous coated)							✓	✓										
CPS uncemented stem		✓	✓	✓	✓							✓				✓	✓	✓
CPT (6 Deg, CoCr) stem											✓							
CPT (6 Deg, SST) stem											✓							
CPT 12/14 CoCr stem	✓	✓	✓	✓	✓							✓				✓	✓	✓
DCM-J Natural-Hip System	✓	✓	✓	✓	✓							✓				✓	✓	✓
DSP / TPP / PDC		✓	✓	✓	✓							✓				✓	✓	✓
Emeraude stem		✓	✓	✓	✓							✓				✓	✓	✓
Epoch stem	✓	✓	✓	✓	✓							✓				✓	✓	✓
ESKA 8/10 stem										✓								
ESKA 12/14 stem								✓										
Etoile cemented stem						✓							✓	✓				✓
EXAFIT Primary stem						✓							✓	✓				✓
EXAFIT Revision stem						✓							✓	✓				✓
Fitmore Hip Stem	✓	✓	✓	✓	✓							✓				✓	✓	✓
FracSure Hip system																		
Heritage stem (XX-8095-XXX-XX)											✓							
Intermedia stem				✓	✓							✓				✓	✓	✓
Lamella stem	✓	✓	✓	✓	✓							✓				✓	✓	✓

Refer to Kinectiv Technology Page

Continued on Page 2

Head and Stem Combinations: Ceramic Femoral Heads

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	Alumina 12/14 Femoral Head (XX-6418-026-XX)	Alumina 12/14 Femoral Head (XX-6418/6428-028-XX)	Alumina 12/14 Femoral Head (XX-6418/6428-032-XX)	BIOLOX-forte head (12.2832.05/07) (12/14 taper)	CERASUL head (12/14 taper)	CERASUL head (8/10 taper)	Continuum Alumina Bearing Head (XX-4102-XXX-XX) (12/14 taper)	Continuum Zirconia Bearing Head (XX-4101-XXX-XX) (12/14 taper)	ESKA BIOLOX Forte head (12/14 taper)	ESKA BIOLOX Forte head (8/10 taper)	NGK Zirconia femoral head (6 degree)	SULOX head (12/14 taper)	SULOX Head (8/10 taper)	Tete Alumine Cone 8/10 D28	Zirconia head (7676-XX-XXX) (12/14 taper)	BIOLOX delta head (00-8775-028/040-xx)(12/14 taper)	BIOLOX OPTION head system (00-8777-028/040-xx)(12/14 taper)	BIOLOX OPTION head system (00-8778-028/040-xx)(8/10 taper)
M/L Taper stem	✓	✓	✓	✓	✓							✓						
M/L Taper stem with Kinectiv Technology	Refer to Kinectiv Technology Page																	
MARBURG stem		✓	✓	✓	✓							✓				✓	✓	
Mayo 12/14 stem	✓	✓	✓	✓	✓											✓	✓	
Mayo 6 Deg stem											✓							
Mercure/F23 cementless stem	✓	✓	✓													✓	✓	
Mercure/F23 cemented stem	✓	✓	✓													✓	✓	
METABLOC stem		✓	✓	✓	✓							✓				✓	✓	
Minima cemented stem						✓							✓	✓				✓
Modular Austin Moore stem										✓								
Modular Calcar Replacement stem										✓								
Moore monoblock stem																		
MOST Options system		✓	✓	✓	✓							✓			✓	✓	✓	
MS-30 stem		✓	✓	✓	✓							✓				✓	✓	
Muller stem (XX-8050-051/056-XX)	✓	✓	✓									✓				✓	✓	
Müller Straight stem PROTASUL-10		✓	✓	✓	✓							✓				✓	✓	
Müller Straight stem PROTASUL-S30		✓	✓	✓	✓							✓				✓	✓	
Müller Straight stem PROTASUL-S30 (monoblock)																		
MultiLock stem											✓							
Natural-Hip LD stem																		
Natural-Hip Short Narrow Neck stem, CoCr				✓	✓							✓			✓	✓	✓	
Natural-Hip Short Narrow Neck stem, porous/titanium		✓	✓	✓	✓							✓			✓	✓	✓	
Natural-Hip System, CoCr												✓			✓	✓	✓	
Natural-Hip System, porous/titanium		✓	✓	✓	✓							✓			✓	✓	✓	
Natural-Hip Calcar Replacement stem												✓			✓	✓	✓	
Omnicasem stem (11.xx.xx.000/004)												✓			✓	✓	✓	
Omnicasem stem (01.060x.xxx)		✓	✓	✓	✓							✓			✓	✓	✓	
OPTAN stem		✓	✓	✓	✓							✓			✓	✓	✓	
Orphee uncemented HA stem												✓			✓	✓	✓	
PF Primary stem (11.xx.00.004/012 or 01.00411.xxx)												✓			✓	✓	✓	
PF Primary stem (01.060xx.xxx)		✓	✓									✓			✓	✓	✓	
PF Revision stem												✓			✓	✓	✓	
PFM-Revision stem		✓	✓	✓	✓							✓			✓	✓	✓	
Precedent Revision Hip system		✓	✓	✓	✓							✓			✓	✓	✓	
PYRAMID stem		✓	✓	✓	✓							✓			✓	✓	✓	
REVITAN PRIMARY stem *		✓	✓	✓	✓							✓			✓	✓	✓	
REVITAN straight/curved stem		✓	✓	✓	✓							✓			✓	✓	✓	
Segmental Proximal Femur	✓	✓	✓	✓	✓							✓			✓	✓	✓	
Segmental TM Proximal Femur	✓	✓	✓	✓	✓							✓			✓	✓	✓	
SF stem	✓	✓	✓												✓	✓	✓	
SKT stem	✓	✓	✓												✓	✓	✓	
SL Müller stem (cemented)		✓	✓	✓	✓							✓			✓	✓	✓	
STELCOR stem																		
TCM stem		✓	✓	✓	✓							✓			✓	✓	✓	
Thompson monoblock stem																		

Continued on Page 3

