

**Zimmer Ceramic Femoral Heads on
Biomet G7 BIOLOX delta Ceramic Liners**

- Indicates "Combination has been assessed as functionally compatible by Zimmer Biomet"
- Indicates "Combination has either not been assessed or is not functionally compatible"

Ceramic Heads All heads listed are compatible with the correspondingly sized liners with no limitations.		G7 BIOLOX delta Ceramic Liners
00-8775-028-01	BioloX Delta Head, 28MM, -3.5MM	✓
00-8775-028-02	BioloX Delta Head, 28MM, +0MM	✓
00-8775-028-03	BioloX Delta Head, 28MM, +3.5MM	✓
00-8775-032-01	BioloX Delta Head, 32MM, -3.5MM	✓
00-8775-032-02	BioloX Delta Head, 32MM, +0MM	✓
00-8775-032-03	BioloX Delta Head, 32MM, +3.5MM	✓
00-8775-032-04	BioloX Delta Head, 32MM, +7MM	✓
00-8775-036-01	BioloX Delta Head, 36MM, -3.5MM	✓
00-8775-036-02	BioloX Delta Head, 36MM, +0MM	✓
00-8775-036-03	BioloX Delta Head, 36MM, +3.5MM	✓
00-8775-036-04	BioloX Delta Head, 36MM, +7MM	✓
00-8775-040-01	BioloX Delta Head, 40MM, -3.5MM	✓
00-8775-040-02	BioloX Delta Head, 40MM, +0MM	✓
00-8775-040-03	BioloX Delta Head, 40MM, +3.5MM	✓
00-8775-040-04	BioloX Delta Head, 40MM, +7MM	✓
00-8777-028-01	BioloX Delta Option Head, 28MM, -3.0MM	✓
00-8777-028-02	BioloX Delta Option Head, 28MM, +0MM	✓
00-8777-028-03	BioloX Delta Option Head, 28MM, +3.5MM	✓
00-8777-028-04	BioloX Delta Option Head, 28MM, +7MM	✓
00-8777-032-01	BioloX Delta Option Head, 32MM, -3.0MM	✓
00-8777-032-02	BioloX Delta Option Head, 32MM, +0MM	✓
00-8777-032-03	BioloX Delta Option Head, 32MM, +3.5MM	✓
00-8777-032-04	BioloX Delta Option Head, 32MM, +7MM	✓
00-8777-036-01	BioloX Delta Option Head, 36MM, -3.0MM	✓
00-8777-036-02	BioloX Delta Option Head, 36MM, +0MM	✓
00-8777-036-03	BioloX Delta Option Head, 36MM, +3.5MM	✓
00-8777-036-04	BioloX Delta Option Head, 36MM, +7MM	✓
00-8777-040-01	BioloX Delta Option Head, 40MM, -3.0MM	✓
00-8777-040-02	BioloX Delta Option Head, 40MM, +0MM	✓
00-8777-040-03	BioloX Delta Option Head, 40MM, +3.5MM	✓
00-8777-040-04	BioloX Delta Option Head, 40MM, +7MM	✓
00-8778-028-01	BioloX Option Hd/Adpt, 8/10, 28 x -3.0	✓
00-8778-028-02	BioloX Option Hd/Adpt, 8/10, 28 x 0	✓
00-8778-028-03	BioloX Option Hd/Adpt, 8/10, 28 x +4.0	✓
00-8778-032-01	BioloX Option Hd/Adpt, 8/10, 32 x -3.0	✓
00-8778-032-02	BioloX Option Hd/Adpt, 8/10, 32 x 0	✓
00-8778-032-03	BioloX Option Hd/Adpt, 8/10, 32 x +4.0	✓
00-8778-036-01	BioloX Option Hd/Adpt, 8/10, 36 x -3.0	✓
00-8778-036-02	BioloX Option Hd/Adpt, 8/10, 36 x 0	✓
00-8778-036-03	BioloX Option Hd/Adpt, 8/10, 36 x +4.0	✓