

Femoral Stems	Longevity & Vivacit-E G7 Non-Constrained Liners	Longevity & Vivacit-E Dual Mobility Bearings
	✓	✓
Alloclassic SL Stem	✓	✓
Alloclassic SLL Revision Stem	✓	✓
Alloclassic Variall SLV Stem	✓	✓
Avenir Cemented Stems	✓	✓
Avenir Muller Cementless Stem	✓	✓
Avenir Complete Cementless Stem	✓	✓
CPT Hip Stem	✓	✓
Fitmore Primary Stem	✓	✓
M/L Taper Stem	✓	✓
M/L Taper Stem with Kinectiv Technology	✓	✓
TM Primary Hip Stem	✓	✓
VerSys Advocate Stem	✓	✓
VerSys Cemented Revision Calcar(CRC) Stem	✓	✓
VerSys Fiber Metal Midcoat Stem	✓	✓
VerSys Fiber Metal Porous Midcoat Stem	✓	✓
VerSys Fiber Metal Taper Stem	✓	✓
VerSys Fiber Metal Porous Taper Stem	✓	✓
VerSys Beaded Fullcoat Stem	✓	✓
VerSys Beaded Fullcoat Calcar Stem	✓	✓
Versys Heritage Stem	✓	✓
Versys Heritage Revision Stem	✓	✓
VerSys LD/Fx Cemented Stem	✓	✓
VerSys LD/Fx Press-Fit Stem	✓	✓
Wagner SL Revision Stem 135 NK Angle	✓	✓
Wagner SL Revision Stem 145 NK Angle	✓	✓
Wagner Revision Cones	✓	✓
CLS Brevius Stem with Kinectiv Technology	✓	✓
CLS Spotorno 125 Degree Stem	✓	✓
CLS Spotorno 135 Degree Stem	✓	✓
CLS Spotorno 145 Degree Stem	✓	✓
CLS Spotorno 145 Degree, HA Stem	✓	✓
Exafit Primary Stem	✓	✓
Exafit Revision Stem	✓	✓
PF Cemented Primary Stem	✓	✓
PF Cemented Revision Stem	✓	✓
PF Cementless Primary Stem	✓	✓
Revitan Revision Bodies and Stems	✓	✓
Weber Fix Cemented Stem	✓	✓
ZMR Revision Stem	✓	✓
Taperloc Complete Bonemaster Stem	✓	✓
Taperloc Complete Stem	✓	✓
Taperloc Bonemaster Reduced Distal Type 1 Stem	✓	✓
Taperloc Complete Bonemaster Reduced Distal 12/14 Stem	✓	✓
Taperloc Complete PPS Full Proximal Type 1 Stem	✓	✓
Taperloc Porous & HA 12/14 Stems	✓	✓

**G7 Liner/Bearing Compatibility with Zimmer Biomet
Femoral Stems and Heads**

Revised Date 12/31/2019
Page 2 of 10

Taperloc Femoral Stems 12/14 Porous - Std & Lat.	✓	✓
Taperloc Bonemaster Type 1 Stems	✓	✓
Taperloc Lateralized 12/14 Stems	✓	
Taperloc Bonemaster Type 1 Stems	✓	✓
Taperloc Standard Porous Type 1 Stems	✓	✓
Taperloc Lateralized Porous Type 1 Stems	✓	✓
Taperloc Bonemaster Porous Type 1 Stems	✓	✓
Taperloc PPS Reduced Distal Type 1 Stems	✓	✓
Taperloc Cemented 12/14 Stems	✓	✓
Taperloc Microplasty Stem	✓	✓
Taperloc HA Type 1 Stems	✓	
Taperloc Type 1 Lateralized RF Stems	✓	
Arcos Modular Revision Stem	✓	✓
Arcos Broach Bodies	✓	
Arcos Calcar Bodies	✓	
Arcos Cone Bodies	✓	
Arcos Interlocking PPS Stem	✓	
Arcos Proximal Taper Distal Stem	✓	✓
Arcos CYL Distal Stem		✓
Mallory/Head Lateralized HA Stem	✓	
Mallory/Head HA Type 1 Stems	✓	✓
Mallory/Head PPS Type 1 Stems	✓	✓
Mallory/ Head Bonemaster Type 1 Stems	✓	
BiMetric CoCr Type 1 Stems	✓	
BiMetric PPS Revision Stems	✓	✓
BiMetric PPS Non-Collared Revision Stems	✓	
BiMetric PPS XR Stems	✓	✓
BiMetric PPS XR Collared Stems	✓	✓
BiMetric HA XR Non-Collared Stems	✓	✓
BiMetric Cementless Porous Coated - Standard	✓	
BiMetric Cementless Porous Coated - Standard & BoneMaster	✓	
BiMetric Cementless Porous Coated - Standard with HA	✓	
BiMetric Cementless Porous Coated - Laterilised	✓	
BiMetric Cementless Porous Coated - Laterilised & Bonemaster	✓	
BiMetric Cementless Porous Coated - Laterilised with HA	✓	
BiMetric CoCr Laterilized Stem	✓	
BiMetric HAP Non-Collared Stem		✓
Bi-Metric CoCr Collared Cemented	✓	
Bi-Metric Collarless CoCr Cemented	✓	
Echo Bi-Metric Stem	✓	✓
Echo FX Stem	✓	✓
Echo PF Stem	✓	✓
Stanmore Stem	✓	✓
GTS Standard Stem	✓	✓
GTS Varized Stem	✓	✓
Integral Lateralized Stems	✓	✓
Integral Reduced Proximal Profile Stems	✓	
Integral Standard Offset Stems	✓	✓
Bio-Moore II Stem		✓
OSS Modular Stem	✓	✓
OSS Elliptical Stem	✓	✓
OSS Low Profile Proximal Stem		✓
Aura 12/14 Cemented Stem	✓	✓
Aura 12/14 PPS HA Stem	✓	✓
Sirius Cemented Stem	✓	✓
Austin Moore Modular Hip Stem	✓	
CMK 12/14 Revision Stem	✓	
Exception Cemented Stems	✓	
Exception Varaised Stems	✓	
Exception Standard Stems	✓	
Primoris Bonemaster Neck/Stem Prosthesis	✓	
Uption Stem	✓	
Metabloc Stem	✓	
MS-30 Hip Stem	✓	
M.E. Muller Stem	✓	

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

**G7 Liner/Bearing Compatibility with Zimmer Biomet
Femoral Stems and Heads**

Revised Date 12/31/2019
Page 3 of 10

CoCr Femoral Heads		Longevity & Vivacit-E G7 Non-Constrained Liners	Longevity & Vivacit-E Dual Mobility Bearings
		✓	✓
00-8018-028-01	12/14 COCR FEMORAL HEAD 28MM -3.5	✓	✓
00-8018-028-02	12/14 COCR FEMORAL HEAD 28MM +0	✓	✓
00-8018-028-03	12/14 COCR FEMORAL HEAD 28MM +3.5	✓	✓
00-8018-028-14	12/14 COCR FEMORAL HEAD 28MM +7	✓	✓
00-8018-028-05	12/14 COCR FEMORAL HEAD 28MM +10.5	✓	
00-8018-032-01	12/14 COCR FEMORAL HEAD 32MM -3.5	✓	
00-8018-032-02	12/14 COCR FEMORAL HEAD 32MM +0	✓	
00-8018-032-03	12/14 COCR FEMORAL HEAD 32MM +3.5	✓	
00-8018-032-14	12/14 COCR FEMORAL HEAD 32MM +7	✓	
00-8018-032-05	12/14 COCR FEMORAL HEAD 32MM +10.5	✓	
00-8018-036-01	12/14 COCR FEMORAL HEAD 36MM -3.5	✓	
00-8018-036-02	12/14 COCR FEMORAL HEAD 36MM +0	✓	
00-8018-036-03	12/14 COCR FEMORAL HEAD 36MM +3.5	✓	
00-8018-036-04	12/14 COCR FEMORAL HEAD 36MM +7	✓	
00-8018-036-05	12/14 COCR FEMORAL HEAD 36MM +10.5	✓	
00-8018-040-01	12/14 COCR FEMORAL HEAD 40MM -3.5	✓	
00-8018-040-02	12/14 COCR FEMORAL HEAD 40MM +0	✓	
00-8018-040-03	12/14 COCR FEMORAL HEAD 40MM +3.5	✓	
00-8018-040-04	12/14 COCR FEMORAL HEAD 40MM +7	✓	
00-8018-040-05	12/14 COCR FEMORAL HEAD 40MM +10.5	✓	
00-9026-028-00	6 DEGREE COCR FEMORAL HEAD 28MM +0 (S)	✓	
00-9026-028-35	6 DEGREE COCR FEMORAL HEAD 28MM +3.5 (S+)	✓	
00-9026-029-00	6 DEGREE COCR FEMORAL HEAD 28MM +7 (M)	✓	
00-9026-029-35	6 DEGREE COCR FEMORAL HEAD 28MM +10.5 (M+)	✓	
00-9026-030-00	6 DEGREE COCR FEMORAL HEAD 28MM +14 (L)	✓	
00-9026-032-00	6 DEGREE COCR FEMORAL HEAD 32MM +0 (S)	✓	
00-9026-032-35	6 DEGREE COCR FEMORAL HEAD 32MM +3.5 (S+)	✓	
00-9026-033-00	6 DEGREE COCR FEMORAL HEAD 32MM +7 (M)	✓	
00-9026-033-35	6 DEGREE COCR FEMORAL HEAD 32MM +10.5 (M+)	✓	
00-9026-034-00	6 DEGREE COCR FEMORAL HEAD 32MM +14 (L)	✓	
00-9026-036-00	6 DEGREE COCR FEMORAL HEAD 36MM +0 (S)	✓	
00-9026-036-03	6 DEGREE COCR FEMORAL HEAD 36MM +3.5 (S+)	✓	
00-9026-036-07	6 DEGREE COCR FEMORAL HEAD 36MM +7 (M)	✓	
00-9026-036-10	6 DEGREE COCR FEMORAL HEAD 36MM +10.5 (M+)	✓	
00-9026-036-14	6 DEGREE COCR FEMORAL HEAD 36MM +14 (L)	✓	
01.01012.288	12/14 COCR FEMORAL HEAD 28MM +8 DURASUL	✓	
01.01012.328	12/14 COCR FEMORAL HEAD 32MM +8 DURASUL	✓	
01.01012.365	12/14 COCR FEMORAL HEAD 36MM -4 DURASUL	✓	
01.01012.366	12/14 COCR FEMORAL HEAD 36MM +0 DURASUL	✓	
01.01012.367	12/14 COCR FEMORAL HEAD 36MM +4 DURASUL	✓	
01.01012.368	12/14 COCR FEMORAL HEAD 36MM +8 DURASUL	✓	
01.01012.444	12/14 COCR FEMORAL HEAD 44MM -8mm DURASUL	✓	
01.01012.445	12/14 COCR FEMORAL HEAD 44MM -4mm DURASUL	✓	
01.01012.446	12/14 COCR FEMORAL HEAD 44MM +0 DURASUL	✓	
01.01012.447	12/14 COCR FEMORAL HEAD 44MM +4mm DURASUL	✓	
01.01012.448	12/14 COCR FEMORAL HEAD 44MM +8mm DURASUL	✓	
14.28.05-20	12/14 COCR FEMORAL HEAD 28MM -4 (S) DURASUL	✓	✓
14.28.06-20	12/14 COCR FEMORAL HEAD 28MM +0 (M) DURASUL	✓	✓
14.28.07-20	12/14 COCR FEMORAL HEAD 28MM +4 (L) DURASUL	✓	✓

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

**G7 Liner/Bearing Compatibility with Zimmer Biomet
Femoral Stems and Heads**

Revised Date 12/31/2019

Page 4 of 10

14.32.05-20	12/14 COCR FEMORAL HEAD 32MM -4 (S) DURASUL	✓	✓
14.32.06-20	12/14 COCR FEMORAL HEAD 32MM +0 (M) DURASUL	✓	✓
14.32.07-20	12/14 COCR FEMORAL HEAD 32MM +4 (L) DURASUL	✓	✓
7210-28-000	12/14 COCR FEMORAL HEAD 28MM +0 DURASUL	✓	✓
7210-28-004	12/14 COCR FEMORAL HEAD 28MM -4 DURASUL	✓	✓
7210-28-400	12/14 COCR FEMORAL HEAD 28MM +4 DURASUL	✓	✓
7210-28-800	12/14 COCR FEMORAL HEAD 28MM +8 DURASUL	✓	✓
7210-32-000	12/14 COCR FEMORAL HEAD 32MM +0 DURASUL	✓	✓
7210-32-004	12/14 COCR FEMORAL HEAD 32MM -4 DURASUL	✓	✓
7210-32-400	12/14 COCR FEMORAL HEAD 32MM +4 DURASUL	✓	✓
7210-32-800	12/14 COCR FEMORAL HEAD 32MM +8 DURASUL	✓	✓
00-8018-022-02	12/14 COCR FEMORAL HEAD 22MM +0	✓	✓
00-8018-022-20	12/14 COCR FEMORAL HEAD 22MM -2	✓	✓
00-8018-022-30	12/14 COCR FEMORAL HEAD 22MM +3	✓	✓
01.01002.225	DURASUL CoCr Head 22.2 (-3.5, 0, +3.5) 8/10	✓	✓
01.01002.226	DURASUL CoCr Head 22.2 (-3.5, 0, +3.5) 8/10	✓	✓
01.01002.227	DURASUL CoCr Head 22.2 (-3.5, 0, +3.5) 8/10	✓	✓
01.01002.285	DURASUL CoCr Head 28 (-4, 0, +4) 8/10	✓	✓
01.01002.286	DURASUL CoCr Head 28 (-4, 0, +4) 8/10	✓	✓
01.01002.287	DURASUL CoCr Head 28 (-4, 0, +4) 8/10	✓	✓
10001007	CoCr Femoral Head 22.2mm Type 1 Std	✓	✓
10001008	CoCr Femoral Head 22.2mm Type 1 -3mm	✓	✓
10001009	CoCr Femoral Head 22.2mm Type 1 -5mm	✓	✓
10001010	COCR FEM HD 28MM TYPE 1 +12MM	✓	✓
10001011	COCR FEM HD 28MM TYPE 1 +9MM	✓	✓
10001012	COCR FEM HD 28MM TYPE 1 +6MM	✓	✓
10001013	COCR FEM HD 28MM TYPE 1 +3MM	✓	✓
10001014	COCR FEM HD 28MM TYPE 1 STD	✓	✓
10001015	COCR FEM HD 28MM TYPE 1 -3MM	✓	✓
10001016	COCR FEM HD 28MM TYPE 1 -6MM	✓	✓
10001017	COCR FEM HD 32MM TYPE 1 +12MM	✓	✓
10001018	COCR FEM HD 32MM TYPE 1 +9MM	✓	✓
10001019	COCR FEM HD 32MM TYPE 1 +6MM	✓	✓
10001020	COCR FEM HD 32MM TYPE 1 +3MM	✓	✓
10001021	COCR FEM HD 32MM TYPE 1 STD	✓	✓
10001022	COCR FEM HD 32MM TYPE 1 -3MM	✓	✓
10001023	COCR FEM HD 32MM TYPE 1 -6MM	✓	✓
10001024	COCR FEM HD 36MM TYPE 1 +12MM	✓	✓
10001025	COCR FEM HD 36MM TYPE 1 +9MM	✓	✓
10001026	COCR FEM HD 36MM TYPE 1 +6MM	✓	✓
10001027	COCR FEM HD 36MM TYPE 1 +3MM	✓	✓
10001028	COCR FEM HD 36MM TYPE 1 STD	✓	✓
10001029	COCR FEM HD 36MM TYPE 1 -3MM	✓	✓
10001030	COCR FEM HD 36MM TYPE 1 -6MM	✓	✓
10001031	COCR FEM HD 40MM TYPE 1 +12MM	✓	✓
10001032	COCR FEM HD 40MM TYPE 1 +9MM	✓	✓
10001033	COCR FEM HD 40MM TYPE 1 +6MM	✓	✓
10001034	COCR FEM HD 40MM TYPE 1 +3MM	✓	✓
10001035	COCR FEM HD 40MM TYPE 1 STD	✓	✓
10001036	COCR FEM HD 40MM TYPE 1 -3MM	✓	✓
10001037	COCR FEM HD 40MM TYPE 1 -6MM	✓	✓
10001038	COCR FEM HD 44MM TYPE 1 +12MM	✓	✓
10001039	COCR FEM HD 44MM TYPE 1 +9MM	✓	✓
10001040	COCR FEM HD 44MM TYPE 1 +6MM	✓	✓
10001041	COCR FEM HD 44MM TYPE 1 +3MM	✓	✓
10001042	COCR FEM HD 44MM TYPE 1 STD	✓	✓
10001043	COCR FEM HD 44MM TYPE 1 -3MM	✓	✓
10001044	COCR FEM HD 44MM TYPE 1 -6MM	✓	✓
10001045	COCR FEM HD 22.2MM TP 12/14 +2	✓	✓
10001046	COCR FEM HD 22.2MM TP 12/14STD	✓	✓
10001047	COCR FEM HD 22.2MM TP 12/14 -2	✓	✓
10001048	COCR FEM HD 28MM TP 12/14+10.5	✓	✓
10001049	COCR FEM HD 28MM TP 12/14 +7MM	✓	✓
10001050	COCR FEM HD 28MM TP 12/14 +3.5	✓	✓
10001051	COCR FEM HD 28MM TP 12/14 STD	✓	✓
10001052	COCR FEM HD 28MM TP 12/14 -3.5	✓	✓

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

**G7 Liner/Bearing Compatibility with Zimmer Biomet
Femoral Stems and Heads**

Revised Date 12/31/2019
Page 5 of 10

10001053	COCR FEM HD 32MM TP 12/14 +12	✓	
10001054	COCR FEM HD 32MM TP 12/14 +8MM	✓	
10001055	COCR FEM HD 32MM TP 12/14 +4MM	✓	
10001056	COCR FEM HD 32MM TP 12/14 STD	✓	
10001057	COCR FEM HD 32MM TP 12/14 -4MM	✓	
10001058	COCR FEM HD 36MM TP 12/14 +12	✓	
10001059	COCR FEM HD 36MM TP 12/14 +8MM	✓	
10001060	COCR FEM HD 36MM TP 12/14 +4MM	✓	
10001061	COCR FEM HD 36MM TP 12/14 STD	✓	
10001062	COCR FEM HD 36MM TP 12/14 -4MM	✓	
10001063	COCR FEM HD 40MM TP 12/14 +12	✓	
10001064	COCR FEM HD 40MM TP 12/14 +8MM	✓	
10001065	COCR FEM HD 40MM TP 12/14 +4MM	✓	
10001066	COCR FEM HD 40MM TP 12/14 STD	✓	
10001067	COCR FEM HD 40MM TP 12/14 -4MM	✓	
10001068	COCR FEM HD 44MM TP 12/14 +12	✓	
10001069	COCR FEM HD 44MM TP 12/14 +8MM	✓	
10001070	COCR FEM HD 44MM TP 12/14 +4MM	✓	
10001071	COCR FEM HD 44MM TP 12/14 STD	✓	
10001072	COCR FEM HD 44MM TP 12/14 -4MM	✓	
11-363660	36MM COCR MOD HD -6MM	✓	
11-363661	36MM COCR MOD HD -3MM	✓	
11-363662	36MM COCR MOD HD STD	✓	
11-363663	36MM COCR MOD HD +3MM	✓	
11-363664	36MM COCR MOD HD +6MM	✓	
11-363665	36MM COCR MOD HD +9MM	✓	
11-363666	36MM COCR MOD HD +12MM	✓	
157444	M2A-MAGNUM MOD HD SZ 44MM	✓	
163220	28MM OFFSET COCR MOD HD -3MMNK	✓	
163221	28MM OFFSET COCR MOD HD STD NK	✓	
163222	28MM OFFSET COCR MOD HD +3 NK	✓	
163223	28MM OFFSET COCR MOD HD +6 NK	✓	
163540	32MM OFFSET COCR MOD HD -6MMNK	✓	
163541	32MM OFFSET COCR MOD HD -3MMNK	✓	
163542	32MM OFFSET COCR MOD HD STD NK	✓	
163543	32MM OFFSET COCR MOD HD +3MMNK	✓	
163544	32MM OFFSET COCR MOD HD +6MMNK	✓	
163637	CO-CR 28MM MODULAR HEAD TAPER TYPE 1 +6MM/SPHERICAL		✓
163638	28MM COCR MOD HD +6MM NO SKIRT	✓	✓
163651	22.2MM DIA COCR MOD HD STD NK		✓
163652	22.2MM DIA COCR MOD HD -5 NK		✓
163653	22.2MM DIA COCR MOD HD -3MM NK		✓
163660	28MM DIA COCR MOD HD -6MM NK	✓	✓
163661	28MM DIA COCR MOD HD -3MM NK	✓	✓
163662	28MM DIA COCR MOD HD STD NK	✓	✓
163663	28MM DIA COCR MOD HD +3MM NK	✓	✓
163664	28MM DIA COCR MOD HD +6MM NK	✓	
163665	28MM DIA COCR MOD HD +9MM NK	✓	
163666	28MM DIA COCR MOD HD +12MM NK	✓	
163667	32MM MOD HEAD COCR -6MM NECK	✓	
163668	32MM MOD HEAD COCR -3MM NECK	✓	
163669	32MM MOD HEAD COCR STD NECK	✓	
163670	32MM MOD HEAD COCR +3MM NECK	✓	
163671	32MM MOD HEAD COCR +6MM NECK	✓	
163672	32MM MOD HEAD COCR +9MM NECK	✓	
163673	32MM MOD HEAD COCR +12MM NECK	✓	
163674	32MM COCR MOD HD +6MM NO SKIRT	✓	
164131	MOD HD 5DEG 42TAPER 22.2MM+2MM		✓
164132	MOD HD 5DEG 42TAPER 22.2MM 0MM		✓
164133	MOD HD 5DEG 42TAPER 22.2MM-2MM		✓
164440	22.2MM DIA COCR MOD HD STD NK		✓
164441	22.2MM DIA COCR MOD HD -5 NK		✓
650-0839	36MM COCR BIOMET FEM HD -6 NK	✓	
650-0840	36MM COCR BIOMET FEM HD -3 NK	✓	
650-0841	36MM COCR BIOMET FEM HD STD NK	✓	
650-0842	36MM COCR BIOMET FEM HD +3 NK	✓	

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

**G7 Liner/Bearing Compatibility with Zimmer Biomet
Femoral Stems and Heads**

Revised Date 12/31/2019

Page 6 of 10

650-0843	36MM COCR BIOMET FEM HD +6 NK	✓	
650-0870	32MM COCR BIOMET FEM HD -6 NK	✓	
650-0871	32MM COCR BIOMET FEM HD -3 NK	✓	
650-0872	32MM COCR BIOMET FEM HD STD NK	✓	
650-0873	32MM COCR BIOMET FEM HD +3 NK	✓	
650-0874	32MM COCR BIOMET FEM HD +6 NK	✓	
650-0877	28MM 12/14 TPR FEM HD -3.5 NK		✓
650-0878	28MM 12/14 TPR FEM HD STD NK		✓
650-0879	28MM 12/14 TPR FEM HD +3.5 NK		✓
650-0887	36MM 12/14 TPR FEM HD -4 NK	✓	
650-0888	36MM 12/14 TPR FEM HD STD NK	✓	
650-0889	36MM 12/14 TPR FEM HD +4 NK	✓	
650-0890	36MM 12/14 TPR FEM HD +8 NK	✓	
650-0891	36MM 12/14 TPR FEM HD +12 NK	✓	
650-1020	EXCEED ABT M2A-40 MOD HD -4MM	✓	
650-1021	EXCEED ABT M2A-40 MOD HD 0MM (	✓	
650-1022	EXCEED ABT M2A-40 MOD HD +4MM	✓	
650-1023	EXCEED ABT M2A-40 MOD HD +8MM	✓	
650-1024	EXCEED ABT M2A-40 MOD HD +12MM	✓	
650-1030	EXCEED ABT M2A-40 MOD HD -6MM	✓	
650-1031	EXCEED ABT M2A-40 MOD HD -3MM	✓	
650-1032	EXCEED ABT M2A-40 MOD HD 0MM (	✓	
650-1033	EXCEED ABT M2A-40 MOD HD +3MM	✓	
650-1034	EXCEED ABT M2A-40 MOD HD +6MM	✓	
650-1035	EXCEED ABT M2A-40 MOD HD +9MM	✓	
650-1036	EXCEED ABT M2A-40 MOD HD +12MM	✓	
P0206C22	COCR HEAD SHORT NECK Ø22.2 / -2		✓
P0206M22	COCR HEAD LONG NECK Ø22.2 / +0		✓
P0206L22	COCR HEAD LONG NECK Ø22.2 / +2		✓
P0206028	CRCO FEM HEAD EXTRA SHORT NECK		✓
P0206C28	COCR HD SHT NK 28 -3.5MM 12/14	✓	✓
P0206M28	COCR HD MD NK 28 STD 12/14	✓	✓
P0206L28	COCR HD LNG NK +3.5MM 12/14	✓	✓
00-8770-028-01	METASUL HEAD 28, 12/14, SIZE S		✓
00-8770-028-02	METASUL HEAD 28, 12/14, SIZE M		✓
00-8770-028-03	METASUL HEAD 28, 12/14, SIZE L		✓
00-8770-028-04	METASUL HEAD 28, 12/14, SIZE XL		✓
19.28.05	METASUL HEAD 28MM "S" 12/14		✓
19.28.06	METASUL HEAD 28MM "M" 12/14		✓
19.28.07	METASUL HEAD 28MM "L" 12/14		✓
11-163676	28MM M2A HI CARBON HD -6MM NK		✓
11-163677	28MM M2A HI CARBON HD -3MM NK		✓
11-163678	28MM M2A HI CARBON HD STD NK		✓
11-163679	28MM M2A HI CARBON HD +3MM NK		✓
11-163680	28MM M2A HI CARBON HD +6MM NK		✓
32-8018-028-01	FEM HD -3.5X28MM ION	✓	
32-8018-028-02	FEM HD +0X28MM ION	✓	
32-8018-028-03	FEM HD +3.5X28MM ION	✓	
32-8018-028-04	FEM HD +7X28MM ION	✓	
32-8018-028-05	FEM HD +10.5X28MM ION	✓	
00-8018-028-04	12/14 FEMORAL HEAD +7X28MM	✓	
650-0863	28MM COCR BIOMET FEM HD -6 NK	✓	✓
650-0864	28MM COCR BIOMET FEM HD -3 NK	✓	✓
650-0865	28MM COCR BIOMET FEM HD STD NK	✓	✓
650-0866	28MM COCR BIOMET FEM HD +3 NK	✓	✓
650-0867	28MM COCR BIOMET FEM HD +6 NK	✓	✓
802202201	UNIVERSAL COCR 22.2MM x -2		✓
802202202	UNIVERSAL COCR 22.2MM x +0		✓
802202203	UNIVERSAL COCR 22.2MM x +3		✓
802202801	UNIVERSAL COCR 28MM x -3.5	✓	✓
802202802	UNIVERSAL COCR 28MM x +0	✓	✓
802202803	UNIVERSAL COCR 28MM x +3.5	✓	✓
802202804	UNIVERSAL COCR 28MM x +7	✓	✓
BH400140	CO CR 40MM MOD HD 12/14 -6	✓	
BH400240	CO CR 40MM MOD HD 12/14 -3	✓	

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

**G7 Liner/Bearing Compatibility with Zimmer Biomet
Femoral Stems and Heads**

Revised Date 12/31/2019

Page 7 of 10

BH400340	CO CR 40MM MOD HD 12/14 STD	✓	✓
BH400440	CO CR 40MM MOD HD 12/14 +3	✓	✓
BH400540	CO CR 40MM MOD HD 12/14 +6	✓	✓
BH400640	CO CR 40MM MOD HD 12/14 +9	✓	✓
P0206E28	COCR HD XL NK 28 +7MM 12/14	✓	✓
P0206228	COCR MOD FEM HEAD Ø28/+10.5MM	✓	✓
P0206C32	COCR HD SHRT NK 32 -4MM 12/14	✓	✓
P0206M32	COCR HD MD NK 32 STD 12/14	✓	✓
P0206L32	COCR HD LNG NK 32 +4MM 12/14	✓	✓
P0206E32	COCR HD XL NK 32 +8MM 12/14	✓	✓
P0206232	COCR MOD FEM HEAD 32/+12	✓	✓
S001140	SELEX/MAGNUM MOD HD 40MM STD	✓	✓
S031140	SELEX/MAGNUM MOD HD 40MM -3	✓	✓
S061140	SELEX/MAGNUM MOD HD 40MM -6	✓	✓
S331140	SELEX/MAGNUM MOD HD 40MM +3	✓	✓
S661140	SELEX/MAGNUM MOD HD 40MM +6	✓	✓
S991140	SELEX/MAGNUM MOD HD 40MM +9	✓	✓

Ceramic Femoral Heads		Longevity & Vivacit-E G7 Non-Constrained Liners	Longevity & Vivacit-E Dual Mobility Bearings
00-8775-028-01	12/14 BIOLOX DELTA FEMORAL HEAD 28MM -3.5MM	✓	✓
00-8775-028-02	12/14 BIOLOX DELTA FEMORAL HEAD 28MM +0	✓	✓
00-8775-028-03	12/14 BIOLOX DELTA FEMORAL HEAD 28MM +3.5	✓	✓
00-8775-032-01	12/14 BIOLOX DELTA FEMORAL HEAD 32MM -3.5	✓	✓
00-8775-032-02	12/14 BIOLOX DELTA FEMORAL HEAD 32MM +0	✓	✓
00-8775-032-03	12/14 BIOLOX DELTA FEMORAL HEAD 32MM +3.5	✓	✓
00-8775-032-04	12/14 BIOLOX DELTA FEMORAL HEAD 32MM +7	✓	✓
00-8775-036-01	12/14 BIOLOX DELTA FEMORAL HEAD 36MM -3.5	✓	✓
00-8775-036-02	12/14 BIOLOX DELTA FEMORAL HEAD 36MM +0	✓	✓
00-8775-036-03	12/14 BIOLOX DELTA FEMORAL HEAD 36MM +3.5	✓	✓
00-8775-036-04	12/14 BIOLOX DELTA FEMORAL HEAD 36MM +7	✓	✓
00-8775-040-01	12/14 BIOLOX DELTA FEMORAL HEAD 40MM -3.5	✓	✓
00-8775-040-02	12/14 BIOLOX DELTA FEMORAL HEAD 40MM +0	✓	✓
00-8775-040-03	12/14 BIOLOX DELTA FEMORAL HEAD 40MM +3.5	✓	✓
00-8775-040-04	12/14 BIOLOX DELTA FEMORAL HEAD 40MM +7	✓	✓
00-8777-028-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM -3.0	✓	✓
00-8777-028-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +0	✓	✓
00-8777-028-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +3.5	✓	✓
00-8777-028-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +7	✓	✓
00-8778-028-01	8/10 BIOLOX OPTION Hd/Adpt 28MM -3.0	✓	✓
00-8778-028-02	8/10 BIOLOX OPTION Hd/Adpt 28MM +0	✓	✓
00-8778-028-03	8/10 BIOLOX OPTION Hd/Adpt 28MM +4.0	✓	✓
00-8777-032-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM -3.0	✓	✓
00-8777-032-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM +0	✓	✓
00-8777-032-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM +3.5	✓	✓
00-8777-032-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM +7	✓	✓
00-8777-036-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM -3.0	✓	✓
00-8777-036-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM +0	✓	✓
00-8777-036-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM +3.5	✓	✓
00-8777-036-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM +7	✓	✓
00-8777-040-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 40MM -3.0	✓	✓
00-8777-040-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 40MM +0	✓	✓

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

**G7 Liner/Bearing Compatibility with Zimmer Biomet
Femoral Stems and Heads**

Revised Date 12/31/2019
Page 8 of 10

00-8777-040-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 40MM +3.5	✓	✓
00-8778-032-03	8/10 BIOLOX OPTION Hd/Adpt 32MM +4.0	✓	✓
00-8778-036-01	8/10 BIOLOX OPTION Hd/Adpt 36MM -3.0	✓	✓
00-8778-036-02	8/10 BIOLOX OPTION Hd/Adpt 36MM +0	✓	✓
00-8778-036-03	8/10 BIOLOX OPTION Hd/Adpt 36MM +4.0	✓	✓
00-8778-040-01	8/10 BIOLOX OPTION Hd/Adpt 40MM -3.0	✓	✓
00-8778-040-02	8/10 BIOLOX OPTION Hd/Adpt 40MM +0	✓	✓
00-8778-040-03	8/10 BIOLOX OPTION Hd/Adpt 40MM +4.0	✓	✓
00-6428-028-01	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM -3.5	✓	✓
00-6428-028-02	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM +0	✓	✓
00-6428-028-03	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM +3.5	✓	✓
00-6428-032-01	12/14 ALUMINA CERAMIC FEMORAL HEAD 32MM -3.5	✓	✓
00-6428-032-02	12/14 ALUMINA CERAMIC FEMORAL HEAD 32MM +0	✓	✓
00-6428-032-03	12/14 ALUMINA CERAMIC FEMORAL HEAD 32MM +3.5	✓	✓
1010-29	CERAMIC ALU HEAD T12/14 28MM S	✓	✓
1010-31	CERAMIC ALU HEAD T12/14 28MM M	✓	✓
1010-33	CERAMIC ALU HEAD T12/14 28MM L	✓	✓
12.28.05	12/14 BIOLOX FORTE FEMORAL HEAD 28MM -3.5	✓	✓
12.28.06	12/14 BIOLOX FORTE FEMORAL HEAD 28MM +0	✓	✓
12.28.07	12/14 BIOLOX FORTE FEMORAL HEAD 28MM +3.5	✓	✓
12.32.05	12/14 BIOLOX FORTE FEMORAL HEAD 32MM -4	✓	✓
12.32.06	12/14 BIOLOX FORTE FEMORAL HEAD 32MM +0	✓	✓
12.32.07	12/14 BIOLOX FORTE FEMORAL HEAD 32MM +4	✓	✓
12-115109	DELTA CERAMIC FEM HD 028/-3MM	✓	✓
12-115110	DELTA CERAMIC FEM HD 028/ 0MM	✓	✓
12-115111	DELTA CERAMIC FEM HD 028/+3MM	✓	✓
12-115114	DELTA CERAMIC FEM HD Ø32/-3MM	✓	✓
12-115115	DELTA CERAMIC FEM HD Ø32/ 0MM	✓	✓
12-115116	DELTA CERAMIC FEM HD Ø32/+3MM	✓	✓
12-115117	DELTA CER FEM HD 32/+6MM T1	✓	✓
12-115120	DELTA CERAMIC FEM HD 036/-3MM	✓	✓
12-115121	DELTA CERAMIC FEM HD 036/0MM	✓	✓
12-115122	DELTA CERAMIC FEM HD 036/+3MM	✓	✓
12-115123	DELTA CERAMIC FEM HD 036/+6MM	✓	✓
P0209C28	ALU FEM HD D28/-3.5/5 42/12-14	✓	✓
P0209L28	ALU FEM HEAD D28/0/5 42/12-14	✓	✓
P0209M28	ALU FEM HD D28/+3.5/5 42/12-14	✓	✓
164135	DELTA CERAMIC FEM HD 028/-3MM	✓	✓
164136	DELTA CERAMIC FEM HD 028/ 0MM	✓	✓
164137	DELTA CERAMIC FEM HD 028/+3MM	✓	✓
164185	DELTA CERAMIC FEM HD Ø32/-3MM	✓	✓
164186	DELTA CERAMIC FEM HD Ø32/ 0MM	✓	✓
164187	DELTA CERAMIC FEM HD Ø32/+3MM	✓	✓
164190	BIOLOX F CER HD 28MM -3.5MM	✓	✓
164191	FORTE CERAMIC FEM HD Ø28/ 0MM	✓	✓
164192	FORTE CERAMIC FEM HD Ø28/+3.5M	✓	✓
650-0660	DELTA CERAMIC FEM HD 036/-3MM	✓	✓
650-0661	DELTA CERAMIC FEM HD 036/0MM	✓	✓
650-0662	DELTA CERAMIC FEM HD 036/+3MM	✓	✓
650-0663	DELTA CERAMIC FEM HD 036/+6MM	✓	✓
650-0667	DELTA CER FM HD 036/+8MM 12/14	✓	✓
650-0830	DELTA CER FM HD 028/-3.5 12/14	✓	✓
650-0831	DELTA CER FM HD 028/0MM 12/14	✓	✓
650-0832	DELTA CER FM HD 028/+3.5 12/14	✓	✓
650-0833	DELTA CER FM HD 032/-4.0 12/14	✓	✓
650-0834	DELTA CER FEM HD 032/0MM 12/14	✓	✓
650-0835	DELTA CER FM HD 032/+4MM 12/14	✓	✓
650-0836	DELTA CER FM HD 036/-4MM 12/14	✓	✓
650-0837	DELTA CER FM HD 036/0MM 12/14	✓	✓
650-0838	DELTA CER FM HD 036/+4MM 12/14	✓	✓
650-1055	CER BIOLOXD OPTION HD 28MM	✓	✓
650-1056	CER BIOLOXD OPTION HD 32MM	✓	✓
650-1057	CER BIOLOXD OPTION HD 36MM	✓	✓
650-1058	CER BIOLOXD OPTION HD 40MM	✓	✓
650-1059	BIOLOX DELTA CERAMIC OPTION HEAD 44MM	✓	✓
650-1157	DELTA CER FEM HD 28/+3MM T1	✓	✓

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

**G7 Liner/Bearing Compatibility with Zimmer Biomet
Femoral Stems and Heads**

Revised Date 12/31/2019
Page 9 of 10

650-1158	DELTA CER FEM HD 28/0MM T1	✓	✓
650-1159	DELTA CER FEM HD 28/-3MM T1	✓	✓
650-1160	DELTA CER FEM HD 32/+6MM T1	✓	✓
650-1161	DELTA CER FEM HD 32/+3MM T1	✓	✓
650-1162	DELTA CER FEM HD 32/0MM T1	✓	✓
650-1163	DELTA CER FEM HD 32/-3MM T1	✓	✓
650-1157	DELTA CER FEM HD 28/+3MM T1	✓	✓
650-1158	DELTA CER FEM HD 28/0MM T1	✓	✓
650-1159	DELTA CER FEM HD 28/-3MM T1	✓	✓
650-1060	CER OPTION 12/14 TPR SLEEVE -3	✓	✓
650-1061	CER OPTION 12/14 TPR SLEEVE +0	✓	✓
650-1062	CER OPTION 12/14 TPR SLEEVE +4	✓	✓
650-1063	CER OPTION 12/14 TPR SLEEVE +7	✓	✓
650-1064	CER OPTION TYPE 1 TPR SLEVE -6	✓	✓
650-1065	CER OPTION TYPE 1 TPR SLEVE -3	✓	✓
650-1066	CER OPT TYPE 1 TPR SLEVE 0MM	✓	✓
650-1067	CER OPTION TYPE 1 TPR SLEVE +3	✓	✓
650-1068	CER OPTION TYPE 1 TPR SLEVE +6	✓	✓
802802801	ZB Type 1 CERAMIC HD 28MM x -3	✓	✓
802802802	ZB Type 1 CERAMIC HD 28MM x +0	✓	✓
802802803	ZB Type 1 CERAMIC HD 28MM x +3	✓	✓
803002801	ZB SELECT CERAMIC HD 28MM	✓	✓
802602202	ZB 12/14 CERAMIC HD 22MM X +0	✓	✓
802602203	ZB 12/14 CERAMIC HD 22MM X +3	✓	✓
802602801	ZB 12/14 CERMC HD 28MM X -3.5	✓	✓
802602802	ZB 12/14 CERAMIC HD 28MM X +0	✓	✓
802602803	ZB 12/14 CERMC HD 28MM X +3.5	✓	✓
802802201	ZB T1 CERAMIC HD 22.2MM X -3	✓	✓
802802202	ZB T1 CERAMIC HD 22.2MM X +0	✓	✓
802802203	ZB T1 CERAMIC HD 22.2MM X +3	✓	✓
802802801	ZB TYPE 1 CERAMIC HD 28MM X -3	✓	✓
802802802	ZB TYPE 1 CERAMIC HD 28MM X +0	✓	✓
802802803	ZB TYPE 1 CERAMIC HD 28MM X +3	✓	✓

Stainless Steel Femoral Heads		Longevity & Vivacti-E G7 Non-Constrained Liners	Longevity & Vivacti-E Dual Mobility Bearings
00-8011-001-28	6 DEGREE CPT SST FEMORAL HEAD 28MM SHT NECK	✓	✓
00-8011-001-32	6 DEGREE CPT SST FEMORAL HEAD 32MM SHT NECK	✓	✓
00-8011-002-28	6 DEGREE CPT SST FEMORAL HEAD 28MM MED NECK	✓	✓
00-8011-002-32	6 DEGREE CPT SST FEMORAL HEAD 32MM MED NECK	✓	✓
00-8011-003-28	6 DEGREE CPT SST FEMORAL HEAD 28MM LNG NECK	✓	✓
00-8011-003-32	6 DEGREE CPT SST FEMORAL HEAD 32MM LNG NECK	✓	✓
30.32.05	12/14 PROTASUL S-30 FEMORAL HEAD 32MM -4 (S)	✓	✓
30.32.06	12/14 PROTASUL S-30 FEMORAL HEAD 32MM +0 (M)	✓	✓
30.32.07	12/14 PROTASUL S-30 FEMORAL HEAD 32MM +4 (L)	✓	✓
30.28.05	12/14 PROTASUL S-30 HEAD 28MM -4 (S)	✓	✓
30.28.06	12/14 PROTASUL S-30 HEAD 28MM +0 (M)	✓	✓
30.28.07	12/14 PROTASUL S-30 HEAD 28MM +4 (L)	✓	✓
01.01102.285	Protasul S30 Head 28 8/10 S	✓	✓
01.01102.286	Protasul S30 Head 28 8/10 M	✓	✓
01.01102.287	Protasul S30 Head 28 8/10 L	✓	✓

**G7 Liner/Bearing Compatibility with Zimmer Biomet
Femoral Stems and Heads**

Revised Date 12/31/2019
Page 10 of 10

P0201C22	STSTL FEM HD Ø22.2/-3.5MM (12/		✓
P0201M22	STSTL FEM HD D22.2/ 0MM 12/14		✓
P0201L22	STSTL FEM HD DIA22.2/+3.5MM 12		✓
P0201C28	STSTL FEM HD Ø28/-3.5MM (12/14	✓	✓
P0201M28	STSTL FEM HD Ø28/ 0MM (12/14)	✓	✓
P0201L28	STSTL FEM HD Ø28/+3.5MM (12/14	✓	✓
P0201C32	STSTL FEM HD Ø32/-4MM (12/14)	✓	
P0201E28	STSTL FEM HD Ø28/+7MM (12/14)	✓	
P0201E32	STAINLESS STEEL MODULAR HEAD 5° 42 - 12/14 DIA 32mm / +8mm	✓	
P0201L32	STSTL FEM HD Ø32/+4MM (12/14)	✓	
P0201M32	STSTL FEM HD Ø32/ 0MM (12/14)	✓	