

## Head and Stem Combinations: Unipolar and Bipolar Femoral Heads

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	<b>Unipolar Heads:</b>						<b>Bipolar Heads<sup>1</sup>:</b>							
	Unipolar Head (4510-XX-XXX) (12/14 taper)	Adapter for Unipolar Head (4510-XX-XXX) (12/14 taper)	Endo Femoral Head (6 degree)	Endo Femoral Head Adapter (6 degree)	Unipolar Head (non-modular) (30.38/60.09) (12/14 taper)	VerSys Endo Femoral Head (12/14 taper)	VerSys Endo Femoral Head Adapter (12/14 taper)	Bipolar Head (AP) (1946/1989)	Bi-polar head (4404-22/26-XXX)	Continuum Bipolar Head	Cupola Mobile	Cupule Mobile Retentive	Durasul Bi-polar head	Modular-bipolar head
ACA stem														
Acora cemented stem														
ALLOCLASSIC SLA and SL-64 stem <sup>2</sup>	✓	✓			✓		✓	✓		✓	✓	✓	✓	✓
ALLOCLASSIC SLL-Revision stem	✓	✓			✓		✓	✓		✓	✓	✓	✓	✓
ALLOCLASSIC stem <sup>2</sup>	✓	✓			✓		✓	✓		✓	✓	✓	✓	✓
ALLOCLASSIC VARIALL SLV stem					✓		✓	✓		✓	✓	✓	✓	✓
ALLOCOR AN stem					✓		✓	✓		✓	✓	✓	✓	✓
Anatomic stem			✓	✓								✓	✓	✓
Anatomic II stem						✓	✓					✓	✓	✓
Anatomique uncemented HA stem					✓		✓	✓		✓	✓	✓	✓	✓
Apollo Hip System	✓	✓			✓		✓	✓		✓	✓	✓	✓	✓
APR Hip System	✓	✓			✓		✓	✓		✓	✓	✓	✓	✓
APS Natural-Hip System						✓	✓					✓	✓	✓
Avenir Müller stem	✓	✓			✓	✓	✓			✓	✓	✓	✓	✓
Avenir Cemented stem	✓	✓			✓	✓	✓			✓	✓	✓	✓	✓
Avenir Complete stem	✓	✓			✓	✓	✓			✓	✓	✓	✓	✓
BEO stem					✓		✓	✓		✓	✓	✓	✓	✓
CDH stem PROTASUL-S30 (monoblock)							✓	✓		✓	✓	✓	✓	✓
CLS SPOTORNO stem	✓	✓			✓		✓	✓		✓	✓	✓	✓	✓
Concept cemented stem							✓	✓		✓	✓	✓	✓	✓
Continuum Cobrex stem									✓					
Continuum F-115 Hip stem									✓					
Continuum F-117 Collarless Hip stem									✓					
Continuum J-Taper stem									✓					
Continuum Proxilock stem (porous and HA coated)									✓					
CPS uncemented stem					✓		✓	✓		✓	✓	✓	✓	✓
CPT (6 Deg, CoCr) stem			✓	✓								✓	✓	✓
CPT (6 Deg, SST) stem								✓						✓
CPT 12/14 CoCr stem						✓	✓					✓	✓	✓
DCM-J Natural - Hip System						✓	✓					✓	✓	✓
DSP / TPP / PDC					✓		✓	✓		✓	✓	✓	✓	✓
Emeraude stem					✓		✓	✓		✓	✓	✓	✓	✓
Epoch stem						✓	✓					✓	✓	✓
ESKA 8/10 stem														
ESKA 12/14 stem														
Etoile cemented stem										✓	✓	✓	✓	✓
EXAFIT Primary stem										✓	✓	✓	✓	✓
EXAFIT Revision stem										✓	✓	✓	✓	✓
Fitmore Hip Stem	✓	✓			✓		✓	✓		✓	✓	✓	✓	✓
FracSure Hip system	✓	✓			✓		✓	✓		✓	✓	✓	✓	✓
Heritage stem (XX-8095-XXX-XX)			✓	✓								✓	✓	✓
Intermedia stem						✓	✓					✓	✓	✓
Lamella stem					✓		✓	✓		✓	✓	✓	✓	✓
M/L Taper stem						✓	✓					✓	✓	✓
M/L Taper stem with Kinectiv Technology														
MARBURG stem					✓		✓	✓		✓	✓	✓	✓	✓
Mayo 12/14 stem						✓	✓					✓	✓	✓
Mayo 6 Deg stem			✓	✓								✓	✓	✓
Mercure/F23 cementless stem						✓	✓					✓	✓	✓

1. See "Articulating Combinations" and "Head and Stem Combinations"  
 2. ALLOCLASSIC SLA and SL-64 (all sizes) and ALLOCLASSIC stems (sizes 01-6) can only be used with Bipolar heads (4404-22/26-XXX, 4405-22/26-XXX) up to size 56.

Continued on Page 2

## Head and Stem Combinations: Unipolar and Bipolar Femoral Heads

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	<b>Unipolar Heads:</b>						<b>Bipolar Heads<sup>1</sup>:</b>							
	Unipolar Head (4510-XX-XXX) (12/14 taper)	Adapter for Unipolar Head (4510-XX-XXX) (12/14 taper)	Endo Femoral Head (6 degree)	Endo Femoral Head Adapter (6 degree)	Unipolar Head (non-modular) (30.38/60.09) (12/14 taper)	VerSys Endo Femoral Head (12/14 taper)	VerSys Endo Femoral Head Adapter (12/14 taper)	Bipolar Head (AP) (1946/1989)	Bi-polar head (4404-22/26-XXX)	Continuum Bipolar Head	Cupola Mobile	Cupule Mobile Retentive	Durasul Bi-polar head	Modular-bipolar head
Mercure/F23 cemented stem														
METABLOC stem														
Minima cemented stem														
Modular Austin Moore stem														
Modular Calcar Replacement stem														
Moore monoblock stem														
MOST Options system														
MS-30 stem														
Muller stem (XX-8050-051/056-XX)														
Müller Straight stem PROTASUL-10														
Müller Straight stem PROTASUL-S30														
Müller Straight stem PROTASUL-S30 (monoblock)														
Multilock stem														
Natural-Hip LD stem														
Natural-Hip Short Narrow Neck stem														
Natural-Hip System														
Omnicares stem														
OPTAN stem														
Orphee uncemented HA stem														
PF Primary stem														
PF Revision stem														
PFM-Revision stem														
Precedent Revision Hip system														
PYRAMID stem														
REVITAN PRIMARY stem														
REVITAN straight/curved stem <sup>3</sup>														
Segmental Proximal Femur														
Segmental TM Proximal Femur														
SF stem														
SKT stem														
SL Müller stem (cemented)														
STELCOR stem														
TCM stem														
Trabecular Metal Primary Hip Stem														
Thompson monoblock stem														
VARIALL stem (except ALLOCLASSIC VARIALL SLV stem)														
VerSys Epoch RNL Stem														
VerSys Epoch Stem														
VerSys System <sup>4</sup>														
VIRTEC stem														
VIRTEC, Long stem														
Wagner Cone prosthesis stem														
Wagner SL Revision stem <sup>5</sup>														
Weber/Stühmer and Weber Curved stem														
ZMR (Body)														

1. See "Articulating Combinations" and "Head and Stem Combinations"
3. REVITAN straight stems in 14mm diameter cannot be used with femoral heads longer than +4mm.
4. VerSys Reduced Metaphyseal (RM) stems cannot be used with femoral heads longer than +7mm.
5. Wagner SL Revision stem in 14mm and 15mm diameters cannot be used with femoral heads longer than +4mm.