

**Zimmer Femoral Stems and Heads on Biomet Active Articulation
Acetabular Polyethylene Liners**

- Indicates "Combination has been assessed as functionally compatible by Zimmer Biomet"
- Indicates "Combination has either not been assessed or is not functionally compatible"

Femoral Stems

	Active Articulation ArCom Hip Bearing	Active Articulation E1 Hip Bearing
Alloclassic SL and SLO	✓	✓
Alloclassic SLL Revision Stem	✓	✓
Alloclassic Variall SLV Stem	✓	✓
Avenir Cemented	✓	✓
Avenir Muller Stem	✓	✓
Avenir Complete Stem	✓	✓
CPT 12/14	✓	✓
CPT 12/14 Long Stem	✓	✓
Fitmore Stem	✓	✓
M/L Taper Stem - Standard and Reduced Neck Length	✓	✓
M/L Taper with Kinectiv Technology	✓	✓
TM Primary Hip Stem	✓	✓
VerSys Advocate Stem	✓	✓
VerSys CEM/REV/CALCAR (CRC)	✓	✓
VerSys Fiber Metal Midcoat	✓	✓
VerSys Fiber Metal Taper	✓	✓
Versys Beaded Fullcoat (6in and 8in)	✓	✓
Versys Heritage stem	✓	✓
VerSys LD/Fx cemented	✓	✓
VerSys LD/Fx Press-Fit	✓	✓
Wagner SL Revision stem	✓	✓
CLS Brevius Stem with Kinectiv Technology	✓	✓
Exafit Stem	✓	✓
Exafit Revision Stem	✓	✓
PF Primary Stem	✓	✓
PF Revision Stem	✓	✓
Revitan Revision System	✓	✓
Weber Fix Cemented Stem	✓	✓
ZMR Revision System	✓	✓

CoCr Femoral Heads

		Active Articulation ArCom Hip Bearing	Active Articulation E1 Hip Bearing
00-8018-028-01	12/14 COCR FEMORAL HEAD 28MM -3.5	✓	✓
00-8018-028-02	12/14 COCR FEMORAL HEAD 28MM +0	✓	✓
00-8018-028-03	12/14 COCR FEMORAL HEAD 28MM +3.5	✓	✓
00-8018-028-14	12/14 COCR FEMORAL HEAD 28MM +7	✓	✓
00-8018-028-05	12/14 COCR FEMORAL HEAD 28MM +10.5	✓	✓
01.01012.288	12/14 COCR FEMORAL HEAD 28MM +8 DURASUL	✓	✓
14.28.05-20	12/14 COCR FEMORAL HEAD 28MM -4 (S) DURASUL	✓	✓
14.28.06-20	12/14 COCR FEMORAL HEAD 28MM +0 (M) DURASUL	✓	✓
14.28.07-20	12/14 COCR FEMORAL HEAD 28MM +4 (L) DURASUL	✓	✓
7210-28-000	12/14 COCR FEMORAL HEAD 28MM +0 DURASUL	✓	✓
7210-28-004	12/14 COCR FEMORAL HEAD 28MM -4 DURASUL	✓	✓
7210-28-400	12/14 COCR FEMORAL HEAD 28MM +4 DURASUL	✓	✓
7210-28-800	12/14 COCR FEMORAL HEAD 28MM +8 DURASUL	✓	✓
00-8018-022-02	12/14 COCR FEMORAL HEAD 22MM +0	✓	✓
00-8018-022-20	12/14 COCR FEMORAL HEAD 22MM -2	✓	✓
00-8018-022-30	12/14 COCR FEMORAL HEAD 22MM +3	✓	✓
01.01002.225	DURASUL CoCr Head 22.2 (-3.5, 0, +3.5) 8/10	✓	✓
01.01002.226	DURASUL CoCr Head 22.2 (-3.5, 0, +3.5) 8/10	✓	✓
01.01002.227	DURASUL CoCr Head 22.2 (-3.5, 0, +3.5) 8/10	✓	✓
01.01002.285	DURASUL CoCr Head 28 (-4, 0, +4) 8/10	✓	✓
01.01002.286	DURASUL CoCr Head 28 (-4, 0, +4) 8/10	✓	✓
01.01002.287	DURASUL CoCr Head 28 (-4, 0, +4) 8/10	✓	✓

**Zimmer Femoral Stems and Heads on Biomet Active Articulation
Acetabular Polyethylene Liners**

		Active Articulation	ArCom	Hip Bearing
Ceramic Femoral Heads				
00-8775-028-01	12/14 BIOLOX DELTA FEMORAL HEAD 28MM -3.5MM	✓		✓
00-8775-028-02	12/14 BIOLOX DELTA FEMORAL HEAD 28MM +0	✓		✓
00-8775-028-03	12/14 BIOLOX DELTA FEMORAL HEAD 28MM +3.5	✓		✓
00-8777-028-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM -3.0	✓		✓
00-8777-028-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +0	✓		✓
00-8777-028-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +3.5	✓		✓
00-8777-028-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +7	✓		✓
00-8778-028-01	8/10 BIOLOX OPTION Hd/Adpt 28MM -3.0	✓		✓
00-8778-028-02	8/10 BIOLOX OPTION Hd/Adpt 28MM +0	✓		✓
00-8778-028-03	8/10 BIOLOX OPTION Hd/Adpt 28MM +4.0	✓		✓
00-6428-028-01	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM -3.5	✓		✓
00-6428-028-02	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM +0	✓		✓
00-6428-028-03	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM +3.5	✓		✓
12.28.05	12/14 BIOLOX FORTE FEMORAL HEAD 28MM -3.5	✓		✓
12.28.06	12/14 BIOLOX FORTE FEMORAL HEAD 28MM +0	✓		✓
12.28.07	12/14 BIOLOX FORTE FEMORAL HEAD 28MM +3.5	✓		✓

		Active Articulation	ArCom	Hip Bearing
Stainless Steel Femoral Heads				
30.28.05	12/14 PROTASUL S-30 HEAD 28MM -4 (S)	✓		✓
30.28.06	12/14 PROTASUL S-30 HEAD 28MM +0 (M)	✓		✓
30.28.07	12/14 PROTASUL S-30 HEAD 28MM +4 (L)	✓		✓
01.01102.285	Protasul S30 Head 28 8/10 S	✓		✓
01.01102.286	Protasul S30 Head 28 8/10 M	✓		✓
01.01102.287	Protasul S30 Head 28 8/10 L	✓		✓
01.01102.225	Protasul S30 Head 22.2 8/10 S	✓		✓
01.01102.226	Protasul S30 Head 22.2 8/10 M	✓		✓
01.01102.227	Protasul S30 Head 22.2 8/10 L	✓		✓