

- Indicates "Combination has been assessed as functionally compatible by Zimmer Biomet"
- Indicates "Combination has either not been assessed or is not functionally compatible"

Femoral Stems	Exceed ABT Cemented
Alloclassic SL and SLO	✓
Alloclassic SLL Revision Stem	✓
Alloclassic Variall SLV Stem	✓
Avenir Cemented	✓
Avenir Muller Stem	✓
CLS Sportorno	✓
CPT 12/14	✓
CPT 12/14 Long Stem	✓
CPT 6 degree	✓
Fitmore Stem	✓
M/L Taper Stem - Standard and Reduced Neck Length	✓
M/L Taper with Kinectiv Technology	✓
Metabloc Cemented	✓
Metabloc Cementless	✓
MS-30 stem	✓
Original M.E. Muller Stem	✓
TM Primary Hip Stem	✓
VerSys Advocate Stem	✓
VerSys CEM/REV/CALCAR (CRC)	✓
VerSys Fiber Metal Midcoat	✓
VerSys Fiber Metal Taper	✓
Versys Beaded Fullcoat (6in 7.5in and 8in)	✓
Versys Heritage stem	✓
VerSys LD/Fx cemented	✓
VerSys LD/Fx Press-Fit	✓
Wagner Cone Prosthesis Stem	✓
Wagner SL Revision stem	✓
CLS Brevius Stem with Kinectiv Technology	✓
Exafit Stem	✓
Exafit Revision Stem	✓
PF Primary Stem	✓
PF Revision Stem	✓
Revitan Revision System	✓
Weber Fix Cemented Stem	✓
ZMR Revision System	✓

CoCr Femoral Heads	Exceed ABT Cemented	
00-8018-022-02	12/14 COCR FEMORAL HEAD 22MM +0	✓
00-8018-022-20	12/14 COCR FEMORAL HEAD 22MM -2	✓
00-8018-022-30	12/14 COCR FEMORAL HEAD 22MM +3	✓
00-8018-028-01	12/14 COCR FEMORAL HEAD 28MM -3.5	✓
00-8018-028-02	12/14 COCR FEMORAL HEAD 28MM +0	✓
00-8018-028-03	12/14 COCR FEMORAL HEAD 28MM +3.5	✓
00-8018-028-14	12/14 COCR FEMORAL HEAD 28MM +7	✓
00-8018-028-05	12/14 COCR FEMORAL HEAD 28MM +10.5	✓
00-8018-032-01	12/14 COCR FEMORAL HEAD 32MM -3.5	✓
00-8018-032-02	12/14 COCR FEMORAL HEAD 32MM +0	✓
00-8018-032-03	12/14 COCR FEMORAL HEAD 32MM +3.5	✓
00-8018-032-14	12/14 COCR FEMORAL HEAD 32MM +7	✓
00-8018-032-05	12/14 COCR FEMORAL HEAD 32MM +10.5	✓
00-8018-036-01	12/14 COCR FEMORAL HEAD 36MM -3.5	✓
00-8018-036-02	12/14 COCR FEMORAL HEAD 36MM +0	✓
00-8018-036-03	12/14 COCR FEMORAL HEAD 36MM +3.5	✓
00-8018-036-04	12/14 COCR FEMORAL HEAD 36MM +7	✓
00-8018-036-05	12/14 COCR FEMORAL HEAD 36MM +10.5	✓
00-8018-040-01	12/14 COCR FEMORAL HEAD 40MM -3.5	✓
00-8018-040-02	12/14 COCR FEMORAL HEAD 40MM +0	✓
00-8018-040-03	12/14 COCR FEMORAL HEAD 40MM +3.5	✓
00-8018-040-04	12/14 COCR FEMORAL HEAD 40MM +7	✓
00-8018-040-05	12/14 COCR FEMORAL HEAD 40MM +10.5	✓
00-9026-021-00	6 DEGREE COCR FEM HEAD 22MM +5	✓
00-9026-022-10	6 DEGREE COCR FEM HEAD 22MM +3.5	✓
00-9026-022-45	6 DEGREE COCR FEM HEAD 22MM +7	✓
00-9026-023-00	6 DEGREE COCR FEM HEAD 22MM +11 (L)	✓
00-9026-028-00	6 DEGREE COCR FEMORAL HEAD 28MM +0 (S)	✓
00-9026-028-35	6 DEGREE COCR FEMORAL HEAD 28MM +3.5 (S+)	✓
00-9026-029-00	6 DEGREE COCR FEMORAL HEAD 28MM +7 (M)	✓
00-9026-029-35	6 DEGREE COCR FEMORAL HEAD 28MM +10.5 (M+)	✓
00-9026-030-00	6 DEGREE COCR FEMORAL HEAD 28MM +14 (L)	✓
00-9026-032-00	6 DEGREE COCR FEMORAL HEAD 32MM +0 (S)	✓
00-9026-032-35	6 DEGREE COCR FEMORAL HEAD 32MM +3.5 (S+)	✓
00-9026-033-00	6 DEGREE COCR FEMORAL HEAD 32MM +7 (M)	✓
00-9026-033-35	6 DEGREE COCR FEMORAL HEAD 32MM +10.5 (M+)	✓
00-9026-034-00	6 DEGREE COCR FEMORAL HEAD 32MM +14 (L)	✓
00-9026-036-00	6 DEGREE COCR FEMORAL HEAD 36MM +0 (S)	✓
00-9026-036-03	6 DEGREE COCR FEMORAL HEAD 36MM +3.5 (S+)	✓
00-9026-036-07	6 DEGREE COCR FEMORAL HEAD 36MM +7 (M)	✓
00-9026-036-10	6 DEGREE COCR FEMORAL HEAD 36MM +10.5 (M+)	✓
00-9026-036-14	6 DEGREE COCR FEMORAL HEAD 36MM +14 (L)	✓
01.01002.225	DURASUL CoCr Head 22.2 (-3.5, 0, +3.5) 8/10	✓
01.01002.226	DURASUL CoCr Head 22.2 (-3.5, 0, +3.5) 8/10	✓
01.01002.227	DURASUL CoCr Head 22.2 (-3.5, 0, +3.5) 8/10	✓
01.01002.285	8/10 COCR FEMORAL HEAD 28MM -4 DURASUL	✓
01.01002.286	8/10 COCR FEMORAL HEAD 28MM +0 DURASUL	✓
01.01002.287	8/10 COCR FEMORAL HEAD 28MM +4 DURASUL	✓
01.01012.288	12/14 COCR FEMORAL HEAD 28MM +8 DURASUL	✓
01.01012.328	12/14 COCR FEMORAL HEAD 32MM +8 DURASUL	✓
01.01012.365	12/14 COCR FEMORAL HEAD 36MM -4 DURASUL	✓
01.01012.366	12/14 COCR FEMORAL HEAD 36MM +0 DURASUL	✓
01.01012.367	12/14 COCR FEMORAL HEAD 36MM +4 DURASUL	✓
01.01012.368	12/14 COCR FEMORAL HEAD 36MM +8 DURASUL	✓
01.01012.405	DURASUL CoCr Head 40, -4	✓
01.01012.406	DURASUL CoCr Head 40, 0	✓
01.01012.407	DURASUL CoCr Head 40, +4	✓
01.01012.408	DURASUL CoCr Head 40, +8	✓

CoCr Femoral Heads (Cont.)		Exceed ABT Cemented
14.28.05-20	12/14 COCR FEMORAL HEAD 28MM -4 (S) DURASUL	✓
14.28.06-20	12/14 COCR FEMORAL HEAD 28MM +0 (M) DURASUL	✓
14.28.07-20	12/14 COCR FEMORAL HEAD 28MM +4 (L) DURASUL	✓
14.32.05-20	12/14 COCR FEMORAL HEAD 32MM -4 (S) DURASUL	✓
14.32.06-20	12/14 COCR FEMORAL HEAD 32MM +0 (M) DURASUL	✓
14.32.07-20	12/14 COCR FEMORAL HEAD 32MM +4 (L) DURASUL	✓
7210-28-000	12/14 COCR FEMORAL HEAD 28MM +0 DURASUL	✓
7210-28-004	12/14 COCR FEMORAL HEAD 28MM -4 DURASUL	✓
7210-28-400	12/14 COCR FEMORAL HEAD 28MM +4 DURASUL	✓
7210-28-800	12/14 COCR FEMORAL HEAD 28MM +8 DURASUL	✓
7210-32-000	12/14 COCR FEMORAL HEAD 32MM +0 DURASUL	✓
7210-32-004	12/14 COCR FEMORAL HEAD 32MM -4 DURASUL	✓
7210-32-400	12/14 COCR FEMORAL HEAD 32MM +4 DURASUL	✓
7210-32-800	12/14 COCR FEMORAL HEAD 32MM +8 DURASUL	✓
05.95001.061	8/10 Durasul CoCr 32mm -4 (S)	✓
05.95001.062	8/10 Durasul CoCr 32mm +0 (M)	✓
05.95001.063	8/10 Durasul CoCr 32mm +4 (L)	✓
Ceramic Femoral Heads		Exceed ABT Cemented
00-8775-028-01	12/14 BIOLOX DELTA FEMORAL HEAD 28MM -3.5MM	✓
00-8775-028-02	12/14 BIOLOX DELTA FEMORAL HEAD 28MM +0	✓
00-8775-028-03	12/14 BIOLOX DELTA FEMORAL HEAD 28MM +3.5	✓
00-8775-032-01	12/14 BIOLOX DELTA FEMORAL HEAD 32MM -3.5	✓
00-8775-032-02	12/14 BIOLOX DELTA FEMORAL HEAD 32MM +0	✓
00-8775-032-03	12/14 BIOLOX DELTA FEMORAL HEAD 32MM +3.5	✓
00-8775-032-04	12/14 BIOLOX DELTA FEMORAL HEAD 32MM +7	✓
00-8775-036-01	12/14 BIOLOX DELTA FEMORAL HEAD 36MM -3.5	✓
00-8775-036-02	12/14 BIOLOX DELTA FEMORAL HEAD 36MM +0	✓
00-8775-036-03	12/14 BIOLOX DELTA FEMORAL HEAD 36MM +3.5	✓
00-8775-036-04	12/14 BIOLOX DELTA FEMORAL HEAD 36MM +7	✓
00-8775-040-01	12/14 BIOLOX DELTA FEMORAL HEAD 40MM -3.5	✓
00-8775-040-02	12/14 BIOLOX DELTA FEMORAL HEAD 40MM +0	✓
00-8775-040-03	12/14 BIOLOX DELTA FEMORAL HEAD 40MM +3.5	✓
00-8775-040-04	12/14 BIOLOX DELTA FEMORAL HEAD 40MM +7	✓
00-8777-028-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM -3.0	✓
00-8777-028-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +0	✓
00-8777-028-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +3.5	✓
00-8777-028-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +7	✓
00-8777-032-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM -3.0	✓
00-8777-032-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM +0	✓
00-8777-032-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM +3.5	✓
00-8777-032-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM +7	✓
00-8777-036-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM -3.0	✓
00-8777-036-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM +0	✓
00-8777-036-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM +3.5	✓
00-8777-036-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM +7	✓
00-8777-040-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 40MM -3.0	✓
00-8777-040-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 40MM +0	✓
00-8777-040-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 40MM +3.5	✓
00-8777-040-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 40MM +7	✓
00-8778-028-01	8/10 BIOLOX OPTION Hd/Adpt 28MM -3.0	✓
00-8778-028-02	8/10 BIOLOX OPTION Hd/Adpt 28MM +0	✓
00-8778-028-03	8/10 BIOLOX OPTION Hd/Adpt 28MM +4.0	✓
00-8778-032-01	8/10 BIOLOX OPTION Hd/Adpt 32MM -3.0	✓
00-8778-032-02	8/10 BIOLOX OPTION Hd/Adpt 32MM +0	✓
00-8778-032-03	8/10 BIOLOX OPTION Hd/Adpt 32MM +4.0	✓
00-8778-036-01	8/10 BIOLOX OPTION Hd/Adpt 36MM -3.0	✓
00-8778-036-02	8/10 BIOLOX OPTION Hd/Adpt 36MM +0	✓
00-8778-036-03	8/10 BIOLOX OPTION Hd/Adpt 36MM +4.0	✓

Ceramic Femoral Heads (Cont.)		Exceed ABT Cemented
00-6428-028-01	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM -3.5	✓
00-6428-028-02	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM +0	✓
00-6428-028-03	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM +3.5	✓
00-6428-032-01	12/14 ALUMINA CERAMIC FEMORAL HEAD 32MM -3.5	✓
00-6428-032-02	12/14 ALUMINA CERAMIC FEMORAL HEAD 32MM +0	✓
00-6428-032-03	12/14 ALUMINA CERAMIC FEMORAL HEAD 32MM +3.5	✓
12.28.05	12/14 BIOLOX FORTE FEMORAL HEAD 28MM -3.5	✓
12.28.06	12/14 BIOLOX FORTE FEMORAL HEAD 28MM +0	✓
12.28.07	12/14 BIOLOX FORTE FEMORAL HEAD 28MM +3.5	✓
12.32.05	12/14 BIOLOX FORTE FEMORAL HEAD 32MM -4	✓
12.32.06	12/14 BIOLOX FORTE FEMORAL HEAD 32MM +0	✓
12.32.07	12/14 BIOLOX FORTE FEMORAL HEAD 32MM +4	✓
17.28.05	12/14 Sulox 28mm -3.5	✓
17.28.06	12/14 Sulox 28mm +0	✓
17.28.07	12/14 Sulox 28mm +3.5	✓
17.32.05	12/14 Sulox 32mm -3.5	✓
17.32.06	12/14 Sulox 32mm +0	✓
17.32.07	12/14 Sulox 32mm +3.5	✓
01.01202.285	8/10 Sulox 28mm -3.5 (S)	✓
01.01202.286	8/10 Sulox 28mm +0 (M)	✓
01.01202.287	8/10 Sulox 28mm +3.5 (L)	✓
Stainless Steel Femoral Heads		Exceed ABT Cemented
00-8011-001-22	CPT MOD HD ST/ST 22MMD SHT NCK	✓
00-8011-002-22	CPT MOD HD ST/ST 22MED SHT NCK	✓
00-8011-003-22	CPT MOD HD ST/ST 22MMD LNG NCK	✓
00-8011-001-28	6 DEGREE CPT SST FEMORAL HEAD 28MM SHT NECK	✓
00-8011-001-32	6 DEGREE CPT SST FEMORAL HEAD 32MM SHT NECK	✓
00-8011-002-28	6 DEGREE CPT SST FEMORAL HEAD 28MM MED NECK	✓
00-8011-002-32	6 DEGREE CPT SST FEMORAL HEAD 32MM MED NECK	✓
00-8011-003-28	6 DEGREE CPT SST FEMORAL HEAD 28MM LNG NECK	✓
00-8011-003-32	6 DEGREE CPT SST FEMORAL HEAD 32MM LNG NECK	✓
01.01102.225	Protasul S30 Head 22.2 8/10 S	✓
01.01102.226	Protasul S30 Head 22.2 8/10 M	✓
01.01102.227	Protasul S30 Head 22.2 8/10 L	✓
01.01102.285	8/10 PROTASUL S-30 FEMORAL HEAD 28MM -4 (S)	✓
01.01102.286	8/10 PROTASUL S-30 FEMORAL HEAD 28MM +0 (M)	✓
01.01102.287	8/10 PROTASUL S-30 FEMORAL HEAD 28MM +4 (L)	✓
30.28.05	12/14 PROTASUL S-30 FEMORAL HEAD 28MM -4 (S)	✓
30.28.06	12/14 PROTASUL S-30 FEMORAL HEAD 28MM +0 (M)	✓
30.28.07	12/14 PROTASUL S-30 FEMORAL HEAD 28MM +4 (L)	✓
30.32.05	12/14 PROTASUL S-30 FEMORAL HEAD 32MM -4 (S)	✓
30.32.06	12/14 PROTASUL S-30 FEMORAL HEAD 32MM +0 (M)	✓
30.32.07	12/14 PROTASUL S-30 FEMORAL HEAD 32MM +4 (L)	✓