

Revised 6/1/2016
Page 1 of 5

*Sidus Stem-Free Shoulder Humeral Head 38-13
*Sidus Stem-Free Shoulder Humeral Head 40-14
*Sidus Stem-Free Shoulder Humeral Head 42-15
*Sidus Stem-Free Shoulder Humeral Head 44-16
*Sidus Stem-Free Shoulder Humeral Head 46-16
*Sidus Stem-Free Shoulder Humeral Head 48-17
*Sidus Stem-Free Shoulder Humeral Head 48-20
*Sidus Stem-Free Shoulder Humeral Head 50-18
*Sidus Stem-Free Shoulder Humeral Head 50-21
*Sidus Stem-Free Shoulder Humeral Head 52-19
*Sidus Stem-Free Shoulder Humeral Head 52-23
Anatomical Shoulder Humeral Head Ø36-12
Anatomical Shoulder Humeral Head Ø36-15
Anatomical Shoulder Humeral Head Ø38-13
Anatomical Shoulder Humeral Head Ø38-16
Anatomical Shoulder Humeral Head Ø40-14
Anatomical Shoulder Humeral Head Ø40-17
Anatomical Shoulder Humeral Head Ø42-15
Anatomical Shoulder Humeral Head Ø42-18
Anatomical Shoulder Humeral Head Ø44-16
Anatomical Shoulder Humeral Head Ø44-19
Anatomical Shoulder Humeral Head Ø46-16
Anatomical Shoulder Humeral Head Ø46-19
Anatomical Shoulder Humeral Head Ø48-17
Anatomical Shoulder Humeral Head Ø48-20
Anatomical Shoulder Humeral Head Ø50-18
Anatomical Shoulder Humeral Head Ø50-21
Anatomical Shoulder Humeral Head Ø52-19
Anatomical Shoulder Humeral Head Ø52-23
Anatomical Shoulder Humeral Head Ø54-20
Anatomical Shoulder Humeral Head Ø54-23
Anatomical Shoulder Humeral Head Ø56-21

| | | | | | | | | | | | | | | | | | | | | | |
|--|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | Comprehensive Modular Hybrid Glenoid Base (SM, MD, LG) |
| | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | Comprehensive Keeled Glenoid 4 mm (SM, MD, LG) |
| | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | Comprehensive Keeled Glenoid 7 mm (SM, MD, LG) |
| | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | Biglani/Flatow (BF) Glenoid w/ Keel 40mm, 46mmx40mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | BF Glenoid w/ Keel 46mm, 40mmx46mm, 52mmx46mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | BF Glenoid w/ Keel 52mm, 46mmx52mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | BF Glenoid w/ Keel 52mmx56mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | BF Glenoid w/ Pegs 40mm, 46mmx40mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | BF Glenoid w/ Pegs 46mm, 40mmx46mm, 52mmx46mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | BF Glenoid w/ Pegs 52mm, 46mmx52mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | BF Glenoid w/ Pegs 52mmx56mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | Trabecular Metal Glenoid 40mm, 46mmx40mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | Trabecular Metal Glenoid 46mm, 40mmx46mm, 52mmx46mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | Trabecular Metal Glenoid 52mm, 46mmx52mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | Trabecular Metal Glenoid 52mmx56mm Articular Surface |
| | | | | | | | | | | | | | | | | | | | | | Anatomical Shoulder Pegged Glenoid (Small) |
| | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | Anatomical Shoulder Pegged Glenoid (Medium) |
| | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | Anatomical Shoulder Pegged Glenoid (Large) |
| | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | Anatomical Shoulder keeled Glenoid (Small) |
| | | | | | | | | | | | | | | | | | | | | | Anatomical Shoulder keeled Glenoid (Medium) |
| | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | Anatomical Shoulder keeled Glenoid (Large) |



This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

Revised 6/1/2016
Page 2 of 5

Anatomical Shoulder Humeral Head Ø56-24
Anatomical Shoulder Removable Head D40 H14
Anatomical Shoulder Removable Head D42 H15
Anatomical Shoulder Removable Head D44 H16
Anatomical Shoulder Removable Head D46 H16
Anatomical Shoulder Removable Head D48 H17
Anatomical Shoulder Removable Head D50 H18
Anatomical Shoulder Removable Head D52 H19
Anatomical Shoulder Removable Head D48 H20
Anatomical Shoulder Removable Head D50 H21
Anatomical Shoulder Removable Head D52 H23
Anatomical Shoulder Fracture Slim head, Left ø36mm
Anatomical Shoulder Fracture Slim head, Right ø36mm
Anatomical Shoulder Fracture Slim head, Left ø40mm
Anatomical Shoulder Fracture Slim head, Right ø40mm
Anatomical Shoulder Fracture head, Left ø40mm
Anatomical Shoulder Fracture head, Right ø40mm
Anatomical Shoulder Fracture head, Left ø44mm
Anatomical Shoulder Fracture head, Right ø44mm
Anatomical Shoulder Fracture head, Left ø48mm
Anatomical Shoulder Fracture head, Right ø48mm
Anatomical Shoulder Fracture head, Left ø52mm
Anatomical Shoulder Fracture head, Right ø52mm
BF Standard Head 15mm x 40mm
BF Standard Head 18mm x 40mm
BF Standard Head 21mm x 40mm
BF Standard Head 24mm x 40mm
BF Standard Head 27mm x 40mm
BF Standard Head 15mm x 46mm
BF Standard Head 18mm x 46mm
BF Standard Head 21mm x 46mm
BF Standard Head 24mm x 46mm

[illegible]

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

Revised 6/1/2016
Page 3 of 5

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

Revised 6/1/2016
Page 4 of 5

Indicates "Combination has been assessed as functionally compatible by Zimmer Biomet"

Indicates "Combination has either not been assessed or is not functionally compatible"

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

Heads to Glenoid Combinations Total Shoulder System

Revised 6/1/2016
Page 5 of 5



Indicates "Combination has been assessed as functionally compatible by Zimmer Biomet"



Indicates "Combination has either not been assessed or is not functionally compatible"

| | Comprehensive Modular Hybrid Glenoid Base (SM, MD, LG) | | | | | | | | | | | | | | | | | |
|--|--|---|---|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
| | Comprehensive Keeled Glenoid 4 mm (SM, MD, LG) | | | | | | | | | | | | | | | | | |
| | Comprehensive Keeled Glenoid 7 mm (SM, MD, LG) | | | | | | | | | | | | | | | | | |
| | Bigliani/Flatow (BF) Glenoid w/ Keel 40mm, 46mmx40mm Articular Surface | | | | | | | | | | | | | | | | | |
| | BF Glenoid w/ Keel 46mm, 40mmx46mm, 52mmx46mm Articular Surface | | | | | | | | | | | | | | | | | |
| | BF Glenoid w/ Keel 52mm, 46mmx52mm Articular Surface | | | | | | | | | | | | | | | | | |
| | BF Glenoid w/ Keel 52mmx56mm Articular Surface | | | | | | | | | | | | | | | | | |
| | BF Glenoid w/ Pegs 40mm, 46mmx40mm Articular Surface | | | | | | | | | | | | | | | | | |
| | BF Glenoid w/ Pegs 46mm, 40mmx46mm, 52mmx46mm Articular Surface | | | | | | | | | | | | | | | | | |
| | BF Glenoid w/ Pegs 52mm, 46mmx52mm Articular Surface | | | | | | | | | | | | | | | | | |
| | BF Glenoid w/ Pegs 52mmx56mm Articular Surface | | | | | | | | | | | | | | | | | |
| | Trabecular Metal Glenoid 40mm, 46mmx40mm Articular Surface | | | | | | | | | | | | | | | | | |
| | Trabecular Metal Glenoid 46mm, 40mmx46mm, 52mmx46mm Articular Surface | | | | | | | | | | | | | | | | | |
| | Trabecular Metal Glenoid 52mm, 46mmx52mm Articular Surface | | | | | | | | | | | | | | | | | |
| | Trabecular Metal Glenoid 52mmx56mm Articular Surface | | | | | | | | | | | | | | | | | |
| | Anatomical Shoulder Pegged Glenoid (Small) | | | | | | | | | | | | | | | | | |
| | Anatomical Shoulder Pegged Glenoid (Medium) | | | | | | | | | | | | | | | | | |
| | Anatomical Shoulder Pegged Glenoid (Large) | | | | | | | | | | | | | | | | | |
| | Anatomical Shoulder keeled Glenoid (Small) | | | | | | | | | | | | | | | | | |
| | Anatomical Shoulder keeled Glenoid (Medium) | | | | | | | | | | | | | | | | | |
| | Anatomical Shoulder keeled Glenoid (Large) | | | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 46 x 18 x 53 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 46 x 21 x 50 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 46 x 24 x 47 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 46 x 27 x 46 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 50 x 21 x 57 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 50 x 24 x 52 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 50 x 27 x 50 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 54 x 21 x 64 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 54 x 24 x 58 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 54 x 27 x 55 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 58 x 24 x 64 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |
| Versa-Dial Titanium Humeral Head 58 x 27 x 61 mm | ✓ | ✓ | ✓ | | | | | | | | | | | | | | | |

* Sidus heads are currently going through clinical trials (IDE) and are not available in the US.