

✔ Indicates that the combination CAN be used  
✘ Indicates that the combination CANNOT be used

Screw Designation	1.5		2.0		2.7		3.5		4.0		4.5		6.5		Washer		4.0	
	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS
Note: This table represents functional compatibility and does not indicate the regulatory approval status. Be sure to check for proper Regulatory approvals in your country prior to using any products found on this website. Some devices may not be currently licensed with Health Canada. Some device compatibilities may not be approved for use by Health Canada. xx or xxx = These are placeholders that represent the item numbers for the range of device sizes. Note: Sterile Plates are export only. Not available in the United States.																		
1.5mm	1.5mm Condylar Plate (Left/Right) 47/00-4915-001/002-19	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	1.5mm H-Plate/Straight Plate 47/00-4915-004-20 47/00-4915-020-02	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	1.5mm T-Plate 47/00-4915-039/049-04	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
2.0mm	2.0mm Condylar Plate (Left/Right) 47/00-4920-001/002-19	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	2.0mm L-Plate (Left/Right) 47/00-4920-002-07	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	2.0mm H-Plate 47/00-4920-004-20	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	2.0mm Compression Plate 47/00-4920-0xx-00	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	2.0mm Straight Plate 47/00-4920-0xx-02	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
2.7mm Plate	2.7mm Condylar Plate (Left/Right) 47/00-4927-001/002-19	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	2.7mm L-Plate (Left/Right) 47/00-4927-003-xx	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	2.7mm Multi-Fragment Plate 47/00-4927-007-18	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	2.7mm Quarter Tubular Plate 47/00-4927-0xx-03	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	2.7mm Compression Plate 47/00-4927-0xx-00	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	2.7mm T-Plate 47/00-4927-003-04	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	2.7mm Straight Reconstruction Plate 47/00-1179-013-xx	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
3.5mm Plate	3.5mm Straight Reconstruction Plate 47/00-1179-005-xx	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	3.5mm Curved Reconstruction Plate 47/00-1179-006-xx	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	3.5mm Calcaneal Plate 47/00-4935-060/070/087-17	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	3.5mm T-Plate 47/00-4935-0xx-04/05	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	3.5mm Cloverleaf Plate 47/00-4935-0xx-16	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	3.5mm One-Third Tubular Plate 47/00-4935-0xx-03	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	3.5mm Contourable Dual Compression Plate 47/00-4935-0xx-07	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	3.5mm Std. Compression Plate 47/00-4935-0xx-00	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
4.5mm Plate	4.5mm Straight Reconstruction Plate 47/00-1179-007-02/16	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	4.5mm Symphyseal Bridge Plate 47/00-1179-002-02,04	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	4.5mm Broad Std Compression Plate 47/00-4945-005-026-01	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	4.5mm Cobra Head Plate 47/00-4945-0xx-14	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	4.5mm L-Butress Plate, 4 Hole 47/00-4945-004-11/12	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	4.5mm Narrow Std. Compression Plate 47/00-4945-0xx-00	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	4.5mm Contourable Dual Compression Plate, Narrow 47/00-4945-0xx-07	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	4.5mm Contourable Dual Compression Plate, Broad 47/00-4945-0xx-08	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	4.5mm Semi-Tubular Plate 47/00-4945-0xx-03	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	4.5mm Spoon Plate 47/00-4945-0xx-13	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	4.5mm T-Butress Plate 47/00-4945-0xx-06	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	4.5mm T-Plate 47/00-4945-0xx-04	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
Washer	2.7mm Spherical Washer 47/00-1179-015-00	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	3.5mm Spherical Washer 47/00-1179-012-00	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	Cortical Nut 47/00-4900-045-51	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	Large Washer 47/00-4900-065-51	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
	Mini Washer 47/00-4900-020-51	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔
Small Washer 47/00-4900-040-51	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	✔	