



**ZIMMER® A.T.S.®
DISPOSABLE
CUFFS**

Your Circle of
Confidence

LOW PROFILE
High Performance



LOW PROFILE *High Performance*

Zimmer Sterile Disposable Cuffs are designed and manufactured to ensure performance while providing single-use convenience. A comprehensive offering provides multiple options to help the clinician meet patient needs. The following features can help you create "Your Circle of Confidence."

Color Coding
For easy size identification.



Protection Sleeves
Specially designed to compliment Zimmer cuffs and to protect the patient from soft tissue damage.



90 Degree Ports
Help keep hoses clear of the surgical site.





Disposable

No cleaning, no cross contamination — just discard after use.

Single/Dual Bladder Options

Offers clinicians convenient choices for cuff style and sizing. Single bladder cuffs available in all sizes. Dual bladder cuffs available in selected sizes for use with intravenous regional anesthesia (IVRA or Bier Block) procedures.



Positive Locking Connectors*

Help prevent unintentional detachment from regulator hoses.

*Cuffs with Luer Lock connectors also available

High Quality Material

Latex-free and durable, manufactured to withstand any procedural demand.

Low Profile Design

Slim and lightweight for enhanced performance and patient comfort. Designed to provide less surgical obstruction.

Single/Dual Port Options

All cuff styles are available in both single and dual port configurations, allowing for use with any A.T.S. tourniquet regulator. Single port cuffs for use with the A.T.S. 750 and A.T.S. 1200 Tourniquet Systems. Dual port cuffs for use with A.T.S. 2000 Tourniquet System.

Contact Closures and Fasteners

Holds the cuff securely in position by virtually eliminating rolling and twisting.



Zimmer offers the most comprehensive line of cuffs on the market. The disposable cuff family provides high quality material without the hassle of cleaning or risk of cross contamination. With a wide variety of sizes, it is easy to meet even the most demanding procedural need.

A.T.S. Disposable Cuffs

- No cross contamination risk
- Durable, high quality material
- Single and dual bladder options
- Single and dual port options

A.T.S. Sterile Disposable Cuffs

DUAL PORT, SINGLE BLADDER [BOX 10]

Without Sleeve	With Sleeve	Size
60-7070-101	60-7070-151	8 in. (20cm)
60-7070-102	60-7070-152	12 in. (30cm)
60-7070-103	60-7070-153	18 in. (46cm)
60-7070-104	60-7070-154	24 in. (61cm)
60-7070-105	60-7070-155	30 in. (76cm)
60-7070-106	60-7070-156	34 in. (86cm)
60-7070-107	60-7070-157	42 in. (107cm)

DUAL PORT, DUAL BLADDER [BOX 10]

Without Sleeve	With Sleeve	Size
60-7080-101	60-7080-151	12 in. (30cm)
60-7080-102	60-7080-152	18 in. (46cm)
60-7080-103	60-7080-153	24 in. (61cm)

SINGLE PORT, SINGLE BLADDER [BOX 10]

Without Sleeve	With Sleeve	Size
60-7075-101	60-7075-151	8 in. (20cm)
60-7075-102	60-7075-152	12 in. (30cm)
60-7075-103	60-7075-153	18 in. (46cm)
60-7075-104	60-7075-154	24 in. (61cm)
60-7075-105	60-7075-155	30 in. (76cm)
60-7075-106	60-7075-156	34 in. (86cm)
60-7075-107	60-7075-157	42 in. (107cm)

SINGLE PORT, DUAL BLADDER [BOX 10]

Without Sleeve	With Sleeve	Size
60-7085-101	60-7085-151	12 in. (30cm)
60-7085-102	60-7085-152	18 in. (46cm)
60-7085-103	60-7085-153	24 in. (61cm)

SINGLE BLADDER CUFFS WITH SLEEVES [BOX 1]

Dual Port	Single Port	Size
60-7070-153-01	60-7075-153-01	18 in. (46cm)
60-7070-154-01	-	24 in. (61cm)
60-7070-155-01	-	30 in. (76cm)
60-7070-156-01	60-7075-156-01	34 in. (86cm)

Limb Protection Sleeves

- Available for all *Zimmer* tourniquet cuffs
- Specific designs for A.T.S. Cylindrical, Contour, and Small Adult Cuffs
- Shown to be more effective than other forms of limb protection
- Helps protect skin from wrinkling, pinching, and shearing
- Helps prevent tourniquet contamination
- Does not clog cuff closure material

Limb Protection Sleeves*

FOR A.T.S. DISPOSABLE AND REUSABLE CYLINDRICAL CUFFS [BOX 10]

Prod. No.	Description
60-8000-008	8 in. (20cm) Cuff Sleeve
60-8000-012	12 in. (30cm) Cuff Sleeve
60-8000-018	18 in. (46cm) Cuff Sleeve
60-8000-024	24 in. (61cm) Cuff Sleeve
60-8000-030	30 in. (76cm) Cuff Sleeve
60-8000-034	34 in. (86cm) Cuff Sleeve
60-8000-042	42 in. (107cm) Cuff Sleeve

*Sold Sterile

Esmark Bandages

LATEX FREE STERILE [BOX 20]

Prod. No.	Description
60-2694-011	3 in. x 9 ft. (8cm x 3m)
60-2694-012	4 in. x 9 ft. (10cm x 3m)
60-2694-013	5.5 in. x 9 ft. (14cm x 3m)
60-2694-014	3 in. x 12 ft. (8cm x 4m)
60-2694-015	4 in. x 12 ft. (10cm x 4m)
60-2694-016	5.5 in. x 12 ft. (14cm x 4m)

STERILE [BOX 20]

Prod. No.	Description
60-2694-001	3 in. x 9 ft. (8cm x 3m)
60-2694-002	4 in. x 9 ft. (10cm x 3m)
60-2694-003	6 in. x 9 ft. (15cm x 3m)
60-2694-004	3 in. x 12 ft. (8cm x 4m)
60-2694-005	4 in. x 12 ft. (10cm x 4m)
60-2694-006	6 in. x 12 ft. (15cm x 4m)

Contact your Zimmer representative or visit us at www.zimmer.com.

