



**ZIMMER®
A.T.S.® REUSABLE
CUFFS**

Your Circle of
Confidence

LOW PROFILE
High Performance



LOW PROFILE *High Performance*

Zimmer Reusable Cuffs are designed and manufactured to ensure superior convenience, comfort and safety in any tourniquet application. With the expanded line of Contour and Small Adult Cuffs, the clinician can feel secure that procedural demands will be met without sacrificing patient comfort and safety. The following features can help you create "Your Circle of Confidence."

90 Degree Ports
Help keep hoses clear of the surgical site.

Contour and Small Adult Cuffs*
Engineered to fit virtually any size or shape of limb, providing the ultimate in patient comfort and clinician confidence.

Protection Sleeves
Specially designed to compliment Zimmer cuffs and to protect the patient from soft tissue damage.





Reusable

Slim design means easy cleaning for comfortable wearing, time and time again.



Single/Dual Bladder Options

Offers clinicians convenient choices for cuff style and sizing. Single bladder cuffs available in all sizes. Dual bladder cuffs available in selected sizes for use with intravenous regional anesthesia (IVRA or Bier Block) procedures.



Positive Locking Connectors

Help prevent unintentional detachment from regulator hoses.

High Quality Material

Latex-free and durable, manufactured to withstand any procedural demand.

Low Profile Design

Slim and lightweight for enhanced performance and patient comfort. Designed to provide less surgical obstruction.



Single/Dual Port Options

All cuff styles are available in both single and dual port configurations, allowing for use with any A.T.S. tourniquet regulator. Single port cuffs for use with the A.T.S. 750 and A.T.S. 1200 Tourniquet Systems. Dual port cuffs for use with A.T.S. 2000 Tourniquet System.

Contact Closures and Fasteners

Holds the cuff securely in position by virtually eliminating rolling and twisting.

Zimmer offers the most comprehensive line of cuffs on the market:

A.T.S. Cylindrical Cuffs – Available in 7 sizes to accommodate a wide range of limb sizes. Individually color-coded for easy identification.

Contour Cuffs* – For use when standard A.T.S. Cylindrical Cuffs do not conform to actual limb shapes. Designed for superior fit and performance on contoured adult limbs. Unique design and adjustable fasteners allow a superior fit on a wide range of limb shapes and sizes.

Small Adult Cuffs* – For use when regular adult tourniquet cuffs are wider or longer than required for surgery. Unique P-shaped design minimizes tourniquet related obstruction of the surgical site. Comes in Extra Small (3 in. [8cm] wide) and Small (3.5 in. [9cm] wide).

* Manufactured by Delfi Medical Innovations, Inc.

Reusable Cuffs

- Available in A.T.S. Cylindrical, Contour or Small Adult styles
- Durable, high-quality material for longer wear
- Single and dual bladder options
- Single and dual port options

A.T.S. Cylindrical Cuffs

DUAL PORT, SINGLE BLADDER

Prod. No.	Description
60-7500-001	8 in. (20cm)
60-7500-002	12 in. (30cm)
60-7500-003	18 in. (46cm)
60-7500-004	24 in. (61cm)
60-7500-005	30 in. (76cm)
60-7500-006	34 in. (86cm)
60-7500-007	42 in. (107cm)

DUAL PORT, DUAL BLADDER

Prod. No.	Description
60-7555-001	12 in. (30cm)
60-7555-002	18 in. (46cm)
60-7555-003	24 in. (61cm)

SINGLE PORT, SINGLE BLADDER

Prod. No.	Description
60-7600-001	8 in. (20cm)
60-7600-002	12 in. (30cm)
60-7600-003	18 in. (46cm)
60-7600-004	24 in. (61cm)
60-7600-005	30 in. (76cm)
60-7600-006	34 in. (86cm)
60-7600-007	42 in. (107cm)

SINGLE PORT, DUAL BLADDER

Prod. No.	Description
60-7666-001	12 in. (30cm)
60-7666-002	18 in. (46cm)
60-7666-003	24 in. (61cm)

Contour Cuffs

DUAL PORT, SINGLE BLADDER

Prod. No.	Description
60-7450-002	Contour Upper Limb Cuff
60-7450-003	Contour Thigh Cuff
60-7450-004	Contour Lower Leg Cuff

DUAL PORT, DUAL BLADDER

Prod. No.	Description
60-7450-011	Contour Upper Limb IVRA Cuff

SINGLE PORT, SINGLE BLADDER

Prod. No.	Description
60-7350-002	Contour Upper Limb Cuff
60-7350-003	Contour Thigh Cuff
60-7350-004	Contour Lower Leg Cuff

SINGLE PORT, DUAL BLADDER

Prod. No.	Description
60-7350-011	Contour Upper Limb IVRA Cuff

Small Adult Cuffs

DUAL PORT, SINGLE BLADDER

Prod. No.	Description
60-8100-001	Extra Small Adult Cuff 3.0 in. (8cm) Wide
60-8100-002	Small Adult Cuff 3.5 in. (9cm) Wide

SINGLE PORT, SINGLE BLADDER

Prod. No.	Description
60-8000-001	Extra Small Adult Cuff 3.0 in. (8cm) Wide
60-8000-002	Small Adult Cuff 3.5 in. (9cm) Wide

Limb Protection Sleeves

- Specific designs for A.T.S. Cylindrical, Contour, and Small Adult Cuffs
- Shown to be more effective than other forms of limb protection
- Helps protect skin from wrinkling, pinching, and shearing
- Helps prevent tourniquet contamination
- Does not clog cuff closure material

Limb Protection Sleeves

FOR A.T.S. DISPOSABLE AND REUSABLE CYLINDRICAL CUFFS [BOX 10]*

Prod. No.	Description
60-8000-008	8 in. (20cm) Cuff Sleeve
60-8000-012	12 in. (30cm) Cuff Sleeve
60-8000-018	18 in. (46cm) Cuff Sleeve
60-8000-024	24 in. (61cm) Cuff Sleeve
60-8000-030	30 in. (76cm) Cuff Sleeve
60-8000-034	34 in. (86cm) Cuff Sleeve
60-8000-042	42 in. (107cm) Cuff Sleeve

* Sold Sterile

FOR CONTOUR CUFFS [BOX 10]**

Prod. No.	Description
60-7920-001	Contour Upper Limb Cuff Sleeve
60-7920-002	Contour Thigh Cuff Sleeve
60-7920-003	Contour Lower Leg Cuff Sleeve

** Sold Non-Sterile

FOR SMALL ADULT CUFFS [BOX 10]**

Prod. No.	Description
60-7910-001	3.0 in. (8cm) Extra Small Adult Cuff Sleeve
60-7910-002	3.5 in. (9cm) Small Adult Cuff Sleeve

** Sold Non-Sterile

Esmark Bandages

LATEX FREE STERILE [BOX 20]

Prod. No.	Description
60-2694-011	3 in. x 9 ft. (8cm x 3m)
60-2694-012	4 in. x 9 ft. (10cm x 3m)
60-2694-013	5.5 in. x 9 ft. (14cm x 3m)
60-2694-014	3 in. x 12 ft. (8cm x 4m)
60-2694-015	4 in. x 12 ft. (10cm x 4m)
60-2694-016	5.5 in. x 12 ft. (14cm x 4m)

STERILE [BOX 20]

Prod. No.	Description
60-2694-001	3 in. x 9 ft. (8cm x 3m)
60-2694-002	4 in. x 9 ft. (10cm x 3m)
60-2694-003	6 in. x 9 ft. (15cm x 3m)
60-2694-004	3 in. x 12 ft. (8cm x 4m)
60-2694-005	4 in. x 12 ft. (10cm x 4m)
60-2694-006	6 in. x 12 ft. (15cm x 4m)

Contact your Zimmer representative or visit us at www.zimmer.com.

