

Persona® OsseoTi® Keel Tibia
Persona® Cemented Tibia
Persona® PPS® Femur

Ordering Information



Table of Contents

Femurs

Cementless

Cruciate Retaining (CR) Porous Plasma Spray (PPS)	2-3
Posterior Stabilized (PS) Porous Plasma Spray (PPS)	3-4

Tibias

Cemented Tibias

Cemented Keel Tibias	7
----------------------------	---

Cementless Tibias

OsseoTi Keel Tibias	8
---------------------------	---

Persona Implant Dimensions

Persona PPS Femur Dimensions	10-11
Persona OsseoTi & Cemented Keel Tibia Dimensions	12

Instrument Ordering Information

Persona 0° Tibia Kit	14-15
Persona 0° Tibia Kit Size J Tray	16
Persona Spike Arm Tibia Resection Tray	17



**PERSONA
FEMORAL
COMPONENTS**



CR/ Left/ Narrow/ Porous Plasma Spray (PPS)

Part Number	Description
42-5080-050-01	Persona CR Femur, Porous Plasma Spray, Narrow, Size 1, Left
42-5080-052-01	Persona CR Femur, Porous Plasma Spray, Narrow, Size 2, Left
42-5080-054-01	Persona CR Femur, Porous Plasma Spray, Narrow, Size 3, Left
42-5080-056-01	Persona CR Femur, Porous Plasma Spray, Narrow, Size 4, Left
42-5080-058-01	Persona CR Femur, Porous Plasma Spray, Narrow, Size 5, Left
42-5080-060-01	Persona CR Femur, Porous Plasma Spray, Narrow, Size 6, Left
42-5080-062-01	Persona CR Femur, Porous Plasma Spray, Narrow, Size 7, Left
42-5080-064-01	Persona CR Femur, Porous Plasma Spray, Narrow, Size 8, Left
42-5080-066-01	Persona CR Femur, Porous Plasma Spray, Narrow, Size 9, Left
42-5080-068-01	Persona CR Femur, Porous Plasma Spray, Narrow, Size 10, Left
42-5080-070-01	Persona CR Femur, Porous Plasma Spray, Narrow, Size 11, Left

CR/ Right/ Narrow/ Porous Plasma Spray (PPS)

Part Number	Description
42-5080-050-02	Persona CR Femur, Porous Plasma Spray, Narrow, Size 1, Right
42-5080-052-02	Persona CR Femur, Porous Plasma Spray, Narrow, Size 2, Right
42-5080-054-02	Persona CR Femur, Porous Plasma Spray, Narrow, Size 3, Right
42-5080-056-02	Persona CR Femur, Porous Plasma Spray, Narrow, Size 4, Right
42-5080-058-02	Persona CR Femur, Porous Plasma Spray, Narrow, Size 5, Right
42-5080-060-02	Persona CR Femur, Porous Plasma Spray, Narrow, Size 6, Right
42-5080-062-02	Persona CR Femur, Porous Plasma Spray, Narrow, Size 7, Right
42-5080-064-02	Persona CR Femur, Porous Plasma Spray, Narrow, Size 8, Right
42-5080-066-02	Persona CR Femur, Porous Plasma Spray, Narrow, Size 9, Right
42-5080-068-02	Persona CR Femur, Porous Plasma Spray, Narrow, Size 10, Right
42-5080-070-02	Persona CR Femur, Porous Plasma Spray, Narrow, Size 11, Right



CR/ Left/ Standard/ Porous Plasma Spray (PPS)

Part Number	Description
42-5086-054-01	Persona CR Femur, Porous Plasma Spray, Standard, Size 3, Left
42-5086-056-01	Persona CR Femur, Porous Plasma Spray, Standard, Size 4, Left
42-5086-058-01	Persona CR Femur, Porous Plasma Spray, Standard, Size 5, Left
42-5086-060-01	Persona CR Femur, Porous Plasma Spray, Standard, Size 6, Left
42-5086-062-01	Persona CR Femur, Porous Plasma Spray, Standard, Size 7, Left
42-5086-064-01	Persona CR Femur, Porous Plasma Spray, Standard, Size 8, Left
42-5086-066-01	Persona CR Femur, Porous Plasma Spray, Standard, Size 9, Left
42-5086-068-01	Persona CR Femur, Porous Plasma Spray, Standard, Size 10, Left
42-5086-070-01	Persona CR Femur, Porous Plasma Spray, Standard, Size 11, Left
42-5086-074-01	Persona CR Femur, Porous Plasma Spray, Standard, Size 12, Left

CR/ Right/ Standard/ Porous Plasma Spray (PPS)

Part Number	Description
42-5086-054-02	Persona CR Femur, Porous Plasma Spray, Standard, Size 3, Right
42-5086-056-02	Persona CR Femur, Porous Plasma Spray, Standard, Size 4, Right
42-5086-058-02	Persona CR Femur, Porous Plasma Spray, Standard, Size 5, Right
42-5086-060-02	Persona CR Femur, Porous Plasma Spray, Standard, Size 6, Right
42-5086-062-02	Persona CR Femur, Porous Plasma Spray, Standard, Size 7, Right
42-5086-064-02	Persona CR Femur, Porous Plasma Spray, Standard, Size 8, Right
42-5086-066-02	Persona CR Femur, Porous Plasma Spray, Standard, Size 9, Right
42-5086-068-02	Persona CR Femur, Porous Plasma Spray, Standard, Size 10, Right
42-5086-070-02	Persona CR Femur, Porous Plasma Spray, Standard, Size 11, Right
42-5086-074-02	Persona CR Femur, Porous Plasma Spray, Standard, Size 12, Right



PS/ Left/ Narrow/ Porous Plasma Spray (PPS)

Part Number	Description
42-5060-050-01	Persona PS Femur, Porous Plasma Spray, Narrow, Size 1, Left
42-5060-052-01	Persona PS Femur, Porous Plasma Spray, Narrow, Size 2, Left
42-5060-054-01	Persona PS Femur, Porous Plasma Spray, Narrow, Size 3, Left
42-5060-056-01	Persona PS Femur, Porous Plasma Spray, Narrow, Size 4, Left
42-5060-058-01	Persona PS Femur, Porous Plasma Spray, Narrow, Size 5, Left
42-5060-060-01	Persona PS Femur, Porous Plasma Spray, Narrow, Size 6, Left
42-5060-062-01	Persona PS Femur, Porous Plasma Spray, Narrow, Size 7, Left
42-5060-064-01	Persona PS Femur, Porous Plasma Spray, Narrow, Size 8, Left
42-5060-066-01	Persona PS Femur, Porous Plasma Spray, Narrow, Size 9, Left
42-5060-068-01	Persona PS Femur, Porous Plasma Spray, Narrow, Size 10, Left
42-5060-070-01	Persona PS Femur, Porous Plasma Spray, Narrow, Size 11, Left

PS/ Right/ Narrow/ Porous Plasma Spray (PPS)

Part Number	Description
42-5060-050-02	Persona PS Femur, Porous Plasma Spray, Narrow, Size 1, Right
42-5060-052-02	Persona PS Femur, Porous Plasma Spray, Narrow, Size 2, Right
42-5060-054-02	Persona PS Femur, Porous Plasma Spray, Narrow, Size 3, Right
42-5060-056-02	Persona PS Femur, Porous Plasma Spray, Narrow, Size 4, Right
42-5060-058-02	Persona PS Femur, Porous Plasma Spray, Narrow, Size 5, Right
42-5060-060-02	Persona PS Femur, Porous Plasma Spray, Narrow, Size 6, Right
42-5060-062-02	Persona PS Femur, Porous Plasma Spray, Narrow, Size 7, Right
42-5060-064-02	Persona PS Femur, Porous Plasma Spray, Narrow, Size 8, Right
42-5060-066-02	Persona PS Femur, Porous Plasma Spray, Narrow, Size 9, Right
42-5060-068-02	Persona PS Femur, Porous Plasma Spray, Narrow, Size 10, Right
42-5060-070-02	Persona PS Femur, Porous Plasma Spray, Narrow, Size 11, Right




PS/ Left/ Standard/ Porous Plasma Spray (PPS)

Part Number	Description
42-5066-054-01	Persona PS Femur, Porous Plasma Spray, Standard, Size 3, Left
42-5066-056-01	Persona PS Femur, Porous Plasma Spray, Standard, Size 4, Left
42-5066-058-01	Persona PS Femur, Porous Plasma Spray, Standard, Size 5, Left
42-5066-060-01	Persona PS Femur, Porous Plasma Spray, Standard, Size 6, Left
42-5066-062-01	Persona PS Femur, Porous Plasma Spray, Standard, Size 7, Left
42-5066-064-01	Persona PS Femur, Porous Plasma Spray, Standard, Size 8, Left
42-5066-066-01	Persona PS Femur, Porous Plasma Spray, Standard, Size 9, Left
42-5066-068-01	Persona PS Femur, Porous Plasma Spray, Standard, Size 10, Left
42-5066-070-01	Persona PS Femur, Porous Plasma Spray, Standard, Size 11, Left
42-5066-074-01	Persona PS Femur, Porous Plasma Spray, Standard, Size 12, Left

PS/Right/ Standard/ Porous Plasma Spray (PPS)

Part Number	Description
42-5066-054-02	Persona PS Femur, Porous Plasma Spray, Standard, Size 3, Right
42-5066-056-02	Persona PS Femur, Porous Plasma Spray, Standard, Size 4, Right
42-5066-058-02	Persona PS Femur, Porous Plasma Spray, Standard, Size 5, Right
42-5066-060-02	Persona PS Femur, Porous Plasma Spray, Standard, Size 6, Right
42-5066-062-02	Persona PS Femur, Porous Plasma Spray, Standard, Size 7, Right
42-5066-064-02	Persona PS Femur, Porous Plasma Spray, Standard, Size 8, Right
42-5066-066-02	Persona PS Femur, Porous Plasma Spray, Standard, Size 9, Right
42-5066-068-02	Persona PS Femur, Porous Plasma Spray, Standard, Size 10, Right
42-5066-070-02	Persona PS Femur, Porous Plasma Spray, Standard, Size 11, Right
42-5066-074-02	Persona PS Femur, Porous Plasma Spray, Standard, Size 12, Right



PERSONA KEEL TIBIA BASEPLATES



Left Persona Cemented Keel Tibia

Part Number	Description
42-5360-058-01	Persona 0° Cemented Keel Tibia, Size A, Left
42-5360-061-01	Persona 0° Cemented Keel Tibia , Size B, Left
42-5360-064-01	Persona 0° Cemented Keel Tibia , Size C, Left
42-5360-067-01	Persona 0° Cemented Keel Tibia , Size D, Left
42-5360-071-01	Persona 0° Cemented Keel Tibia , Size E, Left
42-5360-075-01	Persona 0° Cemented Keel Tibia , Size F, Left
42-5360-079-01	Persona 0° Cemented Keel Tibia , Size G, Left
42-5360-083-01	Persona 0° Cemented Keel Tibia , Size H, Left
42-5360-088-01	Persona 0° Cemented Keel Tibia , Size J, Left

Right Persona Cemented Keel Tibia

Part Number	Description
42-5360-058-02	Persona 0° Cemented Keel Tibia , Size A, Right
42-5360-061-02	Persona 0° Cemented Keel Tibia , Size B, Right
42-5360-064-02	Persona 0° Cemented Keel Tibia , Size C, Right
42-5360-067-02	Persona 0° Cemented Keel Tibia , Size D, Right
42-5360-071-02	Persona 0° Cemented Keel Tibia , Size E, Right
42-5360-075-02	Persona 0° Cemented Keel Tibia , Size F, Right
42-5360-079-02	Persona 0° Cemented Keel Tibia , Size G, Right
42-5360-083-02	Persona 0° Cemented Keel Tibia , Size H, Right
42-5360-088-02	Persona 0° Cemented Keel Tibia , Size J, Right



Left Persona OsseoTi Keel Tibia


Part Number	Description
42-5350-058-01	Persona 0° OsseoTi Keel Tibia, Size A, Left
42-5350-061-01	Persona 0° OsseoTi Keel Tibia, Size B, Left
42-5350-064-01	Persona 0° OsseoTi Keel Tibia, Size C, Left
42-5350-067-01	Persona 0° OsseoTi Keel Tibia, Size D, Left
42-5350-071-01	Persona 0° OsseoTi Keel Tibia, Size E, Left
42-5350-075-01	Persona 0° OsseoTi Keel Tibia, Size F, Left
42-5350-079-01	Persona 0° OsseoTi Keel Tibia, Size G, Left
42-5350-083-01	Persona 0° OsseoTi Keel Tibia, Size H, Left
42-5350-088-01	Persona 0° OsseoTi Keel Tibia, Size J, Left

Right Persona OsseoTi Keel Tibia

Part Number	Description
42-5350-058-02	Persona 0° OsseoTi Keel Tibia, Size A, Right
42-5350-061-02	Persona 0° OsseoTi Keel Tibia, Size B, Right
42-5350-064-02	Persona 0° OsseoTi Keel Tibia, Size C, Right
42-5350-067-02	Persona 0° OsseoTi Keel Tibia, Size D, Right
42-5350-071-02	Persona 0° OsseoTi Keel Tibia, Size E, Right
42-5350-075-02	Persona 0° OsseoTi Keel Tibia, Size F, Right
42-5350-079-02	Persona 0° OsseoTi Keel Tibia, Size G, Right
42-5350-083-02	Persona 0° OsseoTi Keel Tibia, Size H, Right
42-5350-088-02	Persona 0° OsseoTi Keel Tibia, Size J, Right


PERSONA IMPLANT DIMENSIONS

Persona CR Femur (Narrow) Dimensions



Persona CR Femur (Narrow) Dimensions

Size	Overall A/P (mm)	Functional A/P (mm)	Overall M/L (mm)	Distal Thickness (mm)	Condyle Thickness (mm)
1*	48.9	41.0	55.5	8	8
2*	50.7	43.0	57.0	8	8
3	51.9	45.0	58.5	9	9
4	54.0	47.0	60.0	9	9
5	56.0	49.0	61.5	9	9
6	59.0	51.0	63.0	9	9
7	60.1	53.0	64.5	9	9
8	62.1	55.0	66.0	9	9
9	64.6	57.0	67.5	9	9
10	66.6	59.0	69.0	9	9
11	69.3	61.0	70.5	9	9




Persona CR Femur (Standard) Dimensions

Size	Overall A/P (mm)	Functional A/P (mm)	Overall M/L (mm)	Distal Thickness (mm)	Condyle Thickness (mm)
3	53.2	45.0	62.5	9	9
4	55.6	47.0	64.3	9	9
5	57.2	49.0	66.0	9	9
6	59.6	51.0	67.8	9	9
7	62	53.0	69.5	9	9
8	63.8	55.0	71.3	9	9
9	66.2	57.0	73.0	9	9
10	68.5	59.0	74.8	9	9
11	71.1	61.0	76.5	9	9
12	75.2	65.0	77.5	9	9

*Sizes 1-2 not available in Trabecular Metal® option


Cementless

PS Femur Dimensions – Cemented & Trabecular Metal



Persona PS Femur (Narrow) Dimensions

Size	Overall A/P (mm)	Functional A/P (mm)	Overall M/L (mm)	Distal Thickness (mm)	Condyle Thickness (mm)
1*	48.9	42.0	55.5	8	9
2*	51.2	44.0	57.0	8	9
3	52.7	46.0	58.5	9	10
4	55.0	48.0	60.0	9	10
5	57.3	50.0	61.5	9	10
6	59.0	52.0	63.0	9	10
7	60.8	54.0	64.5	9	10
8	63.5	56.0	66.0	9	10
9	65.7	58.0	67.5	9	10
10	67.9	60.0	69.0	9	10
11	70.0	62.0	70.5	9	10



Persona PS Femur (Standard) Dimensions

Size	Overall A/P (mm)	Functional A/P (mm)	Overall M/L (mm)	Distal Thickness (mm)	Condyle Thickness (mm)
3	54.2	46.0	62.5	9	10
4	56.6	48.0	64.3	9	10
5	58.7	50.0	66.0	9	10
6	60.4	52.0	67.8	9	10
7	63.0	54.0	69.5	9	10
8	65.0	56.0	71.3	9	10
9	67.2	58.0	73.0	9	10
10	69.4	60.0	74.8	9	10
11	71.7	62.0	76.5	9	10
12	75.9	66.0	77.5	9	10


*Sizes 1-2 not available in Trabecular Metal option

Persona Keel Tibia Dimensions




Persona OsseoTi Keel and Cemented Keel Tibia

Tibia Size	AB	CD	EF	GH & J
Keel Length	28 mm	40 mm	40 mm	50 mm
Peg Length*	8 mm	8 mm	10 mm	10 mm




Size	Medial A/P (mm)	Lateral A/P (mm)	Overall M/L (mm)
A	40.2	35.1	57.7
B	42.5	37.2	60.8
C	44.9	39.5	63.8
D	47.2	41.8	67.0
E	50.2	44.6	71.0
F	53.3	47.4	75.1
G	56.5	50.2	79.0
H	59.8	53.3	83.0
J	63.5	56.7	88.1

* Only on OsseoTi Keel Version










INSTRUMENTATION


Persona 0° Tibia Kit






Persona 0° Tibia Kit KT-5937-080-10 (PSK500)

Product	Part Number	Description	Bracket No.	Label
	42-5870-007-21	Persona Tibial Impactor Pad	110046271	1
	42-5399-018-21	Persona 0° Keel Tibial Drill	110046269	2
	42-5399-018-31	Persona 0° Keel Tibial Peg Drill	110046270	2
	42-5399-020-20	Persona 0° Keel Tibial Drill Guide	110046276	2
	42-5399-022-21	Persona 0° Keel Tibial Broach, Size A-B	110046264	3
	42-5399-022-23	Persona 0° Keel Tibial Broach, Size C-D	110046265	3
	42-5399-022-25	Persona 0° Keel Tibial Broach, Size E-F	110046266	3
	42-5399-022-27	Persona 0° Keel Tibial Broach, Size G-H	110046267	3
	42-5399-023-20	Persona 0° Keel Tibial Broach Insert/ Extractor Handle	110046261	4
	42-5399-058-21	Persona 0° Keel Tibial Sizing Plate, Size A, Left	110046272	5
	42-5399-058-22	Persona 0° Keel Tibial Sizing Plate, Size A, Right	110046272	5
	42-5399-061-21	Persona 0° Keel Tibial Sizing Plate, Size B, Left	110046272	5
	42-5399-061-22	Persona 0° Keel Tibial Sizing Plate, Size B, Right	110046272	5
	42-5399-064-21	Persona 0° Keel Tibial Sizing Plate, Size C, Left	110046273	5
	42-5399-064-22	Persona 0° Keel Tibial Sizing Plate, Size C, Right	110046273	5
	42-5399-067-21	Persona 0° Keel Tibial Sizing Plate, Size D, Left	110046273	5
	42-5399-067-22	Persona 0° Keel Tibial Sizing Plate, Size D, Right	110046273	5
	42-5399-071-21	Persona 0° Keel Tibial Sizing Plate, Size E, Left	110046274	5
	42-5399-071-22	Persona 0° Keel Tibial Sizing Plate, Size E, Right	110046274	5
	42-5399-075-21	Persona 0° Keel Tibial Sizing Plate, Size F, Left	110046274	5
	42-5399-075-22	Persona 0° Keel Tibial Sizing Plate, Size F, Right	110046274	5
	42-5399-079-21	Persona 0° Keel Tibial Sizing Plate, Size G, Left	110046275	5
	42-5399-079-22	Persona 0° Keel Tibial Sizing Plate, Size G, Right	110046275	5
	42-5399-083-21	Persona 0° Keel Tibial Sizing Plate, Size H, Left	110046275	5
	42-5399-083-22	Persona 0° Keel Tibial Sizing Plate, Size H, Right	110046275	5


Persona 0° Tibia Kit (cont.)




Persona 0° Tibia Kit KT-5937-080-10 (PSK500) (cont.)

Part Number	Part Number	Description	Bracket No.	Label
	42-5398-092-01	Persona Tibial Inserter Size A-D	00-5907-084-05	6
	42-5398-092-05	Persona Tibial Inserter Size E-F	00-5907-084-05	6
	42-5398-092-07	Persona Tibial Inserter Size G-J	00-5907-084-05	6
	110039004	Modular Tray Full Size		
	00-5908-071-00	Stainless Steel Single-Level Shell		
	00-5908-070-00	Lid for Stainless Steel Shell		
	00-5907-007-00	Pre-Lined Tray Side Artwork Panel		
	110038346	Removable Side Panel		
	00-5907-007-01	Persona Tray Side Artwork Panel		


Persona 0° Tibia Kit Size J






Persona 0° Tibia Kit Size J KT-5937-080-20 (PSKJT)

Product	Part Number	Description	Bracket No.	Label
	42-5399-022-29	Persona 0° Keel Tibial Broach, Size J	110046268	1
	42-5399-088-21	Persona 0° Keel Tibial Sizing Plate, Size J, Left	110038078	2
	42-5399-088-22	Persona 0° Keel Tibial Sizing Plate, Size J, Right	110038078	2
	00-5908-070-04	Quarter Size Lid		
	00-5908-071-04	Quarter Size Shell		
	00-5908-073-04	Quarter Size Tray		
	00-5907-007-00	Pre-Lined Tray Side Artwork Panel		
	00-5907-007-01	Persona Tray Side Artwork Panel		

Persona Spike Arm Tibia Cut Guide Tray



Persona Spike Arm Optional KT-5937-081-10 (PSK550)

Product	Part Number	Description	Bracket No.	Label
	42-5399-051-23	Persona Spiked Tibial Cut Guide for 3°	110046262	1
	42-5399-051-25	Persona Spiked Tibial Cut Guide for 5°	110046262	
	42-5399-051-27	Persona Spiked Tibial Cut Guide for 7°	110046262	
	42-5399-051-22	Persona 2° Varus/Valgus Spiked Tibial Cut Guide	110046262	2
	42-5399-001-20	Persona Spiked Tibial Alignment Guide	110046263	3
	110039004	Modular Tray Full Size		
	00-5908-071-00	Stainless Steel Single-Level Shell		
	00-5908-070-00	Lid for Stainless Steel Shell		
	00-5907-007-00	Pre-Lined Tray Side Artwork Panel		
	110038346	Removable Side Panel		
	00-5907-007-01	Persona Tray Side Artwork Panel		

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for health care professionals. Distribution to any other recipient is prohibited.

For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert or contact your local representative; visit www.zimmerbiomet.com for additional product information.

Check for country product clearances and reference product specific instructions for use. Not for distribution in France.

©2022 Zimmer Biomet



 **Legal Manufacturer**
Zimmer, Inc
1800 West Center Street
Warsaw, IN 46581-0708
USA

4037.1-GLBL-en-Issue Date 2022-12

zimmerbiomet.com