

Omni™ Surgical Suite System Release 1.0

1.0.13	Release notes - Omni Zimmer – Version 1.0.13
1.0.13	OMNIZ-734 OMNI HUB - SMS Message content update as approved by ICB
1.0.13	OMNIZ-621 OMNISight-User interface Omnilink is not fully functional
1.0.13	OMNIZ-554 Omnilink - inconsistent issue is getting freeze
1.0.13	OMNIZ-526 OmniLink - Virtual keyboard issue
1.0.13	OMNIZ-525 OmniSight - "Touch and Hold" over Text on OmniSight
1.0.13	OMNIZ-521 Checklist PDF header issue in V1.0.12
1.0.13	OMNIZ-514 Typo in messageconfig table for SMS that goes out to Surgeon
1.0.13	OMNIZ-509 Surgery Freezes when you click on begin button of manual surgery.
1.0.13	OMNIZ-766 OmniLink - updated AI (Object Recognition) label under Advanced setup

Known Issues

- With manual workflow, if the workflow is advanced from anesthesia start to surgery start by mistake and then reset to anesthesia start, the anesthesia timer will not reset. If it is advanced from surgery start to surgery stop by mistake and then reset to surgery start the surgery timer will not be reset.